

Annexure - VI

List of instruments with make and model

Sr No	Name of instrument	Make	Model	INST.ID. NO	Location
1	HPLC	Shimadzu	Prominance	CDTL/LC/02	F-125
2	HPLC	Agilent	1100	CDTL/LC/05	F-125
3	HPLC	Waters	2996	CDTL/LC/06	F-124
4	HPLC	Agilent	1200	CDTL/LC/07	F-125
5	HPLC	Agilent	1200	CDTL/LC/08	F-125
6	HPLC	Perkin Elmer	Flexar	CDTL/LC/09	G-13
7	HPLC	Perkin Elmer	Flexar	CDTL/LC/10	T-302
8	HPLC	Perkin Elmer	Flexar	CDTL/LC/11	T-302
9	HPLC	Perkin Elmer	Flexar	CDTL/LC/12	G-13
10	HPLC	Dionex	3000	CDTL/LC/13	G-13
11	HPLC	Dionex	Ultimate-3000	CDTL/LC/14	G-13
12	HPLC	Dionex	Ultimate-3000	CDTL/LC/15	G-13
13	HPLC	Dionex	Ultimate-3000	CDTL/LC/18	s-202
14	HPLC	Dionex	Ultimate-3000	CDTL/LC/19	s-202
15	HPLC	Dionex	Ultimate-3000	CDTL/LC/20	s-202
16	HPLC	Dionex	Ultimate-3000	CDTL/LC/21	T-302
17	HPLC	Dionex	Ultimate-3000	CDTL/LC/22	T-302
18	HPLC	Dionex	Ultimate-3000	CDTL/LC/23	T-302
19	HPLC	Agilent	1260	CDTL/LC/17	G-13
20	Fast HPLC	Dionex	Ultimate-3000	CDTL/LC/16	G-13
21	GC	ThermoFinnigan	Trace 1100	CDTL/Chem/GC/01	F-121
22	GC	Agilent	7890A	CDTL/Chem/GC/02	F-121
23	GC	ThermoFinnigan	Trace 1100	CDTL/Chem/GC/03	G-13
24	GC Gas Generator(TLC)	PCI	NAG-02M	CDTL//GG/01	F-120
25	GC Gas Generator(GC)	PCI	NAG-01 A	CDTL//GG/02	F-121
26	HPTLC	Camag	Linomat-V	CDTL/HPTLC/02	F-120
27	Linomate	Camag	NA	CDTL/HPTLC-LIN/01	
28	Linomate	Camag	NA	CDTL/HPTLC-LIN/02	
29	TLC(plate Heater)	Camag	NA	CDTL/TPH/01	F-120
30	TLC-Dryer	Camag	NA	CDTL/DR/01	F-120
31	Scanner	Camag	NA	CDTL/HPTLC-SCN/02	F-120
32	Visulizer	Camag	NA	CDTL/HPTLC-VIS/01	F-120
33	TLC-UV Cabinet	Camag	NA	CDTL/UVC/01	F-120
34	TLC-UV Cabinet	Camag	NA	CDTL/UVC/02	F-120
35	Plate Derivatisation	Camag	NA	CDTL/TLC/01	
36	TLC Plate development chamber	Camag	NA	CDTL/TSC/01	
37	Atomic Absorption Spectrophotometer	GBC	Savant AA	CDTL/AAS/2	F-122
38	Autotitrator	Metrohm	877 Titripot(702)	CDTL/AUTO/01	F-104
39	Autotitrator	Metrohm	888 Titrand	CDTL/AUTO/03	G-003
40	Autotitrator	Metrohm	888 Titrand	CDTL/AUTO/04	F-104
41	Autotitrator	Labindia	Titra smart	CDTL/AUTO/05	G00-3
42	Autotitrator	Labindia	Titra smart	CDTL/AUTO/06	G-003
43	KF	Metrohm	890KFT	CDTL/KF/02	F-106
44	KF	Labndia	KAFI smart	CDTL/KF/03	G-003
45	Coulometer	Metrohm	831KFC	CDTL/Coul/01	F-106
46	Autoclave	Osworld	NA	CDTL/AUTOCLAVE/03	S-225
47	Autoclave	Precision	PB-08	CDTL/AUTOCLAVE/04	S-225
48	Autoclave	Osworld	JRIC-39/A/B	CDTL/AUTOCLAVE/01	S-205
49	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/04	F-101
50	Analytical Balance	LCGC Radwag	AS220/C/2	CDTL/BAL/06	S-204
51	Analytical Balance	LCGC	377317/12	CDTL/BAL/07	S-201
52	Analytical Balance	Sartorius	BSA224S-CW	CDTL/BAL/11	S-220
53	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/15	F-118

54	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/16	F-106
55	Analytical Balance	Mettler toledo	ME-204	CDTL/BAL/18	G-003
56	Balance	Citizen	-----	CDTL/BAL/21	F-101
57	Balance	Sartorius	Secura12510IN	CDTL/BAL/22	F-108
58	Balance	Sartorius	Secura12510IN	CDTL/BAL/23	F-108
59	BET Kinetic instrument	Lonza	ELX-808	CDTL/KIN/01	S-217Micro
60	BOD Incubator	Osworld	3979	CDTL/BOD/01	S-223 Micro
61	BOD Incubator	Osworld	3980	CDTL/BOD/02	S-223 Micro
62	BOD Incubator	Metalab	JRIC-10	CDTL/BOD/03	S-223 Micro
63	BOD INCUBATOR	REMI		CDTL/BOD/03	S-223 Micro
64	Incubator	Metalab		CDTL/INC/01	S-223 Micro
65	Incubator	Metalab		CDTL/INC/02	S-223 Micro
66	Biosafety Cabinet			CDTL/BSC/01	S-221 micro
67	Loop sterilizer			CDTL/ELS/01	Micro
68	Centrifuge	Electrolab	CM-6MT	CDTL/CENTRI/01	F-101
69	Centrifuge	Remi motors	R-8C BL	CDTL/CENTRI/04	F-104
70	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/05	S-201
71	Centrifuge	Remi motors	R-8C PLUS	CDTL/CENTRI/06	T-301
72	Dehumidifier	Bry-Air	NA	CDTL/DH/02	F-118
73	Dehumidifier	sharp	NA	CDTL/DH/03	G-003
74	Dehumidifier	sharp		CDTL/DH/04	F-104
75	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/05	F-119
76	Disintegration apparatus	Electrolab	EDI-2SA	CDTL/DT/06	T-304
77	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/01	F-119
78	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/02	F-119
79	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/03	S-204
80	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/04	G-04
81	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/05	G-04
82	Dissolution apparatus	LabIndia	DS8000	CDTL/Disso/06	G-04
83	Dissolution apparatus-Autosampler	LabIndia	DS8000	CDTL/Disso/07	T-305
84	Dissolution apparatus-Autosampler	LabIndia	DS8000	CDTL/Disso/08	G-04
85	Dissolution apparatus-Autosampler	LabIndia	DS8001	CDTL/Disso/09	S-203
86	Dissoution media Prep systm	Lab inida	Media MET+	CDTL/DMP/01	G-04
87	Digital Flow meter for HPLC	Stainmat Scientific		CDTL/FM/02	F-125
88	FTIR	Shimadzu	IR Affinity - 1	CDTL/ /FTIR/02	F-118
89	FTIR (ATR + Transmission module)	Bruker	Alpha-II	CDTL/FTIR/03	F-118
90	Dry box for FTIR	Technosearch eng		CDTL//DB/02	F-118
91	FTIR-Press	Technosearch		CDTL/ /RP/03	F-118
92	Fume Hood	Labguard	ST1500ABP	CDTL/FH/01	F-107
93	Fume Hood	Labguard	ST1500ABP	CDTL/FH/02	S-205
94	Fume Hood	Labguard	ST1500ABP	CDTL/FH/03	S-205
95	Fume Hood	Prama tech	NA	CDTL/FH/04	T-305
96	Funnel shaker for solvent extractor	Jiotech	ST1500ABP	CDTL/SFS/01	T-303
97	Heating Block	Neo lab	NA	CDTL/HB/01	S-217
98	Heating Block	Neo lab	NA	CDTL/HB/02	S-218
99	Heating mentle	Supertech		CDTL/HM/11	Micro
100	Laminar Air Flow	Air borne	NA	CDTL/MIC/LAF/01	S-219 micro
101	Laminar Air Flow	Air borne	NA	CDTL/MIC/LAF/02	S-220 micro
102	Laminar Air Flow	Air borne	NA	CDTL/MIC/LAF/03	S-221 micro
103	Laminar Air Flow	Air borne	NA	CDTL/MIC/LAF/04	S-222 micro
104	Melting point apparatus	LabIndia	MR-VIS	CDTL/MP/01	F-117
105	Melting point apparatus	LabIndia	MR-VIS	CDTL/MP/02	T-303
106	Microscope	Vision2000	NA	CDTL/MIC/01	S-217
107	Microscope	Labomed	NA	CDTL/MIC/03	S-217
108	Muffle furnance	NA	NA	CDTL/MF/01	T-305
109	Muffle furnance	NA	NA	CDTL/MF/02	S-205

110	Oven	PATHAK	NA	CDTL/OV/02	F-104
111	Oven	MIC	NA	CDTL/OV/04	F-102
112	Oven	Pathak	NA	CDTL/OV/05	S-204
113	Oven			CDTL/OV/07	225-MICRO
114	Oven	MODERN	NA	CDTL/OV/08	S-205
115	Oven			CDTL/OV/09	S.M. Scientific
116	Oven			CDTL/OV/10	S.M. Scientific
117	Oil Bath	Metalab	NA	CDTL/OB/01	T-303
118	pH Meter	LabIndia	Pico Plus	CDTL/pH/03	F-117
119	pH Meter	Euotech inst	Ph 700	CDTL/pH/09	F-101
120	pH Meter	Horiba	F72A	CDTL/pH/13	F-117
121	pH- meter	Eu-tech	pH-510	CDTL/pH/04	225 micro
122	pH- Conductometer	Horiba	PC2000	CDTL/pH/15	T-301
123	Polarimeter	Anton-Par		CDTL/POL/02	F-117
124	Polarimeter	Anton-Par		CDTL/POL/03	201
125	Profile projector	SIPCON	SP-300	CDTL/PP/02	T-303
126	Particle counter	Shreedhar Instrumrnts		CDTL/PC/01	micro
127	Refractometer	Atago	RX-7000i	CDTL/RI/02	F-117
128	Rotary shaker	Remi	RS-24BI	CDTL/SH/01	S-203
129	Rotary shaker	Remi		CDTL/SH/02	G-010
130	Sonicattor	MIC		CDTL//Soni/04	T-301
131	Sonicattor	S-clean	NA	CDTL//Soni/11	S-201
132	Sonicattor	Centrofix	SM-100-US	CDTL/SONI/13	F-104
133	Sonicattor			CDTL/SONI/15	G-003
134	Sonicattor	S-clean		CDTL/SONI/12	F-101
135	Sonicattor	Lab devices		CDTL/SONI/08	G-003
136	Chiller Sonicator	S clean	SM200 CH US	CDTL/CSB/01	F-101
137	Chiller Sonicator	S clean	SM200 CH US	CDTL/CSB/02	T-204
138	Tensile compression Device	Mecmesin	AFG-100N	CDTL//TCD/03	T-303
139	Tyvec pouch Burst strength	NA	NA	CDTL/TPBST/01	S- 203
140	TOC Analyzer	shimazu	TOC-LCPH	CDTL/TOC/1	G-02
141	Tablet Friability Tester	lab india		CDTL/TFT/01	G-013
142	Tablet Hardnes Tester	lab india		CDTL/THT/01	G-013
143	UV/VIS SPECTROPHOTOMETER	LabIndia	UV3200	CDTL/ UV/01	G-03
144	UV/VIS SPECTROPHOTOMETER	LabIndia	UV3000+	CDTL/ UV/02	S-203
145	UV/VIS SPECTROPHOTOMETER	Perkin Elmer	UVlamda45	CDTL/ UV/03	F-117
146	UV-Vis Spectorophotometer	Agilent	Carry-60	CDTL/UV/04	T-302
147	Vaccume oven	HMG	NA	CDTL/VO/04	S-204
148	Vaccume oven	Fourtech	NA	CDTL/VO/05	S-205
149	Vaccume oven	META LAB	NA	CDTL/VO/07	G-003
150	Vaccume oven	JEIOTECH	ov-12	CDTL/VO/08	F-105
151	Vaccume oven	JEIOTECH	ov-12	CDTL/VO/09	T-304
152	Water Bath	MIC	NA	CDTL/WB/06	S-225 micro
153	Water Bath 5 lit	borosil	WL S005	CDTL/WB/18	S-205
154	Water Bath	Lab quest	WBC 012	CDTL/WB/17	G-10
155	Chiller Water Bath	Scieticfic sales	NA	CDTL/CWB/01	S-225 micro
156	shaker Water Bath		NA	CDTL/CWB/03	T-305
157	Water purifier	Pure Lab Ultra	ELGA DV-35	CDTL/WP/01	T-304
158	Water purifier	Millipore	Elix	CDTL/WP/05	F-105
159	Water purifier	Thermo	Genpure	CDTL/WP/04	G-10
160	Water purifier	Thermo	Genpure	CDTL/WP/03	S-201
161	Zone Reader			CDTL/A/ZR/01	Micro
162	colony counter			CDTL//DCC/01	Micro
163	Liquid Density meter	Anton Paar	DMA-501	CDTL//DM/01	F-117
164	Rotary Evaporator	Science enterprises	Sciencetech SE 119	CDTL/RE/01	T-304
165	VISCOMETER with Helipath	Anton Paar	visco-QC-300	CDTL/VM/02	T-305

166	Column wahing pump	pereless	WP-2000	CDTL/CWP/01	F-125
167	Vortex	Remi	NA	CDTL/VOR/03	Micro
168	Vortex	SPINIX	NA	CDTL/VOR/04	micro
169	Vortex	Remi	CM101	CDTL/VOR/08	F104
170	Vortex	Labquest	VMSP	CDTL/VOR/10	F101
171	Vortex	Labquest	VMSP	CDTL/VOR/11	T-301
172	Vortex	Labquest	VMSP	CDTL/VOR/12	S-205
173	Vortex	Labquest	VMSP	CDTL/VOR/06	S-204
174	Spot Extractor	Prama		CDTL/SPEX/01	F-104
175	Spot Extractor	Prama		CDTL/SPEX/02	F-104
176	Spot Extractor	Prama		CDTL/SPEX/03	F-101
177	Spot Extractor	Prama		CDTL/SPEX/04	F-101
178	Spot Extractor	Prama		CDTL/SPEX/05	S-202
179	Spot Extractor	Prama		CDTL/SPEX/06	S-203
180	Spot Extractor	Prama		CDTL/SPEX/07	S-204
181	Nail dryer	Labtime		CDTL/COS-ND/01	T-303
182	Pay off measurment	Labtime		CDTL/COS-PO/01	T-303
183	Atomatic scrach heater	Labtime		CDTL/COS-HT/01	T-303
184	Peel tester	Labtime		CDTL/COS-PT/01	T-303
185	Break load tester	Labtime		CDTL/COS-BLT/01	T-303
186	Digital Force Gauge	Chattilon	DFA-010	CDTL/DFG/01	T-303
187	Digital Screw Gauge	Mitutoyo	293-561-30	CDTL/DSG/01	T-303
188	Foam Height Apparatus	NA	NA	CDTL//FHA/01	S-213
189	Needle burner & Syringe distroyer	Hi Media	NA	CDTL/NSD/01	
190	Traceble Digital Calliper	Fisher Scientific	NA	CDTL//TDC/01	T-303
191	Traceble Digital Timer	Fisher Scientific	NA	CDTL//TDT/01	T-303
192	Ultra violet Nail paint Dryer			CDTL/COS-ND/01	T-303



