

Statement Indicating Status on Grant of Approval by CLAA for Manufacture of Notified Medical Devices (January 2017 to December 2017)

S. No.	Name of the Manufacturer	Product Licensed by CLAA	Manufacturing Licence No. & Validity
1.	M/s. Disposafe Health and life care Ltd., Plot No. 1 & 2, Sector-59, Phase II, Ballabgarh-121004, Faridabad , Haryana	1. Intravenous Cannula- Polypropylene, Ploythylene, Polyacetal, Stainless Steel (AISI304)-Sterile-E (Total 2 Items)	Mfg.Lic.No.-714-B(H) valid upto – (26-05-2018) (02-01-2017) Additional Products
2.	M/s. RNS Surgical Pvt. Ltd., Plot No. 199, Sector-16, HSIIDC, Bahadurgarh, Haryana	1. Bone Screws- 316L/ 316LVM/Titanium-Non Sterile-5yrs-E & D 2. Bone Plates - 316L/ 316LVM/Titanium - Non Sterile-5yrs-E & D 3. Intramedullary Nails - 316L/ 316LVM/Titanium - Non Sterile-5yrs-E & D 4. Wires & Pins -316L/ 316LVM/Titanium - Non Sterile-5yrs-E & D 5. Prosthesis -316L/ 316LVM/Titanium - Non Sterile-5yrs-E & D	Mfg.Lic.No.-Nil valid upto -Nil (05-01-2017) Fresh License
3.	M/s. Techno Shine Vie Pvt. Ltd., Plot No. 51-52 C, Sai Industrial Estate, Amroli – Kosad Road, Tal: Dist. Surat-394107, Gujarat	1. Bone Screws (locking & Non Locking)- S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D 2. Bone Plates(locking & Non Locking)- S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D 3. Nails - S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D 4. Pins and Wires - S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D 5. Joint Replacement Prosthesis - S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D 6. Spinal Implants - S.S 216, S.S 316LVM, Titanium/ Titanium Alloy – Non Sterile – E & D	Mfg.Lic.No.-G/28/1583 valid upto – (15-09-2021) (05-01-2017) Fresh License
4.	M/s. Sterimed Medical Devices (P) Ltd., 38 KM. Stone Delhi – Rohtak Highway Bahadurgarh, Haryana	1. Umbilical Catheters – Medical grade PVC- sterile – 5yrs – E & D 2. Rectal Catheters - Medical grade PVC- sterile – 5yrs – E & D 3. Foley Ballon catheter 2 way – Latex, PP, Medical Grade PVC - sterile – 5yrs – E & D 4. Foley Ballon catheter3 way – Latex, PP, Medical grade PVC- sterile – 5yrs – E & D 5. Thoracic Trocar catheters - Medical grade PVC, ABS, PP, SS, ALU Needle - sterile – 5yrs – E & D 6. Thoracic Drainage Catheters - Medical grade PVC ABS- sterile – 5yrs – E & D ABS	Mfg.Lic.No.-661-B(H) valid upto – (22-06-2018) (10-01-2017) Additional Products

		<p>7. Suction catheters - Medical grade PVC- sterile – 5yrs – E & D</p> <p>8. Nelaton catheters - Medical grade PVC- sterile – 5yrs – E & D</p> <p>9. Male External catheters – Natural rubber latex – sterile – 5yrs-E & D</p> <p>10. Oxegen catheters - Medical grade PVC- sterile – 5yrs – E & D</p> <p>11. Redon drainage catheters - Medical grade PVC ABS- sterile – 5yrs – E & D</p> <p>12. Urethral catheters - Medical grade PVC- sterile – 5yrs – E & D</p> <p>13. Ryles Tubes (catheters) - Medical grade PVC- sterile – 5yrs – E & D</p> <p>14. Levin/ somach tubes (catheters) - Medical grade PVC- sterile – 5yrs – E & D</p> <p>15. Infant feeding tubes (catheters)- Medical grade PVC- sterile – 5yrs – E & D</p> <p>16. Salem Sump Tubes (Catheters) - Medical grade PVC, LDPE - sterile – 5yrs – E & D</p> <p>17. Kehr’s “T” Tubes (catheters) - Medical grade PVC- sterile – 5yrs – E & D</p> <p>18. Sterivac (Closed wound suction unit) (catheters)- Medical grade PVC, LDPE, ABS, SS, Needle- sterile – 5yrs – E & D</p> <p>19. Endotracheal tubes - Medical grade PVC, pp- sterile – 5yrs – E & D</p> <p>20. Tracheostomy tube (catheters) - Medical grade PVC, LDPE- sterile – 5yrs – E & D</p> <p>21. Yankaur suction set catheter - Medical grade PVC, ABS- sterile – 5yrs – E & D</p> <p>22. Ureteric catheters- poly urethane – sterile – 5yrs- E & D</p>	
5.	M/s. Multimedics, Plot No. 28, HPSIDC, Phase III, Baddi, Dist. Solan , 173205, Himachal Pradesh (HP)	Biodegradable Polymer based Everolimus Eluting Coronary Stent System – Co-Cr 605 stent PTCA catheter everolimus drug – sterile – 2ys- E & D (total 3 items)	<p>Mfg. Lic.No.- 10-01-2017</p> <p>Valid upto 21-01-2021</p> <p>(10-01-2017)</p> <p>Additional products</p>
6.	M/s. Medi-tech Devices Pvt. Ltd., 24, Gujarat Pharma Techno Park, Opp. Zydus Pharma SEZ, Matoda, Changodar, Ahmedabad – 382213, Gujarat	<p>1. Central Venus Catheter Long Term – silicon – sterile – 5yrs- E & D</p> <p>2. Central Venus Catheter with Port – Silicon _ Sterile -5yrs – E & D</p> <p>3. PICC Central Venous Catheter – PU/ Silicon – 5yrs – E & D</p> <p>4. Peritoneal Dialysis catheter - PU/ Silicon – 5yrs – E & D</p> <p>5. Ureteral Stent (Double J Stent) – Polyurethane – Sterile – 5 yrs – E & D</p> <p>6. Foley’s Catheter – Latex/ silicon –sterile -5yrs- E & D</p> <p>7. Endotracheal Tube – PVC – Sterile – 5yrs – E & D</p> <p>8. Trachostomy Tube –PVC – Sterile -5yrs- E & D</p>	<p>Mfg. Lic.No.- G/28/1468</p> <p>Valid upto – 03-02-2019</p> <p>(10-01-2017)</p> <p>Additional Products</p>
7.	M/s. Aurolab, No. 1,	Hydrophobic foldable acrylic – Hydrophobic foldable acrylic – sterile –	Mfg.Lic.No.-

	Sivagangai Main Road, Veerapanjan, Madurai – 625020	5yrs – E & D	TN 00002501 valid upto 10-08-2019 (10-01-2017) Additional Products
8.	M/s. Vascular Concepts Ltd., Doddaballapur Industrial Area, Doddaballapur, Bangalore – 561203, Karnataka	<ol style="list-style-type: none"> 1. Prolink LP- Coronary Stent System – Stainless Steel Balloon Expandable – Medical Grade 316L SS, Thermoplastic Resin, Platinum Alloy (Pt/Ir 90/10) Hypotube – sterile (Ethylene Oxide Gas) – 2yrs - D 2. Propass-Platinum Activated System - Coronary Stent System – Stainless Steel Balloon Expandable - Medical Grade 316L stainless steel coated with medical grade platinum, Thermoplastic Resin, Platinum Alloy (Pt/Ir 90/10) Hypotube – sterile (Ethylene Oxide Gas) – 2yrs - D 3. Prolink SV - Coronary Stent System – Stainless Steel Balloon Expandable – Medical Grade 316L SS, Thermoplastic Resin, Platinum Alloy (Pt/Ir 90/10) Hypotube – sterile (Ethylene Oxide Gas) – 2yrs - D 4. Speed+- PTCA Catheters – Thermoplastic Resin, Platinum alloy (Pt/Ir 90/10), Hypotube, Polymid compound, Medical Grade AISI 304, multilayer Schlauch (PEBAX/PE), Polyimide - sterile (Ethylene Oxide Gas) – 2yrs – D & E 5. Speed - PTCA Catheters – Thermoplastic Resin, Platinum alloy (Pt/Ir 90/10), Hypotube, Polymid compound, Medical Grade AISI 304, Guide wire tubing, multilayer Schlauch (PEBAX/PE), Polyimide - sterile (Ethylene Oxide Gas) – 2yrs – D 6. Pronova – Sirolimus Eluting Stent System – Cobalt Chromium – L605 cobalt chromium alloy drug, sirolimus polymer, polysulfone catheter body, thermoplastic resin marker bands, platinum alloy (Pt/Ir, 90/10) Hypotube, AISI 304, ballon, PA11 Compound – sterile – 2yrs- D & E 7. Pronova SS- Sirolimus Eluting Stent System – stainless Steel – Medical grade 316SS drug, sirolimus polymer, polysulfone catheter body, thermoplastic resin marker bands, platinum alloy (Pt/Ir 90/10) hypotube, stainless steel balloon, thermoplastic resin – sterile – 2yrs – D 8. Pronova XR – Sirolimus Eluting Stent System –Cobalt Chromium Degradable Matrix- L605 Cobalt Chromium alloy drug, sirolimus polymer, 85/15 DL-Lactide/ Glycolide catheter body, thermoplastic resin marker bands, platinum alloy (Pt/Ir 90/10) Hypotube, stainless steel balloon – 2yrs- sterile – D & E 9. Prostar Plus – balloon expandable peripheral stent system , renal & iliac- AISI 316L stent made of stainless steel marker bands, platinum alloy (Pt/IR 90/10) Hypotube, stainless steel balloon, thermoplastic resin 	Mfg. Lic. No. – KTK/28/359/2006 Valid upto – 08-11-2021 (17-01-2017) (Renewal)

9.	M/s Meril Life Sciences Pvt Ltd., Survey No.135/2a and 135/2/B,174/2(M1-M2, Meril park, Fourth floor), Near G.M Bilakhia Stadium, Mukatanand Marg, Chala vapi, valsad-396191	1. Rapid Exchange Percutaneous Transluminal Coronary Angioplasty balloon Dilatation catheter (Hydrophilic Coated)- Non sterile – 4yrs- E (Total 2 items)	Mfg Lic no.G/28/1273 Valid: 20-06-2013 to 19-06-2018 (19-01-2017) (Additional Products)
10.	M/s Meril Life Sciences Pvt Ltd., Survey No.135/2a and 135/2/B,174/2(M1-M2, Meril park, Fourth floor), Near G.M Bilakhia Stadium, Mukatanand Marg, Chala vapi, valsad-396191	1. Atrial Septal Defect Occluder – Nitinol with polyethylene terephthalate (PET) fabric –sterile – 2yrs – E & D 2. Ventricular septal defect (VSD) Occluder – Nitinol with polyethylene terephthalate (PET) fabric –sterile – 2yrs – E & D 3. Muscular Ventricular septal defect (VSD) Occluder - Nitinol with polyethylene terephthalate (PET) fabric –sterile – 2yrs – E & D 4. Patent Ductus Arteriosus (PDA) Occluder – Nitinol with polyethylene terephthalate (PET) fabric –sterile – 2yrs – E & D 5. ASD Occluder Delivery System – Sterile – 2yrs- E & D 6. VSD Occluder Delivery system – sterile – 2yrs – E &D 7. Muscular VSD Occluder Delivery System – 2 yrs – sterile – E & D 8. PDA Occluder Delivery System – Sterile – 2ys – E & D	Mfg Lic no.G/28/1273 Validity period 20-06-2013 to 19-06-2018 (19-01-2017) (Additional Products)
11.	M/s Meril Life Sciences Pvt Ltd., Survey No.135/2a and 135/2/B,174/2(M1-M2, Meril park, Fourth floor), Near G.M Bilakhia Stadium, Mukatanand Marg, Chala vapi, valsad-396191, Gujarat	1. Self Expanding Nitinol Peripheral Stent System – Nitinol – Sterile – 3yrs – E & D	Mfg Lic no.G/28/1273 Validity period 20-06-2013 to 19-06-2018 (19-01-2017) (Additional Products)
12.	M/s. GPC Medical Ltd., Plot No. A82. Sector A4, UPSIDC, Tronica City, Loni Ghaziabad-UP	1. Cortex Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 2. Cancellous Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 3. Large Cannulated Cancellous Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 4. Small Cannulated Cancellous Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 5. Cannulated Cancellous Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E	Mfg. Lic. No.- ___/SC/P of 2016 Valid : Nil (24-01-2017) Fresh License

		<ol style="list-style-type: none"> 6. Malleolar Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 7. Self-Taping Cortex Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 8. Shaft Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 9. Emergency Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 10. Herbert Cannulated Bone screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 11. Washer for screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 12. Cortical Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 13. D.H.S. /Lag Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 14. D.H.S. Compression Screw/ Blocket Screw- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 15. Paediatric D.H.S. Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 16. Locking Head Screws - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 17. Locking Compression Screw Self Tapping - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 18. Locking Head Cancellous Screw Fully Threaded - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 19. Self Tapping Locking Screw - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 20. Bone Plates – stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 21. Small Fragment Plates- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 22. Large Fragment Plates - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 23. Sherman Bone Plates - stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 24. Locking Plate System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 25. Intramedullary Interlocking Nails - System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E 	
--	--	--	--

		<p>26. Locking Bolts - System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E</p> <p>27. Pins, Wires, Staples - System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E</p> <p>28. Arthroplasty (femoral Heads)- System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E</p> <p>29. Spinal Implants - System- stainless steel & Titanium- sterile & Non sterile – 5 yrs –E</p>	
13.	M/s. Aonelife Medicare Pvt. Ltd., Plot No. D-1/6A, D-1/6B, Sector 13, GIDA, Gorakhpur 273209 UP	<p>1. I.V. Cannula (with injection port & with wing) – FEP/PTFE/PUR, PP, HDPE, LDPE, Poly Acetal, Silicoon Rubber, SS Needle – sterile- 5yrs- E & D</p> <p>2. I.V. Cannula (without injection port & with wing) – FEP/PTFE/PUR, PP, HDPE, LDPE, Poly Acetal, SS Needle – sterile- 5yrs- E & D</p> <p>3. I.V. Cannula (without injection port & without wing) – FEP/PTFE/PUR, PP, HDPE, LDPE, Poly Acetal, Silicoon Rubber, SS Needle – sterile- 5yrs- E & D</p> <p>4. I.V. Cannula (with wing & without injection port Neo) – FEP/PTFE/PUR, PP, HDPE, LDPE, Poly Acetal, SS Needle – sterile- 5yrs- E & D</p>	<p>Mfg. Lic. No.- ___/SC/P of 2016 Valid : Nil (03-02-2017) Fresh License</p>
14.	M/s. Neo Surgical, Plot No. G-306, GIDC Almighty Gate, metoda, Pin-360021	<p>1. DHS Locking Plate-SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>2. DHS Plate D.C. Hole- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>3. DHS Plate Round Hole- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>4. DHS Plate for Femoral Osteotomies - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>5. Dynamic Trochanter Stabilizer Plate used with DHS- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>6. Angle Blade Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>7. DCS Locking Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>8. DCS Plate D.C. Hole- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>9. DCS Plate round Hole - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>10. Proximal Remoral Plate-Locking (Left & Right)- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>11. Distal Femur Locking Plate Left & Right- SS316L, SS316LVM, Titanium-</p>	<p>Mfg. Lic. No. G/28/1585 Valid upto 21-09-2021 (13-02-2017) Additional products</p>

		<p>non sterile-5 yrs- E & D</p> <p>12. Condylar Buttress Plate Left & Right - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>13. Distal Lateral Femoral Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>14. Locking Clavical Plate Left & Right- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>15. Anterior Cervicle Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>16. Posterior Cervicle Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>17. Anterior ThoroColumber Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>18. Occifix Locking Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>19. Cranio Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>20. K wire- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>21. Spacer- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>22. Calceneal Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>23. Dynamic Compression Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>24. Limited Contact Dynamic Compression Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>25. Locking Compression Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>26. Locking Head Screw- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>27. Locking Bolt- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>28. Cortical Screw- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>29. Cortical locking Screw- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>30. Cannulated Cancellous Screw (Short Threaded)- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>31. Cannulated Cancellous Screw (Fully Threaded)- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>32. Cannulated Cancellous Screw (16th Threaded)- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p> <p>33. Cannulated Cancellous Screw (32th Threaded)- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D</p>	
--	--	--	--

34. Locking Medial Distal Tibial Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
35. Lateral Tibial Plate Left & Right- - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
36. "L" Simple Plate Left & Right- - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
37. "T" Buttress Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
38. "L" Buttress Locking Plate & left & Right- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
39. "T" Buttress Locking plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
40. Locking Metaphyseal Plate - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
41. Locking Overleaf Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
42. Locking Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
43. Locking Medial Proximal Tibial Plate Left & Right - SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
44. Low Band Distal Medial Tibia Locking Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
45. Anterper;atera; Distal Tibia Locking Plate Left & Right- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
46. Locking Proximal Tibial Plate Left & Right- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
47. Febula Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
48. Tomo Fix Plate- SS316L, SS316LVM, Titanium-non sterile-5 yrs- E & D
49. Spine Rod- Titanium – Non Sterile – 5 yrs- E & D
50. Staple - Titanium – Non Sterile – 5 yrs- E & D
51. Hook - Titanium – Non Sterile – 5 yrs- E & D
52. Connector- Titanium – Non Sterile – 5 yrs- E & D
53. Washer Small - Titanium – Non Sterile – 5 yrs- E & D
54. Washer Big- Titanium – Non Sterile – 5 yrs- E & D
55. Occipital Connector- Titanium – Non Sterile – 5 yrs- E & D
56. Pre Bent Rod- Titanium – Non Sterile – 5 yrs- E & D
57. Endo Button - Titanium – Non Sterile – 5 yrs- E & D
58. Tibia Nail- SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D
59. Femoral Nail SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D
60. Proximal Femoral Nail - Short SS316L, SS316LVM Titanium – Non

		<p>Sterile – 5 yrs- E & D</p> <p>61. Distal Femoral Nail - SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>62. Radius Ulna Nail SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>63. End Cap for Tibia Nail SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>64. End cap for femoral nail SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>65. End cap for proximal femoral nail-short SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>66. End cap for Distal femoral nail SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>67. ST Pin SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>68. Guide Pin SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p> <p>69. Guide Wire SS316L, SS316LVM Titanium – Non Sterile – 5 yrs- E & D</p>	
15.	M/s. Aeonmed Medical Industries, B-XXIX-536/12-C/4, Singla Cycle Road, mandir Street, Opposite Dhandari Kalan Railway Station, Ludhiana, Punjab	<p>1. Screws- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>2. Bone Plates- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>3. Bone nails - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>4. Wires, Pins and Staples - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>5. Locking Plates - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>6. Locking Plates & Locking Head Screws Small FRG- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>7. Locking Plates & Locking Head Screws large FRG- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>8. Locking Anatomical Bone Plates - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>9. Spinal System- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>10. Anterior Cervical Locking Plate - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>11. Maxillofacial Implants - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>12. Screws (CPI TYPE) – 1.5mm modular Hand Implants -- SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>13. Hip Prosthesis - - SS316LVM, Titanium- Aeonmed, Staymaxx</p> <p>14. Thomson Prosthesis Excek, Standard Stem (Sterile)- Aeonmed, Staymaxx</p> <p>15. Thomson Prosthesis Excek, Standard Stem (Non Sterile)- Aeonmed, Staymaxx</p> <p>16. Thomson Prosthesis Excek, Narrow Stem (Sterile)- Aeonmed, Staymaxx</p>	<p>Mfg. Lic. No.-1992-B Valid upto 26-04-2021</p> <p>(22-03-2017) Fresh License</p>

		<ol style="list-style-type: none"> 17. Thomson Prosthesis Excek, Narrow Stem (Non Sterile)- Aeonmed, Staymaxx 18. Austin Moore – Step ball (Non Sterile) 19. Austin Moore – Step ball (Sterile) 20. Austin Moore Hip Prosthesis– Standard (Sterile), Stainless Steel – Aeonmed , Staymaxx 21. Austin Moore Hip Prosthesis– Standard (Non Sterile), Stainless Steel – Aeonmed , Staymaxx 22. Austin Moore Hip Prosthesis– Narrow Stem(Sterile), Stainless Steel – Aeonmed , Staymaxx 23. Austin Moore Hip Prosthesis– Narrow Stem(Non Sterile), Stainless Steel –Aeonmed , Staymaxx 24. Austin Moore Hip Prosthesis– broad Stem(Sterile), Stainless Steel – Aeonmed , Staymaxx 25. Austin Moore Hip Prosthesis– Broad Stem(Non Sterile), Stainless Steel –Aeonmed , Staymaxx 26. Austin Moore Hip Prosthesis– Long Stem(Sterile), Stainless Steel – Aeonmed , Staymaxx 27. Austin Moore Hip Prosthesis– Long Stem(Non Sterile), Stainless Steel – Aeonmed , Staymaxx 28. Bipolar Hip Prosthesis – (Sterile), Stainless – Aeonmed, Staymaxx 29. Bipolar Hip Prosthesis – (Non Sterile), Stainless – Aeonmed, Staymaxx 30. Bipolar hip Prosthesis, Fenestrated –(Sterile), Stainless Steel – Aeonmed Staymaxx 31. Bipolar hip Prosthesis, Fenestrated –(Non- Sterile), Stainless Steel – Aeonmed Staymaxx 32. THR – C STEM – Aeonmed, staymaxx 33. THR – C CUP – Aeonmed, Staymaxx 34. C – Bipolar Cup – Aeonmed, Staymaxx 35. THR – M Stem – Aeonmed, Staymaxx 36. Thr – M head, Stainless Steel – Aeonmed, Staymaxx 37. THR-M Cup (Acetabular Cup) – Aeonmed, Staymaxx 38. THR – M Bipolar Cup – Aeonmed, Staymaxx 	
16.	M/s. Romsons Scientific & Surgical Industries Pvt. Ltd., 63-64, industrial Estate, nunhai, Agra- 282006 (U.P.)	<ol style="list-style-type: none"> 1. I.V. Cannula with injection port, non-return valve & wing – FEP, S.S, Sili cone, P.P, POM & M.B – sterile -5yrs- E & D 2. I.V. Cannula with wing but without injection port- FEP, S.S, M.B. Silicone, P.P, POM – sterile -5yrs- E & D 3. I.V. Cannula with integrated 3 way stop cock & wing – FEP, S.S, M.B., Silicone, P.P, PC, POM –Sterile – 5ys – E & D 	<p>Mfg. Lic. No.- Drug/3164/ Valid upto 17-03-2018</p> <p>28-02-2017 (Additional Products)</p>

4. I.V. Cannula without injection port & without wings – FEP, S.s, P.P, POM, MB – 5ys- E & D
5. Paediatric I.V. Cannula with wings but without injection port – FEP, S.S, P.P, POM, MB- sterile – 5ys – E & D
6. Infant I.V. Cannula with wings but without injection port – FEP, S.S, P.P, POM, MB – sterile – 5ys- E & D
7. Paediatric I.V. Cannula with wings but without injection port – FEP, S.S, P.P, POM, MB – sterile – 5ys- E & D
8. I.V. Cannula with injection port & with wings – FEP, S.S, MB, Silicone, P.P, POM - sterile – 5ys-E & D
9. I.V. Cannula with Injection Port – FEP, S.S, MB, Silicone, P.P - sterile – 5ys-E & D
10. I.V. Cannula with wing- Paediatric use – FEP, P.P, POM, MB- sterile – 5YS- e & d
11. Safety cannula with Luer Lock – FEP, S.S, Silicone, P.P, POM & MB – 5ys –sterile – E & D
12. Safety cannula with injection port, wing & Luer Lock Cap – FEP, S.S, P.P, POM, MB – Sterile -5ys- E & D
13. I.V. Cannula – TEFLONSS, MB, PP, AC ETAL, COP OLYMER, SILICONE –sterile- 5ys- E& D
14. Pediatric I.V. Cannula – Teflon, SS, MB, PP, Acetal Copolym ER – sterile- 5ys- E & D
15. Infant I.V. Cannula - Teflon, SS, MB, PP, Acetal Copolym ER – sterile- 5ys- E & D
16. Scalp Vein Set – PVC, MB, PE, SS- sterile – 5ys- E & D
17. Penile Sheath / External Catheter – latex, Sterile – 3ys- E & D
18. Nelaton catheter for male- PVC, MB – Sterile – 5ys- E & D
19. Nelaton catheter for female- PVC, MB – Sterile – 5ys- E & D
20. Nelaton catheter Hydrophyllic Coated- PVC, MB – Sterile – 5ys- E & D
21. Foley Balloon Catheter, 100% Silicone – Silicone – sterile- 5ys – E & D
22. Umbilical catheter (Arterial & venous catheter) – PVC sterile 5ys- E & D
23. Foley Balloon Catheter, Natural rubber latex silicone elastomer bonded – PP, ABS, POM, MB, SS – sterile – 3ys – E & D
24. Urethral Catheter- PP, MB, sterile-5ys- E & D
25. Ureteric catheter (open Tip) – PU, PE, Polycarbonate – sterile- 5yrs- E & D
26. Ureteric catheter (Close Tip) – PU, PE, Polycarbonate – sterile- 5yrs- E & D
27. Kit-Double”j” stent, Pusher, Guide – wire & cath clamp (Open & Close

- end) – SS – sterile – 5ys-E & D
28. Double “J” Stent (Open & Close End) – PP, PE PC, POM, SS, MB- 5ys-sterile- E & D
 29. Suction Catheter – Plain- PVC – sterile – 5ys – E & D
 30. Suction catheter- TC (Thumb Control Suction Catheter) – PVC – sterile – 5ys- E & D
 31. Suction Catheter FT (Finger Tip Control Suction Catheter) – PVC, MB-sterile- 5ys- E & D
 32. Suction Catheter Suckath (Finger Tip Control, Graduated, Suction Catheter) – PVC, MB – sterile -5ys- E & D-
 33. Suction catheter integrated ring type flower tip with round eyes with thumb control – pvc, mb, 5ys- sterile- E & D
 34. Suction catheter – cee tee cath- PVC, MB –sterile – 5ys- E & D
 35. Ryle’s Tube Naso Gastric Tube – PVC, MB, SS, sterile- 5ys- E & D
 36. Romolene – Ryle’s tube with Funnel & Luer Connector – PVC, MB-sterile – 5ys- E & D
 37. Levin’s Tube naso Gastric Tube - PVC, MB- sterile – 5ys- E & D
 38. Salem Sump Tube – PVC, MB, PP – Sterile -5ys- E & D
 39. Infant Feeding Tube – PVC, MB – sterile -5ys- E & D
 40. 100% Silicone Enteral Feeding Tube for Long Term use – silicone – 5ys- E & D
 41. Chest drainage catheter – straight –PVC, MB, PP –sterile – 5ys- E & D
 42. Extra-Soft Thoracic Drainage catheter – straight – PVC, MB- sterile – E & D
 43. Extra-Soft Thoracic Drainage catheter – Right Angle – PVC, MB- sterile – E & D
 44. Thoracic Drainage catheter with Trocar – PVC, MB- sterile-5ys – E & D
 45. Epidural Catheter – PEBAX, P- sterile – 5ys- E & D
 46. Central Venous Catheter – Kit with Soft “J” Tip Guidewire with Dispenser – PU, ASB, PP, PC, SS, PE, MB- Sterile – 5ys- E & D
 47. Kehr’s “T” Tube- PVC – sterile – 5ys- E & D
 48. Peritoneal Dialysis catheter set – SS, ABS, PVC, PP, MB, PVC Film, Isoprene –sterile-5ys- E & D
 49. Peritoneal Dialysis Transfusion set – PE, ABD, PP, PVC, Isoprene Acrylic, Polyamide-sterile-5ys – E & D
 50. Supra Public balloon Catheter & Trocar – PE, SS PVC, MB, ABS – Sterile – 5ys- E & D
 51. Under water seal drainage system- K Resin, PVC, PE, PP, SS, MB, Adhesive paper – sterile – 5ys- E & D

- | | | | |
|--|--|--|--|
| | | <p>52. Closed wound suction unit-PP, PE, PVC, ABS, SS, MB –sterile- 5ys- E & D</p> <p>53. Closed wound suction unit-mini -PP, PE, PVC, ABS, SS, MB –sterile- 5ys- E & D</p> <p>54. Closed wound drainage system – PP, PE, PVC, ABS, SS, MB-sterile – 5ys- E & D</p> <p>55. Abdominal Drainage kit – PVC, MB-sterile – 5ys- E & D</p> <p>56. Oro-pharyngeal Airways –LDPE, PE, ABS, MB-sterile – 5ys- E & D</p> <p>57. Nasopharyngeal Airways-PVC-Sterile- 5ys- E & D</p> <p>58. Enteral Feeding Bag with Ice Poucj & Transfer set with drip chamber & flow controller & tapered step connector – PVC, PP, ABS, MB – Sterile - 5ys- E & D</p> <p>59. Sump Suction handle with suction tube, enables modification as per suitability; accessory of catheter – ABS, PVC, MB- sterile – 5ys- E & D</p> <p>60. Epidural catheter- soft tip – PEBAX, P- sterile – 5ys- E & D</p> <p>61. Kit-Epidural catheter-soft Tip, Epidural Needle & Accessories –PEBAX, ABS, PC, PE, SS, MB – Sterile -5ys- E & D</p> <p>62. Sterile tube for suction- PE, PVC, MB-sterile-5ys-E & D</p> <p>63. Mist Blower for improved visibility at anastomotic of other surgical sites – PE, ABS,MB,PC, SS, Polymide –Sterile – 5ys- E & D</p> <p>64. Suction handle with Tip,Handle provided wih filter,accessory of catheter- PVC,ABS – Sterile – 5 – e & d</p> <p>65. Suction tube with connectors at both end to connect to Suction Handle/Catheter;Accessory of Catheter- ABS,PVC,MB – Sterile – 5 - E & D</p> <p>66. Suction Handle (standard) with suction tube with flexible connectors at ends;Accessory of Catheters – PP,PVC,MB – Sterile – 5 - E & D</p> <p>67. Suction tube with Flexible connectors at both ends & provided with spare step connectors for connection to catheter;Accessory to Catheter – ABS,PVC,MB - -Sterile – 5 - E & D</p> <p>68. Suction Handle (Standard & Crown Tip) with Suction Tube;Accessory of Catheter – PVC,ABS,MB – Sterile – 5 - E & D</p> <p>69. Closed Suction system with suction Catheter – ABS, PVC, MB, PE - Sterile – 5 - E & D</p> <p>70. Intermittent Urethral Catheterization Procedure Kit – ABS,PVC,MB - Sterile – 5 - E & D</p> <p>71. Kit-Urine Collecting Bag with Foley Balloon Catheter – ABS,PVC,MB,PP,NRL - Sterile – 5 - E & D</p> <p>72. Guidewire-Straight and “J” Tip – S.S.P.E. - Sterile – 5 - E & D</p> | |
|--|--|--|--|

73. 100 % Silicone Supra Pubic balloon Catheter with trocar – Silicone,S.S.P.P,ABS,PVC,MB - Sterile – 5 - E & D
74. Foley balloon Catheter 3-way, Natural Rubber Latex,Silicone Elastomer Bonded – P.P,ABS,POM,MB,SS - Sterile – 3 - E & D
75. 2-way & 3-way Foley Balloon Catheter,Silicone coated – P.P,ABS,POM,MB,SS - Sterile – 3 - E & D
76. Suction Catheter with flower tip with Thumb Control – PVC,M.B - Sterile – 5 - E & D
77. Suction Catheter- Plain with flower tip – PVC,MB - Sterile – 5 - E & D
78. Catheter with Mucus Extractor with/without filter – P.P,PVC, PE,MB - Sterile – 5 - E & D
79. Catheter with Mucus Extractor for attachment to Bronchoscope – PE,PS,P.P,PVC,MB - Sterile – 5 - E & D
80. Catheter –left ventricular Sump vent Catheter – PVC, M.B - Sterile – 5 - E & D
81. Catheters-Ante grade coronary perfusion Cannula – Silicone,ABS,MB - Sterile – 5 - E & D
82. Kit-Catheter Naso-Gastric,Syringe 20 ml & Collection Bag (2000 ml) – PVC,M.B. - Sterile – 5 - E & D
83. Polyurethane Naso-gastric,FeedingTube with Radio-Opaque Line and Optional Guidewire – Silicone - Sterile – 5 - E & D
84. Cardioplegia Cannula – Silicone, ABS,MB,SS - Sterile – 5 - E & D
85. Oral/Nasal Tracheal Tube without cuff & with Murphy eye & Radio-opaque Line – PVC,P.P,M.B - Sterile – 5 - E & D
86. Oral/Nasal Tracheal Tube with cuff,Murphy eye & Radio-opaque Line,Pilot Balloon – PVC, P.P,M.B - Sterile – 5 - E & D
87. Oral/Nasal Tracheal Tube Cuff, Tube with inner wire inside wall for reinforcement,kink resistant,Murphy eye radio-opaque line, Pilot Balloon – PVC, P.P, M.B - Sterile – 5 - E & D
88. Laryngeal mask,airways Catheter – PVC - Sterile – 5 - E & D
89. Laryngeal mask, Airways Catheter-100% Silicone – PVC - Sterile – 5 - E & D
90. Tracheostomy Tube Catheter-Plain, without cuff with Stylet, with radio-opaque line – P.P, PVC,PE,Polyamide M.B - Sterile – 5 - E & D
91. Tracheostomy Tube Catheter-Plain,with pilot balloon with Stylet, with radio opaque line – P.P, PVC, PE, Polyamide M.B - Sterile – 5 - E & D
92. 100 % Silicone Tracheostomy Tube Catheter with Cuff & Pilot Balloon & radio-opaque line – P.P,PVC, PE, Polyamide - Sterile – 5 - E & D
93. Enteral Feeding bag with Ice pouch & transfer Set with-Drip Chamber &

		<p>Flow Controller & Tapered Step Connector – PVC, PP, ABS, MB - Sterile – 5 - E & D</p> <p>94. Polyurethane (PUR) Naso-gastric Tube, Closed Tip, Radio Opaque Line & Tip.Connector compatible with Catheter Mount Syringe, Lucer Slip Syringe & Feeding Bag tapered Connector.Provided with Lateral Eye – PUR,PVC, M.B - Sterile – 5 - E & D</p> <p>95. Kit for Pleural Drainage.Kit consists of 32 sets of Extra Soft Thoracic Drainage Catheter & 3 sets of Drainage bottles with individual hanger & stand - Sterile – 5 - E & D</p> <p>96. External Catheter (Penile Sheath) – Latex - Sterile – 5 - E & D</p> <p>97. Foley Ballooon Catheter-2 Way-Paediatric,Natural Rubber Latex, Silicone Coated – P.P,ABS,POM,MB,SS - Sterile – 3 - E & D</p> <p>98. Foley Balloon Cathetar-2 Way- Adult, Natural Rubber Latex, Silicone Coated – P.P,ABS,POM,MB,SS - Sterile – 3 - E & D</p> <p>99. Foley Balloon Cathetar-3 Way- Adult, Natural Rubber Latex, Silicone Coated – P/P, ABS,POM, MB, SS - Sterile – 3 - E & D</p> <p>100.Suction Cathetar with funnel Connector (Endo Bronchial Suction Cathetar) (Suction Cathetar with Funnel Connector, Open Tip & Lateral Eye) – PVC- Sterile – 5 - E & D</p> <p>101. Suction Cathetar Thumb Control valve with Stopper(Endo Bronchial Suction Cathetar), (Open tip & Lateral Eye) – PVC -Sterile – 5 - E & D</p> <p>102. Guedel Airways- Oropharyngeal Airways – LDPE,PE,ABS,M.B - Sterile – 5 - E & D</p> <p>103. Ryle’s Tube-Naso Gastric Tube for Enternal Feeding and also for Aspiration of gastric secretions,with Radio-opaque line & Radio-opaque Tip& closed Tip – PVC, M.B, S.S - Sterile – 5 - E & D</p> <p>104. Infant Feeding Tube-Naso Gastric Tube for Enteral feeding and Aspiration of gastric secretions with Radio-opaque line and closed tip.Provided with Lateral eyes.Tube marked at 20cm distance from Tip – PVC, M.B - Sterile – 5 - E & D</p> <p>105. Chest Drainage Cathetar-Right Angle Provided with large lateral eyes.Open tipped-Cathetar with radio-opaque line,Spare chest connector provided – PVC,M.B,P.P - Sterile – 5 - E & D</p> <p>106. Chest Drainage Cathetar-Right Angle Provided with large Lateral Eyes.oPEN Tipped-Cathetar with radio –opaque line,spare chest connector provided – PVC, M.B, P.P - Sterile – 5 - E & D</p> <p>107. Closed wound suction unit- for drainage of surgical debris post-operatively under negative pressure.Consists of perforated catheter with Radio-Opaque line, connecting tube with clamp & collection Jar. (Bellow</p>	
--	--	---	--

		<p>shaped) – P.P,PE,PVC,ABS,S.S,M.B - Sterile – 5 - E & D</p> <p>108. Yankaur suction Set-Suction Handle (Standard & Crown Tip) with suction tube – PVC - Sterile – 5 - E & D</p> <p>109. (Silicone, Male Condom Cathetar) (Penile Sheath) – Silicone - Sterile – 3 - E & D</p> <p>110. Percutaneous Tracheostomy Tube (Suitable for Percutaneous Procedure) – PVC, PP, MB, ABS, SS - Sterile – 5 - E & D</p> <p>111. Percutaneous Tracheostomy Kit (kit with Guidewire Dilating Forceps) – PVC, PP, MB, ABS, Polyamide, Teflon - Sterile – 5 - E & D</p> <p>112. Percutaneous Tracheostomy Kit (Kit without dilating forceps) – PVC, PP, MB, ABS, Polymide, Teflon - Sterile – 5 - E & D</p> <p>113. Endobronchial Tube-Left with Stylet, 15mm Connector, Circuit Adaptor with in-built Sampling Port & a pair of suction Catheters- PVC, MB, ABS, SS, PP, Aluminium - Sterile – 5 - E & D</p> <p>114. Endobronchial Tube-Right with Stylet, 15mm connector, Circuit Adaptor with in-built Sampling Port & a pair of Suction Catheters – PVC, MB, ABS, SS, PP, Aluminium - Sterile – 5 - E & D</p> <p>115. Laryngeal Mask Airway- Double Lumen – PVC, SS, MB, ABS, POM - Sterile – 5 - E & D</p> <p>116. Under water Seal drainage system with excess negative pressure control valve – K-Resin, SS, PP, PVC, PE, MB, TPE, POM - Sterile – 5 - E & D</p> <p>117. 100% Silicone thoracic drainage Cathetar-Straight with Radio-Opaque Line – 100 % Silicone Radio Opaque, MB, PP - Sterile – 5 - E & D</p> <p>118. 100% Silicone thoracic drainage Cathetar-Angles with Radio-opaque Line – 100% Silicone Radio Opaque, MB, PP - 100% Silicone thoracic drainage</p> <p>119. Manual Vacuum Aspiration Cannula – PP, MN- Sterile – 5 - E & D</p> <p>120. High Capacity Chest Drainage System – PP, MB, ABS, PVC, TPR, PE - Sterile – 5 - E & D</p> <p>121. High Capacity Chest Drainage System – PP, MB, ABS, PVC, TPR, PE - Sterile – 5 - E & D</p>	
17.	<p>M/s Meril Healthcare Pvt.Ltd., Survey No.135/2/B & 174/2, First Floor, H1-H3, Meril Park, Muktanand Marg, Chala, Vapi-369191</p>	<p>1. Total High Replacement (Sterile) – Ti-6Al-4V ELI – Coating material :- CPTi+HA coating, CPTi coating, HA Coating – Sterilised by Gamma Radiation – 5 – E & D</p> <p>2. Total High Replacement (Non Sterile) - Ti-6Al-4V ELI – Coating material :- CPTi+HA coating, CPTi coating, HA Coating – Sterilised by Gamma Radiation – 5 – E & D</p> <p>3. Additional Components for “Latitud” - Total Hip Replacement</p>	<p>Mfg. Lic. No.- G/28/1465 Valid upto 31.01.2019</p> <p>28-02-2017 (Additional Products)</p>

		System(Sterile) – Modular Femoral Stem, Modular Femoral Neck, Proximally coated Uncemented Stem (Ti+HA Coating), (Ti Coating), (HA Coating),Sealing Plug 4. Additional Components for “Latitud”-Total Hip Replacement System (Non Sterile) – Modular Femoral Stem, Modular Femoral Neck, Proximally coated Uncemented Stem (Ti+HA Coating), (Ti Coating), (HA Coating),Sealing Plug	
18.	M/s Meril Life Sciences Pvt.Ltd., Survey No.135/2a and 135/2/B & 174/2, (M1-M2,Meril Park, Fourth Floor), Near G.M.Bilakhia Stadium,Muktanand Marg,Chala,Vapi,Valsad-369191,Gujarat,India	1.Balloon Sinus Dilatation System – Pebax – Sterile – 2 - E & D 2. Mesire Sphenoid – Polyurathene & Stainless Steel, Stainless Steel – Sterile – 2 - E & D	Mfg. Lic. No.- G/28/1273 Valid upto 19.06.2018 28-02-2017 (Additional Products)
19.	M/s Wellmed International Industries Pvt.Ltd.,A-176 & 177,Sector-63,Noida-201301,Uttar Pradesh	IV Cannula with Wings & 3 Way – FEP/PTFE/PUR,PP,PC,HDPE,LDPE,POLY,ACETAL,SILICON,RUBBER,SS NEEDLE – Sterile – 5 – E & D	Mfg. Lic. No.- 07/SC/P of 2008 Valid upto 28-02-2017 (Additional Products)
20.	M/s Meril Life Sciences Pvt.Ltd., Survey No.135/2a and 135/2/B,174/2 (M1-M2,Meril Park,Fourth floor),Near G.M.Bilakhia Stadium,Muktanand Marg,Chala,Vapi,Valsad-396191,Gujarat,india	Sirolimus Eluting BioResorbable Vascular Scaffold System – PLLA & PDLA Polymers – Sterile – 1 – E & D	Mfg. Lic. No.- G/28/1273 Valid upto 19.06-.2018 07-03-2017 (Additional Products)
21.	M/s Harsoria Healthcare Pvt.Ltd.,110-111,Phase-IV,Udyog Vihar,Gurgaon,Haryana-122015	1. Healfon (Sterile Disposable IV Cannula with Catheter and Injection Valve for single use) – Polypropylene,HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 Silicone tube – Sterile – 5 – E & D 2. Healvan (Sterile Disposable IV Cannula with Catheter & and Without Injection Valve & with small wings for single use only) - Polypropylene,HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS	Mfg. Lic. No.- DML:-637-B (H) Valid upto 18-09-2021 08-03-2017 (Renewal of License)

Needle in AISI304 – Sterile – 5 – E & D

3. Healcathtm (Sterile disposable IV Cannula with catheter and without injection Valve & without wings for single use only)- Polypropylene,HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 – Sterile – 5 – E & D
4. Healvantm & Autosaftm (Sterile Disposable IV Cannula for single use only) - Polypropylene,HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 Silicone Tube – Sterile – 5 – E & D
5. HealcathTM & AutosafTM (Sterile disposable Safety IV Cannula with Catheter for single use only) - Polypropylene,Polycarbonate ,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 SS spring AISI302 – Sterile – 5 – E & D
6. Canusoft (Sterile Disposable IV Cannula for single use only) - Polypropylene,HDPE,ABS,PolyAcetal,Cathetar (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube – Sterile – 5 – E & D
7. Healvan SafloTM (Sterile Disposable IV Cannula with Catheter & Without Injection Valve with small wings with in line one valve for single use only) – Polypropylene ,HDPE,ABS,Polyacetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
8. Inflo (Sterile Disposable IV Cannula with catheter for single use only) – Polypropylene HDPE,ABS,Poly Acetal,Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
9. Sensecure (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE,ABS,Poly Acetal,Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 – E & D
10. Healflon Plus (Sterile Disposable IV Cannula for single use only) – Polyprpylene HDPE,ABS,Poly Acetal,Cathetar (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube – Sterile – 5 – E & D
11. Healthcath Plus (Sterile Disposable IV Cannula for single use only) – Polypropylene HDPE,ABS,Poly Acetal,Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
12. Healvan Plus (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE,ABS,Poly Acetal,Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
13. Nubenco Port (Sterile Disposable IV Cannula for single use only) – Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
14. Nubenco Wings (Sterile Disposable IV Cannula for single use only) –

		<p>Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <ol style="list-style-type: none"> 15. JCM MED (Sterile Disposable IV Cannula with Catheter and with injection valave for single use only) – Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D 16. JCM MED (Sterile Disposable IV Cannula with catheter and without injection valve and with small wings for single use only (Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D 17. JCM MED (Sterile Disposable IV Cannula with Catheter and Without injection Valve and without wings for single use only) - (Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D 18. JCM MED I.V.SAFE (Sterile Disposable safety IV Cannula with Catheter and injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in Slicione Tube - Sterile – 5 - E & D 19. JCM MED I.V Safe (Sterile disposable safety IV Cannula with Catheter and without injection valve and with small wings for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 - E & D 20. JCM MED I.V Safe (Sterile disposable safety IV Cannula with catheter and without injection valve and without wings for single use only) -Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D 21. NEOTOC (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D 22. NEOTEC-1 (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D 23. NEOTEC-1W (Sterile disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D 24. VENOFIX (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D 25. Venofix Plus (Sterile Disposable IV Cannula for single use only) - 	
--	--	---	--

Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D

26. Euroflow (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D
27. Vivon (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D
28. Vivon Van (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D
29. Vivon Cath (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D
30. Nubenco (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D
31. Medcateter (Sterile Disposable IV Cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D
32. Vasocath (Sterile Disposable IV Cannuls for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D
33. Teflo (Sterile Disposable IV Cannula with Catheter and injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D
34. Neotef (Sterile Disposable IV Cannula with Catheter and without injection valve and with small wings for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D
35. Veincath (Sterile Disposable IV Cannula with Catheter for single use only)-Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D
36. Medin + (Sterile Disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D

- | | | | |
|--|--|--|--|
| | | <p>37. Steracath (Sterile disposable IV Cannula with catheter for single use only) – Polypropylene, Poly Carbonate, ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304, SS Spring AISI302) - Sterile – 5 - E & D</p> <p>38. Paedifix (Sterile Disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D</p> <p>39. Intrafix (Sterile Disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D</p> <p>40. Inflo Plus (Sterile Disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D</p> <p>41. Medin (Sterile Disposable Safety IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D</p> <p>42. Medcatheter Plus (Sterile Disposable safety IV Cannula with catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 - E & D</p> <p>43. Euroflow (Sterile disposable safety IV Cannula with catheter and with injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 - E & D</p> <p>44. Euroflow-1 (Sterile Disposable Safety IV Cannula with catheter and without injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 - E & D</p> <p>45. Neotec (Sterile Safety IV Cannula with Catheter and with injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 - E & D</p> <p>46. Neotec-1 (Sterile Safety IV Cannula with Catheter and without injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) – Sterile – 5 – E & D</p> <p>47. Canuvein (Sterile disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5</p> | |
|--|--|--|--|

– E & D

48. Galaxy (Sterile disposable IV Cannula with Catheter without port and without wings for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
49. Orocath Healflon™ (Sterile Disposable IV Cannula with Catheter and Injection Valve with small wings with in line one way valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 – E & D
50. Unoflon (Sterile Disposable Safety IV Cannula with Catheter and with injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 – E & D
51. Unique-In (Sterile disposable IV Cannula with Catheter for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D
52. Canuvein Safety (Sterile Disposable safety IV Cannula with catheter and with injection valve for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 – E & D
53. Healflon Safety (With injection valve and with wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D
54. Healvan Safety (Safety IV Cannula with Catheter and without injection valve and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D
55. Healcath Safety (Safety IV Cannula with Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D
56. Quciksafe (Safety IV Cannula with Catheter and Injection Valve with wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D
57. Easyfix (IV Cannula with Catheter with injection valve and with suture wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D

- | | | | |
|--|--|---|--|
| | | <p>58. Neofix (IV Cannula with Catheter with injection valve and with suture wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>59. Venfic Safety WP (Safety IV Cannula with Catheter and with injection Valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>60. Venfic Pen (IV Cannula with Catheter and without injection valve and without wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>61. Venfic WP (IV Cannula with Catheter with injection valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>62. Venfic (IV Cannula with Catheter and without injection valve and without wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>63. JCM MED SV CONCEPT (Sterile disposable IV Cannula with Catheter and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>64. Harson (IV Cannula with Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>65. Nextech Safety (Safety IV Cannula with Catheter with injection valve with wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D</p> <p>66. Safnule (Safety IV Cannula with Catheter and without injection Valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302 Porous plug) - Sterile – 5 – E & D</p> <p>67. Clear Safe Safety IV Cather (Safety IV Cannula with Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) - Sterile – 5 – E & D</p> <p>68. Kamemep (Sterile Disposable IV cannula for single use only) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone tube) - Sterile – 5 – E & D</p> <p>69. Kamemep (Sterile Disposable IV Cannula for single use only) -</p> | |
|--|--|---|--|

		<p>Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone tube) - Sterile – 5 – E & D</p> <p>70. Nextech (IV Cannula with Catheter and without injection valve and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>71. Nextech (IV Cannula with Catheter and without injection valve and without wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>72. Softmed (IV Cannula with Catheter with injection valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone tube) - Sterile – 5 – E & D</p> <p>73. Medsource IV Safety Catheter (Sterile Disposable Safety IV Cannula with Catheter and without injection Valve and without small wings for single use only(Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>74. Usafe (Safety IV Cnnula with Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip in AISI302) - Sterile – 5 – E & D</p> <p>75. Nervtec (IV Cannula with Catheter and With Injection Valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone tube) - Sterile – 5 – E & D</p> <p>76. Eli-Cath Hemodialysis Catheter Kit Dual Lumen - Polyurethane Nitinol & SS PP SS Needle ABS Non woven Fabric - Sterile – 3 – Dosmetic</p> <p>77. Eli-Cath Hemodialysis Catheter Kit Triple Lumen - Polyurethane Nitinol and SS PP SS Needle ABS Non woven fabric - Sterile – 3 – Dosmetic</p> <p>78. Harsoria triple Lumen Hemodialysis Catheter Kit – Polyurethane Nitinol & SS PP SS Needle ABS Non woven Fabric - Sterile – 3 – Dosmetic</p> <p>79. Single Lumen, Central Vanous Catheter – Polyurethane Silicone – Sterile – 3 - E& D</p> <p>80. Double Lumen, Central Venous Catheter - Polyurethane Silicone – Sterile – 3 - E& D</p> <p>81. Triple Lumen, Central Venous Catheter - Polyurethane Silicone – Sterile – 3 - E& D</p> <p>82. Clear Safe Comfort, Safety IV Catheter (Without Injection Valve and without small wings) – Polypropylene HDPE, ABS, Poly Acetal</p>	
--	--	--	--

		<p>Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip AISI302) – 5 - E & D</p> <p>83. Harsoria Two Lumen Hemodialysis Catheter Kit – Polyurethane Nitinol and SS PP SS Needle ABS Non woven Fabric - Sterile – 3 - E & D</p> <p>84. Harsoria Single Lumen Hemodialysis Catheter kit - Polyurethane Nitinol and SS PP SS Needle ABS Non woven Fabric - Sterile – 3 - E & D</p> <p>85. True Safe Comfort Safety I.V. Catheter – Sterile – 5 – E & D</p> <p>86. True Safe Safety I.V Catheter – Sterile – 5- E & D</p> <p>87. Global Care (IV Cannula with Catheter and without injection valve and without wings) – Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304) - Sterile – 5 – E & D</p> <p>88. Cadon Single Lumen Central Venous Catheter – Polyurethane Silicone – Sterile – 3 – E & D</p> <p>89. Cadon Double Lumen(Central Venous Catheter) – Polyurethane Silicone – Sterile – 3 – E&D</p> <p>90. Cadon Triple Lumen (Central Venous Catheter) - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>91. Cadon + Single Lumen Cental Venous Catheter -Polyurethane Silicone – Sterile – 3 – E&D</p> <p>92. Cadon + Double Lumen (Central Venous Catheter) - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>93. Cadon + Triple Lumen (Central Venous Catheter) - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>94. Harsoria Single Lumen (Central Venous Catheter) - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>95. Harsoria Double Lumen (Central Venous Catheter) - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>96. Harsoria Triple Lumen Central Venous Catheter - Polyurethane Silicone – Sterile – 3 – E&D</p> <p>97. Venfic Passive Safety 1SW(Safety IV Cannula with Catheter and without injection valve and with small wings - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 locking strip AISI302) - Sterile – 5 – E & D</p> <p>98. Venfic Safety II (Safety IV Cannula with Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 locking strip AISI302 Porous plug) - Sterile – 5 – E & D</p>	
--	--	---	--

- | | | | |
|--|--|--|--|
| | | <p>99. Truflon (IV Cannula with Catheter and with injection valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <p>100. Primefix (IV Cannula with Catheter and with injection valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <p>101. Truflon Safety (Safety IV Cannula with Catheter with injection Valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip in AISI302 - Sterile – 5 – E & D</p> <p>102. Healfon Insta - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube - Sterile – 5 – E & D</p> <p>103. Healfon Insta Plus - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone Tube) - Sterile – 5 – E & D</p> <p>104. Venfic SW (IV Cannula with Catheter and without injection valve and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <p>105. Vanfic Safety (Safety IV Cannula with Catheter and without injection valve and without wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <p>106. Truflon Insta (IV Cannula with injection Valve) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Silicone tube - Sterile – 5 – E & D</p> <p>107. Truflon Neo (IV Cannula with Catheter without injection valve and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 - Sterile – 5 – E & D</p> <p>108. Harsoriaa Guide Wire (Unocated) – Nitinol Core wire SS Spring – Sterile – 3 – E & D</p> <p>109. Harsoria Guide Wire (Coated) – SS Core Wire SS Spring with PTFE Coating) – Sterile – 3 – E & D</p> <p>110. Caress Safety (With Catheter and injection valve and with wings) – Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip in AISI302 – Sterile – 5 – E & D</p> <p>111. Caress Safety (With Catheter ad Without injection valve and with small wings) - Polypropylene HDPE ABS Poly Acetal Catheter</p> | |
|--|--|--|--|

		(PTFE/PUR/FEP SS Needle in AISI304 Locking strip in AISI302 – Sterile – 5 – E & D 112. Caress Safety (With Catheter and without injection valve and without small wings) - Polypropylene HDPE ABS Poly Acetal Catheter (PTFE/PUR/FEP SS Needle in AISI304 Locking strip in AISI302 – Sterile – 5 – E & D	
22.	M/s MIV Therapeutics (I) Pvt.Ltd., 136B, Surat Special Economic Zone,Daimond Parl,GIDC,Sachin Surat,Gujarat (India)	1.Yalimus Sirolimus Eluting Coronary Stent System (eight Brand) – L605 Cobalt Chromium, Biodegradable Polymer Based Sirolimus Eluting Coronary Stent	Mfg. Lic. No.- G/28/1261 Valid upto 05.10.2017 22-03-2017 (Additional Products)
23.	M/s Biotech Vision Care Pvt.Ltd.,Plot No.555,556,557,Khatraj-Vadsar Road,At & Post Khatraj,Tal.Kalol,Distt.Ga ndhi Nagar	1.YHSAS600,YHSAS 4X4,YHSAS600ROH – Natural Yellow Hydrophilic Acrylic Foldable Buttons/Sheet – Sterile – Export 2.PCK120NH,PKC125NH,PKC130NH – Hydrophilic Acrylic Foldable Buttons/Sheet- Sterile – Export 3.YHSA600,YHSAS 4X4,YHSAS600ROH – Natural Yellow Hydrophilic Acrylic Foldable Buttons/Sheet – Sterile – Domestic 4.PKC120NH,PKC125NH,PKC130NH – Hydrophilic Acrylic Foldable Buttons/Sheet – Sterile – Domestic 5. SQ605 – PMMA Sheet – Sterile – Domestic	Mfg. Lic. No.- G/1298 Valid upto 22-03-2017 (Additional Products)
24.	M/s Aurolab,No.1,Sivagangai Main Road,Veerapanjan,Madurai-625020	1.Aurobil – Poly hydroxyethyl Methacrylate – Sterile – 4 – E & D	Mfg. Lic. No.- D.DIS NO.969/D/13/2017(TN0000 2501) Valid upto 10-08-2019 22-03-2017 (Additional Products)
25.	M/s TTK Healthcare Limited (Ortho Division), No.290,SIDCO Industrial Estate, Ambattur,Chennai-98,Tamil Nadu	1.Buechel-Pappas High Flex Total Knee Replacement System – Co-Cr Alloy – Sterile – 3 – Domestic 2.Buechel Pappas High Flex Total Knee Replacement System – UHMWPE - Sterile – 3 – Domestic 3.Hip replacement System – PMMA- Sterile – 3 – Domestic	Mfg. Lic. No.- D.DIS NO.22009 (TN00002456)01/03/2016 Valid upto 12/07/2019 22-03-2017 (Additional Products)

26.	M/s La-Med Healthcare Pvt.LTd., (previously M/s lakhani Medicare Pvt.Ltd.),Plot No.136,Sector-24,Faridabad-121005	Alfa-med (I.V Cannula) – PP,LDPE,HDPE,FEP,PTFE,SS304,Silicon – Sterile – 5 – Export	Mfg. Lic. No654-B (H) Valid upto 31/07/2019 22-03-2017 (Additional Products)
27.	M/s Excel Optics Pvt.Ltd., No.15 A &B, Loganathan Nagar Second Street, Choolaimedu,Chennai-94	1.Excelens – Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 2.TrueLens - Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 3.Novalens - Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 4. Solaris - Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 5. Freedom - Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 6. Excelens - Poly Methyl Methacrylate (PMMA)– Sterile – 5 – E & D 7.Novaflex – Hydroxy Ethyl Methacrylate – Sterile – 3 – E & D 8. Truefold - Hydroxy Ethyl Methacrylate – Sterile – 3 – E & D 9.Truefold - Hydroxy Ethyl Methacrylate – Sterile – 3 – E & D	Mfg. Lic. No.- D.DIS NO.17493D1/01/2016 (TN00002843) Valid upto 30.09.2020 15-03-2017 (Additional Products)
28.	M/s Disposafe Health and Life Care Ltd.,Plot No.1 & 2,Sector-59,Phase-II,Ballabgarh-121004,Faridabad,Haryana,India	Samin – Polypropylene Polyethylene Polycetal Stainless Steel (AISI304) – Sterile – 5 - Export	Mfg. Lic. No.- 714-B(H) Valid upto 26-05-2016 22-03-2017 (Additional Products)
29	M/s Poly Medicure Ltd., Plot No.115-116,HSIDC,Sector-59,Ballabgarh,Distt.Faridabad,Haryana	1.KNOW-VEN I.V. Cannula with Injection port and wings – PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 2. KNOW-VEN I.V. Cannula with wings and without port – PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 3.FREBOR I.V Cannula with injection port with wings - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 4.Cat.Catheter (for veterinary Use only) - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 5.Dog Catheter (for veterinary use only) - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 6.Endoscope Flushing Catheter - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 7.Sensimedical I.V. Cannula / Catheter Intravenoso - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 8. Tro Venocath Flex - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 9. Safety I.V. Cannula with extension tube - PTFE/FEP/PUR PE,PP,SS,PA – Sterile – 5 – Export 10.Epidural Catheter Kit – PUR,PP,SS - Sterile – 5 – Export	Mfg. Lic. No.- 627-B (H) Valid upto 05-03-20121 15-03-2017 (Additional Products)

		<p>11. Alfasafe Catheter Intravenoso - PTFE/FEP/PUR PE,PP,SS,PA – Sterile for single used only – 5 – Export</p> <p>12. Alfasafe Catheter Intravenoso de Seguridad - PTFE/FEP/PUR PE,PP,SS,PA – Sterile for single used only – 5 – Export</p> <p>13. Epicent – PUR,PP,SS –Sterile for single used only – 5 – E & D</p> <p>14. Hospitech Plus - PTFE/FEP/PUR PE,PP,SS,PA – Sterile for single used – 5 – Export</p> <p>15. Novocent Protekt-Mono – PVC,ABS,PC,SS,PP – Sterile for single use only – 5 – E & D</p> <p>16. Novocent Protekt-Duo – PVC,ABS,PC,SS,PP – Sterile for single use only – 5 – E & D</p> <p>17. Novocent Protekt-Trio – PVC,ABS,PC,SS,PP – Sterile for single use only – 5 – E & D</p> <p>18. Novocent Protekt-Quadruple – PVC,ABS,PC,SS,PP – Sterile for single use only – 5 – E & D</p> <p>19. ART-CAN - PTFE/FEP/PUR PE,ABS,LDPE,SS,PC – Sterile for single used – 5 – E & D</p> <p>20. Instaflon Plus – PTFE/FEP/PUR,PE,PP,SS,PA Sterile for single used – 5 – Export</p> <p>21. Jelco Seriva - PTFE/FEP/PUR,PE,PP,SS,PA Sterile for single used – 5 – Export</p>	
30	M/s Biorad Medisys Pvt.Ltd.,48/3 &48/7,Sus-Pashan Road, Pune-411021	<p>1. Intramedullary Nails – SS316L/Titanium – Non Sterile – Domestic</p> <p>2. Bone Plate - SS316L/Titanium – Non Sterile – 5 – Domestic</p> <p>3. Bone Screw - SS316L/Titanium – Non Sterile – Domestic</p> <p>4. Nut & Bolt - SS316L- Non Sterile – Domestic</p> <p>5. Hip and Sholder Fixetion Prosthesis – SS316L- Sterile/Non Sterile -5 – Domestic</p> <p>6. Spine Fixetion Implants - SS316L/Titanium – Non Sterile – Domestic</p> <p>7. Locking System - SS316L/Titanium – Non Sterile – Domestic</p> <p>8. Mandible Plants 1.5/2.0/2.5mm - SS316L/Titanium – Non Sterile – Domestic</p> <p>9. Mendibal Screw - SS316L/Titanium – Non Sterile – Domestic</p> <p>10. Pins-Screws- SS316L- Non Sterile – Domestic</p> <p>11. Interlocking Bolt and Fragment Screws - SS316L/Titanium- Non Sterile – Domestic</p> <p>12. Guide Pin - SS316L- Non Sterile – Domestic</p> <p>13. Wires-Staples SS316L/Titanium - Non Sterile – Domestic</p>	<p>Mfg. Lic. No.- PD/124 Valid upto 30-07-2020</p> <p>15-03-2017 (Additional Products)</p>
31	M/s HCM Orthocare	1.HCM(Bone Screw) – SS316L/Titanium-5 – Non Sterile – E & D	Mfg. Lic. No.- G/28/1571

<p>Pvt.Ltd., 20 Jagannath Industrial Estate.Opp.Karma Gas Agency,Nr.Gangotri Estate,Odhav,Ahmedabad-382415</p>	<ol style="list-style-type: none"> 2.HCM (Cortex Screw) - SS316L/Titanium-5 – Non Sterile – E & D 3.HCM (DHS LEG Screw) - SS316L/Titanium-5 – Non Sterile – E & D 4.HCM (Interloc king screws) - SS316L/Titanium-5 – Non Sterile – E & D 5. HCM (PFN Leg Screw) - SS316L/Titanium-5 – Non Sterile – E & D 6.HCM (Cancellous Screw) - SS316L/Titanium-5 – Non Sterile – E & D 7.HCM (Cannulated Screw) - SS316L/Titanium-5 – Non Sterile – E & D 8.HCM (Locking Head Screw) - SS316L/Titanium-5 – Non Sterile – E & D 9.HCM HCM (Cancellous Screw with locking head) - SS316L/Titanium-5 – Non Sterile – E & D 10. HCM (Bone plates) - SS316L/Titanium-5 – Non Sterile – E & D 11. HCM (Dynamic Compression Plate) - SS316L/Titanium-5 – Non Sterile – E & D 12. HCM (Dyna Hip and Condylar) - SS316L/Titanium-5 – Non Sterile – E & D 13. HCM (Locking Compression Plate) - SS316L/Titanium-5 – Non Sterile – E & D 14.HCM (DHS LCP) - SS316L/Titanium-5 – Non Sterile – E & D 15. HCM (DCS Locking Plate) - SS316L/Titanium-5 – Non Sterile – E & D 16. HCM (Anatomical locking plates) - SS316L/Titanium-5 – Non Sterile – E & D 17. HCM (Distal Demural) - SS316L/Titanium-5 – Non Sterile – E & D 18. HCM (Distal Medial Tibia) - SS316L/Titanium-5 – Non Sterile – E & D 19. HCM (Proximal Lateral Tibia) - SS316L/Titanium-5 – Non Sterile – E & D 20. HCM (Raft Plate) - SS316L/Titanium-5 – Non Sterile – E & D 21. HCM (Clavical LCP) - SS316L/Titanium-5 – Non Sterile – E & D 22. HCM (Distal Humerous Dosola tral without support LCP) - SS316L/Titanium-5 – Non Sterile – E & D 23. HCM (Distal Humerous Dosola tral with support LCP) - SS316L/Titanium-5 – Non Sterile – E & D 24. HCM (Distal Medial Huemros) - SS316L/Titanium-5 – Non Sterile – E & D 25. HCM (Distal Humerous extra Articular) - SS316L/Titanium-5 – Non Sterile – E & D 26. HCM (Calcanium) - SS316L/Titanium-5 – Non Sterile – E & D 27. HCM (2.7 MM Distal Radius) - SS316L/Titanium-5 – Non Sterile – E & D 28. HCM (Philos LCP) - SS316L/Titanium-5 – Non Sterile – E & D 29. HCM (Olecron on LCP) - SS316L/Titanium-5 – Non Sterile – E & D 30.HCM (Reconstruction LCP) - SS316L/Titanium-5 – Non Sterile – E & D 	<p>Valid upto 24/05/2021</p> <p>16-03-2017 (Fresh License)</p>
--	--	--

- | | | |
|--|--|--|
| | <p>31. HCM (One-third LCP) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>32. HCM (Posterior Medical Tibia LCP) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>33. HCM (Proximal Tibia LCP) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>34. HCM (Tomofix LCP) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>35. HCM (Nails) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>36. HCM (Tibia Nails) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>37. HCM (Femur Inter Locking nails) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>38. HCM (Femoral proximal inter locking nail R & L) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>39. HCM (Femoral proximal inter locking nail PFN Short) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>40. HCM (Ten Nail) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>41. HCM (Rush Nail) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>42. HCM (Enders Nail) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>43. HCM (Syrus Tibia-Femur Nail) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>44. HCM (Expert Tibia Nails) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>45. HCM (Wire and Pins) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>46. HCM (K-Wire) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>47. HCM (ST-Pin) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>48. HCM (Schnews Screw) - SS316L/Titanium-5 – Non Sterile – E & D</p> <p>49. HCM (Spine) - Titanium-5 – Non Sterile – E & D</p> <p>50. HCM (Pedical Screw) - Titanium-5 – Non Sterile – E & D</p> <p>51. HCM (Mono Axial Screw) - Titanium-5 – Non Sterile – E & D</p> <p>52. HCM (Poly Axial Screw) - Titanium-5 – Non Sterile – E & D</p> <p>53. HCM (Single lock reduction screw) - Titanium-5 – Non Sterile – E & D</p> <p>54. HCM (Cervical Plate) - Titanium-5 – Non Sterile – E & D</p> <p>55. HCM (Cage System) - Titanium-5 – Non Sterile – E & D</p> <p>56. HCM (Tlif Cage) - Titanium-5 – Non Sterile – E & D</p> <p>57. HCM (Peek Spine Fusion Cage) - Titanium-5 – Non Sterile – E & D</p> <p>58. HCM (Rods) - Titanium-5 – Non Sterile – E & D</p> <p>59. HCM (Cervical Screw System) - Titanium-5 – Non Sterile – E & D</p> <p>60. HCM (Arthoscopic implants) - Titanium-5 – Non Sterile – E & D</p> <p>61. HCM (HIP Prosthesis) - Titanium-5 – Non Sterile – E & D</p> <p>62. HCM (BI-Polar) - Titanium-5 – Non Sterile – E & D</p> | |
|--|--|--|

32.	M/s Relisys Medical Devices Limited, Sy.No.312,Pocharam Road, Magalpally Village, Ibrahimpatnam Mandal,Ranga Reddy,Distt.Telangana	<p>1.Legend – Stainless Steel 316 LVM – Sterile – 5 – Export</p> <p>2.Stainless Steel Coronary Stent System – Stent (Stainless Steel 316LVM),Stent System (nylon) – Sterile – 3 – E & D</p> <p>3.Cobal + C^{VIINI} – Stent (Cobalt Chromium L 605), Stent System (Nylon) – Sterile – 3 –Export</p> <p>4. Cobal + C – Stent (Cobalt Chromium L 605), Stent System (Nylon) – Sterile – 3 –Export</p> <p>5.Paclitaxel Eluting Coronary Stent System - Stent (Stainless Steel 316 LVM), Stent System (Nylon), Coating Material (Polylactic Acid) – Sterile – 1-E& D</p> <p>6.Cobalt Chromium Coronary Stent Delivery System – Stent (Cobalt Chromium L 605) ,Stent System (Nylon) – Sterile – 3 - E & D</p>	<p>Mfg. Lic. No.- 19/RR/AP/2006/MD/R Valid upto 28-11-2021</p> <p>24-03-2017 (Renewal License)</p>

		<p>7.PTCA Catheter stent delivery device – Nylon – 3 - E & D 8. PTCA Catheter Dilatation – Nylon – 3 - E & D 9.Sirolimus Eluting Coronary Stent System - Stent (Cobalt Chromium L 605), Stent System (Nylon), Coating matrix (PLA Poly Lactic Acid) – Sterile – 2 – E & D 10.Release (R) - Stent (Cobalt Chromium L 605), Stent System (Nylon), Coating matrix (PLA Poly Lactic Acid) – Sterile – 2 – E & D 11.Everolimus Eluting Cobalt Chromium Coronary Stent System - Stent (Cobalt Chromium L 605), Stent System (Nylon), Coating matrix (PLGA Poly D,L-Lactide-Co-Glycolide) – Sterile – 2 – E & D 12.Pristine - Stent (Cobalt Chromium L 605), Stent System (Nylon), Coating matrix (PLGA Poly D,L-Lactide-Co-Glycolide) – Sterile – 2 – Export 13.Aspiration Catheter/Embolectomy Catheter – Pebax – Sterile – 3 – E & D 14.Coabl Chromium Stent (Co-Cr Stent) – Cobalt Chromium L 605 – Non Sterile – 5 – Export 15.Judkins Left,Judkins Right,Pig Tail Degree,Tiger Straight,Tiger Angles,Jacky,Amplatz Right,Amplatz Left – Braided Shaft (Pebax + 30% BaSO₄) – Sterile – 3 – E & D 16. Judkins Left,Judkins Right,Pig Tail Degree,Tiger,jacky,Amplatz - Braided Shaft (Pebax + 30% BaSO₄) – Sterile – 3 – Export 17.Guding Catheter - Braided Shaft (PTFE AND PEBAX) – Sterile – 3 – E & D 18.Pre-Dilatation Catheter – Nylon – Sterile – 3 - E & D 19.Post-Dilatation Catheter – Nylon – Sterile – 3 - E & D 20.Cantral Venous Catheter – Polyurethane - Sterile – 3 - E & D</p>	
33.	M/s Orthocare and Cure (I) Pvt.Ltd.,Plot No.22,S.No.327,B/1,Opp. Vedh Textile,At Post Mudshingi,Tal-Karvir,Dist-Kolhapure-416119	<p>1.Bone Screws – SS 316L/Titanium – Non Sterile – E & D 2.Bone Plates - SS 316L/Titanium – Non Sterile – E & D 3.Nails - SS 316L/Titanium – Non Sterile – E & D 4.Wire & Pins - SS 316L – Non Sterile – E & D</p>	<p>Mfg. Lic. No.- MH/102003 Valid upto 08.01.2022 28-03-2017 (Fresh License)</p>

34.	M/s General Surgical Company India Pvt.Ltd., No,43F Tass Industrial Estate,Amabattur,Chennai-98	<ol style="list-style-type: none"> 1.Pass Poserior and Anterior Spinal System – Titanium – Sterile – 5 – E& D 2. Pass Poserior and Anterior Spinal System – Stainless Steel – Sterile – 5 – E& D 3. Pass Poserior and Anterior Spinal System Alpha – Titanium – Sterile – 5 – E& D 4. Pass Poserior and Anterior Spinal System Alpha Plus – Titanium – Sterile – 5 – E& D 5.Caps Cervical Anterior plating system - Titanium – Sterile – 5 – E& D 6. Caps Plus Cervical Anterior plating system - Titanium – Sterile – 5 – E& D 7.K-Fixator multi-segmental Internal Fixator - Titanium – Sterile – 5 – E& D 8.Repoflx cervical and dorso lumbar - Titanium – Sterile – 5 – E& D 9.CVJ Loop-Grand Vertebral Junction Loop - Titanium – Sterile – 5 – E& D 10.Ravi`s Rectangle - Titanium – Sterile – 5 – E& D 11.Wires Kirschner - Titanium – Sterile – 5 – E& D 12. Wires Kirschner – Stainless Steel – Sterile – 5 – E& D 13.Posterior Occipital Cervico-thoracic fixation - Titanium – Sterile – 5 – E& D 14. Posterior Occipital Cervico-thoracic fixation (contd.) - Titanium – Sterile – 5 – E& D 15.Micro plates for cranial oral maxilla facial - Titanium – Sterile – 5 – E& D 16.Micro plates for cranial oral maxilla facial – Stainless Steel – Sterile – 5 – E& D 17.Bone Plates – Titanium - Sterile – 5 – E& D 18. Bone Plates (Contd.) – Titanium - Sterile – 5 – E& D 19. Bone Plates – Stainless Steel - Sterile – 5 – E& D 20. Bone Plates (Contd.)– Stainless Steel - Sterile – 5 – E& D 21.Bone Screws - Titanium - Sterile – 5 – E& D 22. Bone Screws (Contd.) - Titanium - Sterile – 5 – E& D 23. Bone Screws - Stainless Steel – 5 – E& D 24. Bone Screws (Contd.)- Stainless Steel – 5 – E& D 25.Cervical Disc – Titanium - Sterile – 5 – E& D 26.Lumbar disc - Titanium - Sterile – 5 – E& D 27.Expa Cage - Titanium - Sterile – 5 – E& D 28.Octa Cage Dorso/Lumbar - Titanium - Sterile – 5 – E& D 29.Octa cage cervical - Titanium - Sterile – 5 – E& D 30.Ethan-cervical anterior spacer - Titanium - Sterile – 5 – E& D 31.Ethan dorso lumbar anterior spacer - Titanium - Sterile – 5 – E& D 32.Fuzer cervical threaded cage - Titanium - Sterile – 5 – E& D 	<p>Mfg. Lic. No.- TN00002539 Valid upto 29.12.2018</p> <p>23-03-2017 (Renewal License)</p>
-----	---	---	--

- | | | |
|--|---|--|
| | <p>33.Fuzer dorso/lumbar threaded fusion cage - Titanium - Sterile – 5 – E& D</p> <p>34. Recage cervical - Titanium - Sterile – 5 – E& D</p> <p>35.Recage droso Lumbar - Titanium - Sterile – 5 – E& D</p> <p>36. Das-cervical - Titanium - Sterile – 5 – E& D</p> <p>37.Das Lumbar - Titanium - Sterile – 5 – E& D</p> <p>38. Bean Cage - Titanium - Sterile – 5 – E& D</p> <p>39. COFS Cannulated odontoid fixation system - Titanium - Sterile – 5 – E& D</p> <p>40.C-Fix Anterior disc cage with screw - Titanium - Sterile – 5 – E& D</p> <p>41. C-Fix Plus Anterior disc cage with screw - Titanium - Sterile – 5 – E& D</p> <p>42. GFS-Goel Facet Spacer - Titanium - Sterile – 5 – E& D</p> <p>43.GFS Goel Facet Prosthesis - Titanium - Sterile – 5 – E& D</p> <p>44. Gaas Goel atlanto axial system - Titanium - Sterile – 5 – E& D</p> <p>45.Atlanto axial prosthesis - Titanium - Sterile – 5 – E& D</p> <p>46. Das Ultra Distractable anterior spacer with holding device - Titanium - Sterile – 5 – E& D</p> <p>47.Alpha Fen Fenestrate D Pedicle system and cannula - Titanium - Sterile – 5 – E& D</p> <p>48.Alpha Ostero Ha coated pedicle screw system - Titanium - Sterile – 5 – E& D</p> <p>49.Rigid Pedicle Screw System - Titanium - Sterile – 5 – E& D</p> <p>50. Rigid Pedicle Screw System (contd.) - Titanium - Sterile – 5 – E& D</p> <p>51.RigidFen Fenestrated Pedicle System and cannula - Titanium - Sterile – 5 – E& D</p> <p>52. Rigid Osteo HA Coated pedicle screw system - Titanium - Sterile – 5 – E& D</p> <p>53 Cerlock Cervical Anterior Plating System - Titanium - Sterile – 5 – E& D</p> <p>54. Posterior/anterior C1/C2 Implants - Titanium - Sterile – 5 – E& D</p> <p>55.Uno Lock Pedicle Screw System (Contd.) - Titanium - Sterile – 5 – E& D</p> <p>56.Uno Plus Pedicle screw system - Titanium - Sterile – 5 – E& D</p> <p>57. Uno Plus Pedicle screw system(Contd.) - Titanium - Sterile – 5 – E& D</p> <p>58.Uno lock pedicle screw system – Stainless steel – 5 – E & D</p> <p>59. Uno lock pedicle screw system(Contd.) – Stainless steel – 5 – E & D</p> <p>60.Bullcage – Titanium – Sterile – 5 – E & D</p> <p>61.Plif Cage – Titanium – Sterile – 5 - E & D</p> <p>62.Tlif Cage - Titanium – Sterile – 5 - E & D</p> <p>63.Bean Ultra - Titanium – Sterile – 5 E & D</p> <p>64.Anterior spinal Implants - Titanium – Sterile – 5 - E & D</p> | |
|--|---|--|

65. Posterior Spinal Implants - Titanium – Sterile – 5 - E & D
66. Posterior Spinal Implants (Contd.) - Titanium – Sterile – 5 - E & D
67. Magnum minimal invasive spinal system - Titanium – Sterile – 5 - E & D
68. Magnum Plus minimal invasive spinal system - Titanium – Sterile – 5 - E & D
69. GSGM Ganga Spinal growth modulator - Titanium – Sterile – 5 - E & D
70. Infio CMF Cranial and maxilla facial/plastic/reconstructive/hand/surgery implants - Titanium – Sterile – 5 - E & D
71. Knee, Elbow, Hand, Cranial - Titanium – Sterile – 5 - E & D
72. Knee, Elbow, Hand, Cranial – Stainless steel – Sterile – 5 - E & D
73. Halder Humeral Nail - Stainless steel – Sterile – 5 - E & D
74. Halder Humeral Nail - Titanium – Sterile – 5 - E & D
75. Halifax Naul Inter Trochanteric Nail - Stainless steel – Sterile – 5 - E & D
76. Halifax Naul Inter Trochanteric Nail – Titanium – Sterile – 5 - E & D
77. Halifax Naul IntraTrochanteric Nail – Titanium – Sterile – 5 - E & D
78. Halifax Naul IntraTrochanteric Nail – Stainless steel – Sterile – 5 - E & D
79. Halder Tibial Nail - Stainless steel - Sterile – 5 - E & D
80. Halder Tibial Nail (Contd.)- Stainless steel - Sterile – 5 - E & D
81. Halder Tibial Nail - Titanium - Sterile – 5 - E & D
82. Halder Femoral Nail - Stainless steel - Sterile – 5 - E & D
83. Halder Femoral Nail - Titanium - Sterile – 5 - E & D
84. Halder Femoral Nail(contd.) - Titanium - Sterile – 5 - E & D
85. Plates - Titanium - Sterile – 5 - E & D
86. Screws - Titanium - Sterile – 5 - E & D
87. Cygnet Pediatric Scoliosis correction and fixation system - Titanium - Sterile – 5 - E & D
88. Cygnet (Contd.) - Titanium - Sterile – 5 - E & D
89. Direct lateral interbody cage DLIF - Titanium - Sterile – 5 - E & D
90. Alif Anteriorlumbar Fusion Cage - Titanium - Sterile – 5 - E & D
91. Inter Spinous Process - Titanium - Sterile – 5 - E & D
92. K-Cer Cervical Anterior Plating system – Titanium – Sterile – 5 – E & D
93. Laminofix Laminoplasty implant - Titanium - Sterile – 5 - E & D
94. Wire - Titanium - Sterile – 5 - E & D
95. Wire – Stainless Steel – Sterile – 5 – E & D
96. Denique Pedicle Screw System – Titanium – Sterile – 5 – E & D
97. Denique Pedicle Screw System (Contd.) – Titanium – Sterile – 5 – E & D
98. Facet Link System for facet Stabilization - Titanium – Sterile – 5 – E & D
99. Dynamic Facet Cage - Titanium – Sterile – 5 – E & D
100. Static Facet Cage - Titanium – Sterile – 5 – E & D

		<p>101.Sacrolitic Fixation System - Titanium – Sterile – 5 – E & D</p> <p>102.Lateral Lumbar cage Interbody LLIF - Titanium – Sterile – 5 – E & D</p> <p>103.Nail – Stainless Steel - Sterile – 5 – E & D</p> <p>104.Nail - Titanium – Sterile – 5 – E & D</p> <p>105.Expandable Interbody Cage-Cresco – Titanium – Sterile – 5 – E & D</p> <p>106.Thoracic Lumbar Plate with Screw - Titanium – Sterile – 5 – E & D</p> <p>107. Cervical Anterior Plate –LP - Titanium – Sterile – 5 – E & D</p> <p>108. Cervical Anterior Distractable Plate - Titanium – Sterile – 5 – E & D</p> <p>109. Posterior Occipital Cervico Thoracic Fixation - Titanium – Sterile – 5 – E & D</p>	
35.	M.s Angiplast Private Limited,Plot No.4803 & 4804 Phase –IV, G.I.D.C. Vatva,Ahmedabad-382445,Gujarat	<p>1.Sterile Disposable Single use Haemodialysis Blood Tubing Set :Angeltouch – TT -12A “ – Discharge Hook(PP), Rotating Lucer Lock(ABS),C Clamp (HDPE), Injection Port(ABS),Stopper(Latex Fre),Pump Connector(PVC),Pump tube (PVC),Female lucer lock(ABS),Cap(PP),Chamber Cap(PVC),Chamber(PVC),Dialyser Connector(ABS), Cap(PP),Tubing(PVC) – Sterile –Five years from date of manufacturing</p> <p>2. Sterile Disposable Single use Haemodialysis Blood Tubing Set “ ANGELTOUCH-TT-12B” – The sterile Haemodialysis blod line is used for dialysis of blood in patients whose kidney fails to purify flood - Sterile –Five years from date of manufacturing</p> <p>3. Sterile Disposable Single use IV Cannula with wings and injection port “ANGELTOUCH-IT-31” – The Sterile IV Cannula is used infusion or transfusion on continous /intermittent basis and short term peripheral catheterization of vein to give access for drug delivery - Sterile –Five years from date of manufacturing</p>	<p>Mfg. Lic. No.- G/28/1357</p> <p>Valid upto 24-11-2020</p> <p>30-03-2017</p> <p>(Renewal License)</p>
36.	M/s Surgimade having manufacturing premises at Plot No.459/15/B,GIDC,Porra mangamdi, At:Por,Dist.Vadodara	<p>1.Humerus nail Cannulated – SS316L/Titanium/SS316LVM – Non Sterile – E & D</p> <p>2.Humerus Nail Solid - SS316L/Titanium/SS316LVM – Non Sterile – E & D</p> <p>3.Tibia Nail Cannulated - SS316L/Titanium/SS316LVM – Non Sterile – E & D</p> <p>4.Tibia Naild Solid - SS316L/Titanium/SS316LVM – Non Sterile – E & D</p> <p>Femur Nail Cannulated - SS316L/Titanium/SS316LVM – Non Sterile – E & D</p>	<p>Mfg. Lic. No.- G/28/1591</p> <p>Valid upto 03-01-2022</p> <p>30-03-2017</p> <p>(Fresh License)</p>
37.	M/s Unisys Ortho,51,52,53,Bhagirath -4 Industrial Estate,Opp.	<p>1.Self Compression Plate 3.5mm – Stainless Steel/Titanium – Non Sterile – E & D</p> <p>2.Limited Conatct Self Compression Plate 3.5mm - Stainless Steel/Titanium</p>	<p>Mfg. Lic. No.- G/28/1319</p> <p>Valid upto 17-11-2019</p>

	Jawaharnagar, Near Nageshwar Estate, Amaraiwade, Ahmedabad-380026, Gujarat	– Non Sterile – E & D 3.Small T Plate 3.5mm - Stainless Steel/Titanium – Non Sterile – E & D 4.Small Oblique T Plate 3.5mm Left and Right - Stainless Steel/Titanium – Non Sterile – E & D 5.Cloverleaf Plate 3.5mm - Stainless Steel/Titanium – Non Sterile – E & D 6.Y Plate 3.5mm - Stainless Steel/Titanium – Non Sterile – E & D 7. Reconstruction Plate 3.5mm - Stainless Steel/Titanium – Non Sterile – E & D 8.One Third Tubular Plate 3.5mm - Stainless Steel/Titanium – Non Sterile – E & D 9.One third Locking Plate - Stainless Steel/Titanium – Non Sterile – E & D	30-03-2017 (Renewal License)
--	--	---	---------------------------------

S. No.	Name of the Manufacturer	Product Licensed by CLAA	File No.	Manufacturing Licence No. & Validity
38.	M/s Surgimedik Healthcare, Surgimedikk House, Plot No.223, Surat Special Economic Zone, Sachin, Surat-394230	Ureteral Stent- Polyurethane-Sterile-5- E& D	CLAA/MD/GUJ/09/2009-DC	Mfg Licno.G/28/1384 validity period 17/08/2016 to (05-04-2017) (Approval for additional Product)
39.	M/s Neo Surgical.Plot No.G-306,GIDC Almighty Gate, Metoda,Pin-360021,Rajkot	<ol style="list-style-type: none"> 1. Malleolar Screw – SS316L,S316LVM,Titanium-Non Sterile-5-E& D 2. Herbert Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 3. PFN Leg Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 4. DHS Leg Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 5. Top Screw used with DHS Leg Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 6. MonoAxial Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 7. Polyaxial Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 8. Mono Reduction Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 9. Poly Reduction Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 10. Multiaxial Cannulated Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 11. Maspoly Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 12. Bone Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 13. Occipital Bone Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 14. ACL Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 15. PCL Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 16. Cranio Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 17. Emergency Screw- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 18. Semi Plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 19. One Third Tebular Plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 20. Reconstruction plate-round hole- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 21. Reconsturction locking plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 22. “Y”Reconstruction ion plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 23. “T” Small locking plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D 	CLAA/MD/GUJ/01/2016-DC	Mfg Licno.G/28/1585384 validity period 22/09/2016 to 21/09/2021 (05-04-2017) (Approval for additional Product)

		<p>24. "T"Olique locking plate left and right- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>25. "T"Small- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>26. "T" Olique- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>27. Distal Radius Plate-Extra Atricular- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>28. Distal Radius Plate-Straight- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>29. Distal Radius Angeled- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>30. Distal Radius Plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>31. Distal Radius Plate-L Olique- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>32. Distal Radius Plate-T- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>33. Distal Radius Plate-Buttress- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>34. Distal Radius Plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>35. Humerous Locking Plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>36. Locking distal humerus plate(dorsolateral, with support, without support, medial)- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>37. Locking Proximal humerus plate (Philos)plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>38. Olacranon plate left and right - SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>39. Extra Atricular distal humerus locking plate left and rught- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p> <p>40. Metaphyseal distal medial humerus plate- SS316L,S316LVM,Titanium-Non Sterile-5-E& D</p>		
40.	M/s HAROSIA Healthcare Pvt.Ltd.,110-111,Phase-IV, Udyog Vihar,Gurgaon-122015,Haryana,India	<p>1.Amieflow-Polyurethane Silicone- Sterile- 3- E & D</p> <p>2.Amieflow (Double Lumen Central Venous Catehter)- Polyurethane Silicone - Sterile- 3- E & D</p> <p>3. Amieflow (Triple Lumen Central Venous Catheter) - Polyurethane Silicone - Sterile- 3- E & D</p>	CLAA/MD/HR/04/2006-DC (Pt-2)	<p>Mfg Licno/DML No.,637-B (H) validity period 09-2016 to 18-09-2021</p> <p>(05-04-2017) (Approval for additional Product)</p>
41.	M/s MediPlus India Ltd.,1261,1262 MIE,Bahadurgarh,Dist.Jhajjar,Haryana-124507	Neotec-1 – Polypropylene,Polythylene,Polyacteal,Stainless Steel (AISI304)- Sterile- 5- E &D	CLAA/MD/HR/06/2006-DC	Mfg Licno No.,646-B (H) validity period 01/06/2012 to 31/05/2017

				(05-04-2017) (Approval for additional Product)
42.	M/s Mais India Medical Devices Pvt.Ltd.,Plot No.525P, Sector-37, Pace City-II, Gurgaon-122001, Haryana	1.Liberton- Plastic raw material,Catheter,SS Needle- Sterile – 5- E 2.Healthhease- Plastic raw material,Catheter,SS Needle- Sterile – 5- E	CLAA/MD/HR/03/2013-DC	Mfg Licno No.,726- B (H) validity period 15-05-2014 to 14-05-2019 (05-04-2017) (Approval for additional Product)
43.	M/s Meril Life Sciences Pvt.ltd., Survey No.135/2a and 1352/B,174/2 (M1-M2,Meril Park, Fourth Floor),Near G.M.Bilakhia Stadium, Muktanand Marg, Chala,Vapi,Valsad-396191,Gujarat,India	Meriskin™(PTFE Covered Stent graft)- L605 Cobalt Chromium stent covered with PTFE (Polytetrafluoroethylene) – Sterile- 2 – E & D	CLAA/MD/GUJ/01/2008-DC (Pt.IV)	Mfg Licno No.,G/28/1273 validity period 20-06-2013 to 19-06-2018 (05-04-2017) (Approval for additional Product)
44.	M/s Bonetech Medisys Pvt.ltd.,50,51,52,53 Sarovar Industrial Park,Opp.Steel Konnect India Pvt.ltd., S.P.Ring Road,Hathijan-Vatva,Ahmedabad	1.Bonetech (Medisys Pvt.ltd) – Stainless Steel & Titanium – Non Sterile- E & D	CLAA/MD/GUJ/08/2015-DC	Mfg Licno No.,G/28/1273 validity period 06-05-2016 to 05-05-2021 (11-04-2017) (Approval for additional Product)
45.	M/s Adler Mediequip Pvt.Ltd.,Plot No.A-1,MIDC Sadavali.Tal.Sangameshwar,Dist:Ratnagiri-415804,Maharsashtra	1.Adler (Bone Screws)-Stainless Steel (ISO 5832-1/9), Titanium (ISO 832-2), Titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile &Non-Sterile- 5 – E & D 2. Adler (Bone Plates)-Stainless Steel (ISO 5832-1/9), Titanium (ISO 832-2), Titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile &Non-Sterile- 5 – E & D 3. Adler (Pins and Wires)-Stainless Steel (ISO 5832-1/9), Titanium (ISO 832-2), Titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile &Non-Sterile- 5 – E & D 4. Adler (BoHemi Endoprosthesis)-Stainless Steel (ISO 5832-1/9), Titanium (ISO 832-2), Titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile &Non-Sterile- 5 – E & D 5. Adler (Joint Replacement Prosthesis)-Titanium ALLOY (ISO 5832-3/11,ASTM F1108),Cobalt Chrome (ISO 5832-4-8/12),Stainless Steel (AISI 316L Gr X02Cr17Ni12Mo2-IS 6911/IS 6528/ISO 5832-1/9),Hydroxyapatite (ISO 13/79),Peek (ASTM F2026),Tantalum (ISO	CLAA/MD/MS/02/2016-DC	Mfg Licno No.,KD-649 validity period 17-09-2016 to 16-09-2021 (17-04-2017) (License Form 28)

13782) - Sterile & Non-Sterile- 5 – E & D

6.Sushrut/S`RUT(Bone Screws)- Stainless Steel (AISI 316L GrX02Cr11Ni12Mo2-IS 6911/IS 6528)- Sterile & Non-Sterile- 5 – E & D

7.Sushrut/S`RUT(Pins and Wires)- Stainless Steel (AISI 316L GrX02Cr11Ni12Mo2-IS 6911/IS 6528)- Sterile & Non-Sterile- 5 – E & D

8.Sushrut/S`RUT(IM Nails)- Stainless Steel (AISI 316L GrX02Cr11Ni12Mo2-IS 6911/IS 6528)- Sterile & Non-Sterile- 5 – E & D

9.Sushrut/S`RUT(Hemi Endoprosthesis)- Stainless Steel (AISI 316L GrX02Cr11Ni12Mo2-IS 6911/IS 6528)- Sterile & Non-Sterile- 5 – E & D

10.UMEX(Pins and Wires)- Stainless Steel (ISO 5832-1/9)- Sterile and Non Sterile- 5- E & D

11.Eazee (Bone Screws)- Stainless Steel (ISO 5832-1), Titanium (ISO 5832-2),Titanium Alloy ISO 5832-3/ISO 5832-11)- Sterile and Non Sterile -5- E & D

12.Eazee (IM Nails)- Stainless Steel (ISO 5832-1), Titanium (ISO 5832-2),Titanium Alloy ISO 5832-3/ISO 5832-11)- Sterile and Non Sterile -5- E & D

13.Zeta (Spine Implants)- Titanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),PeeK (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9),Cobalt Chrome Allo (ISO 5832-4-8/12),Tantalum (ISO 13782) - Sterile and Non Sterile -5- E & D

14.Onelock (Bone Screws) - Titanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),PeeK (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9),Cobalt Chrome Allo (ISO 5832-4-8/12),Tantalum (ISO 13782) - Sterile and Non Sterile -5- E & D

15.Onelock Spine Implants) - Titanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),PeeK (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9),Cobalt Chrome Allo (ISO 5832-4-8/12),Tantalum (ISO 13782) - Sterile and Non Sterile -5- E & D

16.Endofix (Bone Plates)- Totanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),PeeK (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9)

17.Atlas (Cone Screws)- Stainless Steel (AISI 316L Gr X02Cr17Ni2Mo2-IS 6911/IS 6528/ISO 5832-1/9),Titanium Alloy (ISO 5832-3/11),Titanium (ISO 5832-2)-Sterile and Non Sterile – 5- E & D

18.Atlas (IM Nails)- Stainless Steel (AISI 316L Gr X02Cr17Ni2Mo2-IS 6911/IS 6528/ISO 5832-1/9),Titanium Alloy (ISO 5832-3/11),Titanium (ISO 5832-2)-Sterile and Non Sterile – 5- E & D

19.ResTOR (Joint Replacement Prosthesis)- titanium Alloy (ISO 5832-3/11,ASTM F1108), Cobalt Chrome (ISO 5832-4-8/12),Stainless Steel (AISI 316L GrX02Cr17Ni12Mo2-is 6911/is 6528/iso 5832-1/9),Hydroxyapatite (ISO 13779)- Sterile and Non Sterile – 5- E & D

		<p>20.Accis (IM Nails)- Stainless Steel (AISI 316L GrX02Cr17Ni12Mo2-IS 6911/IS 6528/ISO 5832-1/9),Titanium Alloy (ISO 5832-3/11),Titanium (ISO 5832-2)- Sterile and Non Sterile- 5- E & D</p> <p>21.Axilock (bone Screws) – Stainless Steel (ISO 5832-1/9),Titanium (ISO 5832-2),titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile and Non Sterile – 5- E & D</p> <p>22.Axilock (Bone Plates) – Stainless Steel (ISO 5832-1/9),Titanium (ISO 5832-2),titanium Alloy (ISO 5832-3/ISO 5832-11)- Sterile and Non Sterile – 5- E & D</p> <p>23.Cerviedge (Spinal Plates)-Titanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),Peek (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9),Cobalt Chrome Alloy (ISO 5832-4-8/12)- Sterile and Non Sterile – 5- E &D</p> <p>24.CerviTwin (Spinal Plates)- Titanium (ISO 5832-2),Titanium Alloy (ISO 5832-3/ISO 5832-11),Peek (ASTM F2026),Stainless Steel (ISO 5832-1/ISO 5832-9),Cobalt Chrome Alloy (ISO 5832-4-8/12)- Sterile and Non Sterile – 5- E &D</p> <p>25.Moduloc (Hemi Endoprosthesis)-Stainless Steel (AISI 316L GrX02Cr17Ni12M02-IS 6911/IS 6528/ISO 5832-1/9),Cobalt Chrome Alloy (ISO 5832-4-8/12),UHMWPE (ISO 5834-1/2)- Sterile and Non Sterile – 5- E & D</p> <p>26.Adler Legend (Joint Replacement Prosthesis)- Titanium Alloy (ISO 5832-3/11), Hydroxyapatite (ISO 13779)- Sterile and Non sterile – 5- E& D</p> <p>27.Adler Legend (Joint Replacement Prosthesis)- Stainles Steel (ISO 5832-1/9),Cobalt Chrome Alloy (ISO 5832-4-8/12) Sterile and Non sterile – 5- E& D</p> <p>28.Atlas Proxima (IM Nails) – Stainless Steel (AISI 316L GrX02Cr17Ni12Mo2-IS 6911/IS 6528/ISO 5832-1/9),Titanium Alloy (ISO 5832-3/11),Titanium (ISO 5832-2)- Sterile and Non sterile – 5- E& D</p> <p>.</p>		
46.	M/s Miraculus Orthotech Pvt.Plot No.1112,Phase-3, GIDC, Vatva,Ahmedabad-382445,Gujarat (india)	<p>1.Locking Screw- SS 316/LVM/Titanium- Non Sterile – E & D</p> <p>2.Extra-Articular Distal Humerus Plate- SS 316/LVM/Titanium- Non Sterile – E & D</p> <p>3.ML-Tibia-Titanium- Non Sterile – E & D</p>	CLAA/MD/GUJ/03/2013-DC	Mfg Licno No.,.G/28/1449 validity period 28-06-2013 to

- 4.ML-Titbia End Cap - Titanium- Non Sterile – E & D
- 5.Headless Comprassion Screw- Titanium- Non Sterile – E & D
6. ASLP 3.5 Supreior Anterior Clavicle Plates - SS 316/LVM/Titanium- Non Sterile – E & D
7. 3.5mm ASLP Superior Anterior Clavicle plate with lateral extention - SS 316/LVM/Titanium- Non Sterile – E & D
- 8.ASLP-3.5 Periarticular Proximal humerus plates -SS 316/LVM/Titanium- Non Sterile – E & D
- 9.ASLP- 3.5 Posteromedial Promximal Tibia Plates - SS 316/LVM/Titanium- Non Sterile – E & D
10. ASLP- 3.5/4.5/5.0 Metaphyseal Plates - SS 316/LVM/Titanium- Non Sterile – E & D
11. ASLP- 3.5 Proximal medical tibia plates - SS 316/LVM/Titanium- Non Sterile – E & D
- 12.Pelvic Curve Plate- SS 316/LVM/Titanium- Non Sterile – E & D
13. Pelvic Straight Plate- SS 316/LVM/Titanium- Non Sterile – E & D
14. Pelvic J Plate - SS 316/LVM/Titanium- Non Sterile – E & D
15. CTL Mesh Cage – Titanium – Non Sterile – E & D
- 16.ACC- Stand alone Cage – Titanium + Peek – Non Sterile – E & D
17. 3mm Locking screws – Tiatnium – Non sterile – E & D
18. Simple femur nail – Tatanium – Non Sterile – E & D
19. Variable Angle, Extra Articular – Titanium – Non Sterile – E & D
20. Variable Angle,Two Colum-
- 21.MS-2 Poly asixal Screws Dual threaded – Titanium – Non Sterile – E & D
22. Posterior Cervical Screws – Titanium – Non Sterile – E & D
23. Posterior cervical screws, Rod – Titanium – Non Sterile – E & D
24. Phenistrated SCREWS – Titanium – Non Sterile – E & D
- 25.M.I.S. Rod – Titanium – Non Sterile – E & D
26. Posterior Cervical Inner – Titanium – Non Sterile – E & D
- 27.ASLP-2.7/3.5,Distal Humerus,Posterolateral with lateral support plates – SS 316 LVM /Titanium- Non Sterile – E & D
- 28.ASLP-3.5,Distal humerus Medical Plates- SS 316 LVM/Titanium- Non Sterile – E & D
- 29.ASLP-Metaphyseal Plates 3.5 for Distal,Medical Humerus- SS 316 LVM/Titanium- Non Sterile – E & D
- 30.ASLP-3.5 Olecranon Plates - SS 316 LVM/Titanium- Non Sterile – E & D
31. ASLP-4.5/5.0, DHS Plates with long barrel (38 mm)- SS 316 LVM/Titanium- Non Sterile – E & D
32. ASLP-3.5,Proximal Latreal Tibia Plates - SS 316 LVM/Titanium- Non Sterile – E & D

(18-04-2017)
 (Approval for additional Product)

		<p>33. ASLP-3.5 Distal Medial Tibia plates - SS 316 LVM/Titanium- Non Sterile – E & D</p> <p>34. ASLP 3.5 Distal Anterolateral tibia plates - SS 316 LVM/Titanium- Non Sterile – E & D</p> <p>35. ASLP 3.5 Proximal medial tibia plates - SS 316 LVM/Titanium- Non Sterile – E & D</p> <p>36. 3.5, One third tubular plates (non locking)- SS 316 LVM/Titanium- Non Sterile – E & D</p> <p>37. ASLP 3.5 Reconstruction plates - SS 316 LVM/Titanium- Non Sterile – E & D</p> <p>38. SL Tibia- Titanium- Non Sterile – E & D</p> <p>39. Miraculus PFA Nail - Titanium- Non Sterile – E & D</p> <p>40. PFA Blade Nail - Titanium- Non Sterile – E & D</p> <p>41. 5mm Locking Bolt- Titanium- Non Sterile – E & D</p> <p>42. End Cap for PFA nail - Titanium- Non Sterile – E & D</p> <p>43. 4.5mm Cancellous locking bolt- Titanium- Non Sterile – E & D</p> <p>44. End Cap for tibia II Nail - Titanium- Non Sterile – E & D</p> <p>45. 4.5 mm Cannulated screws, partially/fully threaded- SS 316 lvm/Titanium – Non Sterile – E & D</p> <p>46. 6.5mm cannulated screws, partially/fully threaded - SS 316 lvm/Titanium – Non Sterile – E & D</p>		
47.	M/s Meril Helathcare Pvt.ltd., Survey o.135/2/B & 74/2, 1 st floor, H1-H3, Meril Park, Muktanand Marg, Chala, Vapi-396191	<p>1. Bone Plates (KET)- SS316LVM – Non Sterile – E & D</p> <p>2. Bone Plates (DHUM)- SS316L – Non Sterile – E & D</p> <p>3. Bone Screws (KET)- SS316LVM- Non Sterile – E & D</p> <p>4. Bone Screws (DHUM) – SS316L – Non Sterile – E & D</p>	CLAA/MD/GUJ/05/2013-DC	<p>Mfg Licno No., G/28/1465 validity period 01.02.2014 to 31.01.2019</p> <p>(18-04-2017) (Approval for additional Product)</p>
48.	M/s Bio-Med Helathcare Care Product Pvt.ltd., Plot No.62-65, Sector-6A, SIDCUL, Ranipur, Haridwar-249403, Uttarakhand	<p>1. Ritemed- Polypropylene, stainless steel (14-26G) – Sterile – Export</p> <p>2. Manaflo- Polypropylene, stainless steel (14- 26G) – Sterile – Export</p>	CLAA/MD/UA/02/2007-DC	<p>Mfg Licno No/D.L.No.123/UA/SC/P-2007</p> <p>(24-04-2017) (Approval for additional Product)</p>

49.	M/s Poly Medicure Ltd., Plot No.115-116,HSIDC, Sector-59, Ballagarh, Distt.Faridabad-121004,Haryana,India	<ol style="list-style-type: none"> 1.Polysafety-BC- PTFE/FEP/PUR, PE/PP/SS/PA (14-27G)- Sterile – 5- E & D 2.ADVA Catheter Para infusao Intravenosa De Uso Periferico Sobre Agulha - PTFE/FEP/PUR, PE/PP/SS/PA (14-27G)- Sterile – 5- E & D 3. ADVA Dispositivo Seguranca Cateter Para Infusao Intracenososa De Uso Periferico Sobre Agulha com Dispositivo De Seguranca-PTFE/FEP/PUR, PE/PP/SS/PA (14-27G)- Sterile – 5- Export 4. Polyflon-Pro- PTFE/FEP/PUR/PuMEX/PE/PP/SS/PA (14-27G)- Sterile – 5- E & D 5.PTFE/FEP/PUR/PuMEX/PE/PP/SS/PA (14-27G)- Sterile – 5 – E & D 	CLAA/MD/HR/01/2006-DC	Mfg Licno No.,627-B (H) (24-04-2017) (Approval for additional Product)
50.	M/s Zedex Orthocare Pvt.ltd., 1 st Floor, Shop No.1 to 11, Krishna Corporation Coumpund, M/s Associates Chemical, Near Municipal Labour Quarters,Naroda,Ahmedabad-380025	<ol style="list-style-type: none"> 1.Bone Screw- Zedex- SS316L/Titanium LVM- Non Sterile – E& D 2.Bone Plate- Zedex- SS316L/Titanium LVM- Non Sterile – E& D 3.Interlocking Nail- Zedex-SS316L /Titanium LVM- Non Sterile – E& D 4.Locking Plate- Zedex-SS316L /Titanium LVM- Non Sterile – E& D 5.Wire and Pins- Zedex-SS316L /Titanium LVM- Non Sterile – E& D 6.Maxillofacial Implants - Zedex-SS316L /Titanium LVM- Non Sterile – E& D 	CLAA/MD/GUJ/06/2016-DC	Mfg Licno No. G/28/1580 validity period from 09.08.2016 to 08.08.2021 (25-04-2017) (Approval for additional Product)
51.	M/s Suhradam Ortho, Survey No.183, B/H, Satyam Mall,Nr.Cadila Pharmaceutical Ltd., 1 st floor, Shop No.1 to 11, Krishna Corporation Compound, Sakhej- Bavla Highway- Changdor-3882213,Tal:Sanand,Dist:A hmedabad	1.Bone Screws- SS316L/Titanium/SS316LVM- Non Sterile – E & D	CLAA/MD/GUJ/03/2017-DC	Mfg Licno No. G/28/1606 validity period from (01-05-2017) (Approval for additional Product)
52.	M/s Meril Life Sciences Pvt.ltd., Survey no.135/2a and 135/2b,174/2 (M1-M2,Meril Park, Fourth Floor),near G.M.Bilakhia Stadium,Muktanand Marg,Chala,Vapi,Valsad-396191,Gujarat	1. Rapid Exchange Percutaneous Transluminal Coronary Angioplasty Balloon Dilatation Catheter (Hydrophilic Coated)- Xpedient- Non Sterile – 4- E & D	CLAA/MS/GUJ/01/2008-DC (Pt.IV)	Mfg Licno No.G/28/1273 validity from 20.06.2013 to 19.06.2018 (04-05-2017) (Approval for additional Product)
53.	M/s Ribeel International Limited, 20 th Mile,Jatheri	1.Harsoria Foley Catheter- Natural Latex-Sterile/Non-Sterile- 5- E &D	CLAA/MD/HR/03/2014-DC	Mfg Licno No.740-B (H)

	Road, P.O.Rai,Sonepat			(04-05-2017) (Approval for additional Product)
54.	M/s iLife Medical Devices Pvt.Ltd., #86,Sector-56,Phase-IV, I.E.Kundli,Sonepat-131028,Haryana	1.VEX-Catheter,Polypropylene,HDPE,LDPE &Stainless Steel- Sterile -5-dOMESTIC	CLAA/MD/HR/03/2009-DC	Mfg Licno No.673-B (H) (04-05-2017) (Approval for additional Product)
55.	M/s Harsoria Healthcare Pvt.ltd., 110-111,Phase-Iv, Udyog Vihar, Gurgaon-122015,Haryana	<p>1.Elite Harsoria (Hemodialysis catheter kit Single Lumen)- Polyurethane Nitinol and SS PP SS Needle ABS Non Woven Fabric- Sterile -3- E & D</p> <p>2.Elite Harsoria (Hemodialysis catheter kit Dual Lumen)- Polyurethane Nitinol and SS PP SS Needle ABS Non Woven Fabric- Sterile -3- E & D</p> <p>3.Elite Harsoria (Hemodialysis catheter kit TRIPLE Lumen)- Polyurethane Nitinol and SS PP SS Needle ABS Non Woven Fabric- Sterile -3- E & D</p> <p>4.Elite Harsoria (Uncoated)- Ntinol core wire SS Spring- Sterile -3- E & D</p> <p>5.Elite Harsoria (PTFE Coated)- SS Core Wire SS Spring with PTFE Coating – Sterile – 3- E & D</p> <p>6.IV Cannula with catheter with injection Valve – Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p> <p>7.IV Cannula with catheter without injection Valve – Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p> <p>8.IV Cannula with catheter & without injection valve and without small wings – Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p> <p>9.IV Cannula with catheter & with injection valve– Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p> <p>10.IV Cannula with catheter & without injection valve and with small wings – Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p> <p>11.IV Cannula with catheter & without injection valve and without small wings – Polypropylene HDPE,ABS, Poly Acteal ,Catheter (PTFE/PUR/FEP/ETFE/PUR/SS Needle in AISI304 Silicone tube – Sterile -5- E& D</p>	CLAA/MD/HR/04/2006-DC (Pt-2)	Mfg Licno No.637-B (H) Validity from 19-09-2016 to 18-09-2021 (04-05-2017) (Approval for additional Product)

56.	M/s Excel Optics Pvt.ltd., No.15 A & B, Loganathan Nagar Second Street, Choolaimedu,Chennai-94	<p>1.Truevue (Hydrophobic Foldable Arcylic Diffractive multifocal intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p> <p>2.Novaflex(Hydrophobic Foldable Arcylic Diffractive multifocal intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p> <p>3.Truevue (Hydrophobic Foldable Arcylic Toric multifocal intraocular Lens)- Hydrophobix foldable arylc – Sterile -5- E &D</p> <p>4. Novaflex (Hydrophobic Foldable Arcylic Toric multifocal intraocular Lens)- Hydrophobix foldable arylc – Sterile -5- E &D</p> <p>5. Excellena (intraocular lens)- Poly methyl Methacrylate (PMMA)- Sterile -5 – E & D</p> <p>6.Solars (intraocular Lens)- Poly methyl Methacrylate (PMMA)- Sterile -5 – E & D</p> <p>7.Truevue (Hydrophobic Foldable Arcylic intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p> <p>8. Novaflex (Hydrophobic Foldable Arcylic intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p> <p>9. Truevue (Hydrophobic Foldable Arcylic Multifocal intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p> <p>10. Novaflex (Hydrophobic Foldable Arcylic Multifocal intraocular Lens)- Hydrophobic foldable arylc – Sterile -5- E &D</p>	CLAA/MD/TN/0 3/2010-DC	<p>Mfg Licno No.TN0002843 Validity from 01.10.2010 to 30.09.2020</p> <p>(05-05-2017) (Approval for additional Product)</p>
57.	M/s Dalvi and Sons, G- 10,Additional MIDC, Kudal,Dist.Sindhudurg,Mah arashtra	1.Locking Plate and Screw- 316L (Steel)- Non Sterile – Domestic	CLAA/MD/MS/0 1/2015-DC	<p>Mfg Licno No.MH/10512 Validity from 23-10-2015</p> <p>(05-05-2017) (Approval for additional Product)</p>

58.	M/s B.D.Surgical Industries, Plot No.2082,MIE, Part-II,Bahadurgarh,Jhajjar,Haryana-124507	<p>1. Mini Cortex Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>2. Mini Cortex Screw (Self Tapping)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>3. Cortex Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>4. Cortex Screw (Self tapping)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>5. Cancellous screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>6. Cannulated cancellous screw (self tappng/self drilling)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>7. Maxillofacial Screws (single slotted/cross headed slotted)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>8. Lag Screw (Maxillofacial single slotted/cross head slotted)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>9. Cannulated herber screw (Self tapping/self drilling)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>10. Malleolar screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>11. Shaft Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>12. Hip Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>13. Compression Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>14. DHS/DCS Compression Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>15. End Caps for interlocking nails - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>16. Locking Screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>Plug Screw for interlocking nail - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>17. Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>18. End cap for screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>19. Locking screws (Self tapping/self drilling)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>20. Plug Screw for screws - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>21. Set Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p>	CLAA/MD/HR/05/2016-DC	<p>Mfg Licno No.</p> <p>(05-05-2017)</p> <p>(Approval for additional Product)</p>
-----	---	--	-----------------------	---

22.Locking Bolt- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

23.Cancellous traction screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

24.ACL Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

25.Anchor Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

26.Orthopaedic Screws,Bone- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

Cannulated conical screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

27.Lag Screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

28.Locking Spacer- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

29.Pediatric Lag screw- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

30.Small washers- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

31.Spike washers- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

32.Button- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

33.Suture Washer- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

34.Self Tapping screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

35.Self tapping emergency screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

36.Self Drilling screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

37.BMR Screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

38.IMF Screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

39.Orthodontic Screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

40.T Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

41.DCP Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

42. L Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

43.H Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

44.Reconstruction plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

45.Straight Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

46.Dorsal Plates(Straight/Right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

47.Volar Plates (Straight/left/right) - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

48.Volar Locking Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

49.Wrist Plates (Straight/Right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

50.Wrist Locking plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

51.Ulna Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

52.Ulna Locking Plates(Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

53.Fibular Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

54.Fibular locking Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

55.Elbow Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

56.Elbow Locking Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

57.Humerus Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

58.Humerus locking Plates -(Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

59.Pelvis Plate ((Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

60.Pelvis Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

61.Femoral Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

62.Femoral Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

63.Tibial Plates (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

64.Tibial Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

65.Distal Radius Dorsal Plate Locking (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

66.Distal Radius Volar Plate Locking (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

67.Distal Humerus Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

68.Medial Distal tibia Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

69.Distal Humerus Medial Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

70.Small Dynamic Compression locking plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

71.Dynamic Compression plate iwth LC Under Cuts- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

72.Orthopaedic Plate Bone - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

73.Small T Locking Plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

74.T locking Plate (straight)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

75.Proximal tibia locking plate (Straight/right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

76.Anterolateral distal tibia locking plate (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

77.Cloverleaf locking plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

78.Metaphyseal locking plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

79.Metaphyseal Locking Plate Distal Medial Humerus (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

80.Metaphyseal locking Plate Distal Medial tibia (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

81.Pilion Lokcing plate (Straight/Right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

82.Calcaneal Plate (Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

83.Narrow Locking plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

84. Broad locking plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

85. Broadly curved locking plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

86. Clavicle locking plate system (straight/right/left/curved)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

87. Clavicle Plates (STRAIGHT/RIGHT/LEFT/CURVED)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

88. Distal femoral locking plate (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

89. Locking Plate,Bone- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

90. Lateral Tibia LockinG plate (Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

91. Proximal Femoral locking plate (Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

92. Condylar Femoral Locking Plate (straight/Left/Right)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

93. Trochanter locking plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

94. Trochanter Plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

95. One Third Tubular palte- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

96. Small T Plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

97. Lecestre plate (Straight/Right/Left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

98. Anatomical Plates 9Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

99. Distal Fibular Plate (straight/Right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

100. Proximal Ulna Plate (Striaght/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

101. Y-reconstruction plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

102. Reconsrtuction plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

103. Semi Tubular Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

104. Y-Fork Plates(straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

105. Distal Lateral Humeral Plate (Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

106. Pheelos locking plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

107. dhs LOCKING Plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

108. DCS Locking Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

109. Broad Synamic Compression Plate - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

110. Lengthening plate narrow - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

111. Lengthening palte broad- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

112. DCS Plate (ALL Angles)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

113. DCS Plates (All Angles)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

114. Condylar Buttress plate (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

115. Lateral tibia plate (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

116. Angles blades plate (all angles and blade lengths)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

117. Roy Camille plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

118. Sherman Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

119. Proximal palte lateral side (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

120. Pediatric condylar plate (all nagles and blade lengths)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

121. Distal Humeral Plate (Straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

122. Cobra palte- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

123. Hook plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

124. Spoon plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

125. Broad dynamic compression plate with LC under cuts- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

126.Narrow dynamic compression plate with LC under cuts- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

127.Jewett nail plate (all angles and blade lengths)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

128.L Plate (straight/right/left)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

129.T Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

130. T Buttress Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

131.Pediatric DHS Plate (all angles)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

132.DHS Plates (all angles)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

133.Condylar balde plate 9all angles and blade lengths)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

134.Angles blade palte for interochanteric femoral osteotomies)- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

135.Col De cygnet plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

136.Pediatric DCS Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

137.Orbital plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

138.Continuous iiole plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

139.X Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

140.Z Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

141.Double Y Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

142.Mesh Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

143.Square Frame plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

144.Rectangular fram plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

145.TMJ Condyle plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

146.TMJ Condyle plate with reconstruction- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

147.Chin plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

148.Angular plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

149.Cranioo maxillofacial plates- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

150.Splitfix plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

151.Curved plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

152.Strut Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

153.Mandible plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

154.Mandible reconstruction plates- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

155.Double I Plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

156.Straight plate with bridge- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

157.Straight plate SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

158.Straight plate with Gap - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

159.Orbital plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

160. "L" shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

161."L" Shape plate with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

162."I" SHAPE PLATE 110- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

163."T"shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

164."T" shape plate with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

165.Double "T" shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

166.Double "T"shape plate with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

167."H" Shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

168."Y" Shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

169."Y" shape plate with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

170. "X" shape plate- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

171."X" shape plate iwth gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

172."Z" shape palte with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

173.CAD plate with triangular- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

174.CAD plate square- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

175.CAD plate rectangular- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

176.CAD Plate rectangular strengthened - SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

177.Mesh plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

178.Lefort plate with gap 110- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

179.Chin plate standard- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

180.Chin plate bull- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

181.Mandibular plate- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

182.Mandibular plate curved- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

183.Mandibular plate straight- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

184.BMR plate straight- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

185.BMR Plate with gap- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

186.BMR plate angled- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

187.BURR hole cover plate- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D

188.Neuro gap plate- SS316LVM, SS316L, Titanium,Chormium Cobalt-Non Sterile – 5- E &D

189. Interlocking femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

190. Interlocking tibial nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

191. Femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

192. Tibial nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

193. Multi fix femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

194. Cannulated tibial nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

195. Unreamed tibial nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

196. Unreamed femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

197. Unreamed multifix femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

198. Proximal femoral nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

199. Orthopaedic nails n bones- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

200. Gamma nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

201. Intramedullary supra condylar nails (right/left)- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

202. Recon Nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

203. Unreamed Humeral nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

204. Cannulated humeral nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

205. Elastic nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

206. Bone, nails, orthopaedic- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

207. Kuntscher nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

208. Rush Nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

209. Enders Nail- SS316LVM, SS316L, Titanium, Chormium Cobalt- Non Sterile – 5- E & D

		<p>210.Square Nail- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>211.Bend Kuntscher nail- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>212.Interlocking nail- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>213.V nails- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>214.Knowel pins- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>215.Steinmann pin- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>216.Schanz Screws- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>217.Centrally threaded schanz pins- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>218.Kirschner wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>219.Guide wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>220.Drumpers wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>221.Wire,orthopaedics- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>222.Wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>223.Cercalage wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p> <p>224.Threaded Krischner wire- SS316LVM, SS316L, Titanium,Chormium Cobalt- Non Sterile – 5- E &D</p>		
59.	M/s Meril Healthcare Pvt.ltd., Survey No.135/2/B & 174/2,1 st floor, H1- H3,meri lpark, muktanand	<p>1.Total Knee System-Destiknee- CO-Cr-Mo Alloy-ETO Gas- 5- E & D</p> <p>2.Total Knee System-Freedom Altimo-Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E & D</p> <p>3.Total Knee System-Freedom Avon – Co-Cr-Mo Alloy/HXLPE- ETO Gas-</p>	CLAA/MD/GUJ/05/2013-DC	Mfg Licno No.G/28/1465 Validity from 01-02-2014 to 31-01-2019

	marg, chala,vapi-396191	<p>5- E& D</p> <p>4. Total Knee System-Freedom Extend – Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E& D</p> <p>5.Total Knee System-Freedom Hi-Flex – Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E& D</p> <p>6. Total Knee System-Freedom 360 – Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E& D</p> <p>7. Total Knee System-Freedom Prime – Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E& D</p> <p>8. Total Knee System-E-CIMA – Co-Cr-Mo Alloy/HXLPE- ETO Gas- 5- E& D</p> <p>9. Total Knee System-Freedom XP – Co-Cr-Mo Ally/HXLPE- ETO Gas- 5- E& D</p> <p>10. Total Knee System-Freedom Excel – Co-Cr-Mo Ally/HXLPE- ETO Gas- 5- E& D</p> <p>11. Total Knee System-Freedom i2 – Co-Cr-Mo Ally/HXLPE- ETO Gas- 5- E& D</p>		(30-05-2017) (Approval for additional Product)
60.	M/s Disposafe Health and Life Care Ltd., Plot No.1 & 2, Sector-59,Phase-II,Ballabgarh-121004,Faridabad (Haryana)	<p>1.Intravenous/Cannula Catheter (14-26G)- 3 Med-Polypropylene Polyethylene Polyacetal Stainless Steel (AISI304)- Sterile – 5 – Export</p> <p>2. Intravenous/Cannula (14-26G)- Nervetec- Polypropylene Polyethylene Polyacetal Stainless Steel (AISI304)- Sterile – 5- E & D</p> <p>3. Intravenous/Cannula (14-26G)- Winner- Polypropylene Polyethylene Polyacetal Stainless Steel (AISI304)- Sterile – 5- E & D</p>	CLAA/MD/HR/01/2013-DC	<p>Mfg Licno No.714-B (H) Validity from 27.05.2013 to 26.05.2018</p> <p>(30-05-2017) (Approval for additional Product)</p>
61.	M/s Sutures India Pvt.ltd., No.472-D, 13 th cross, 4 th phase, Peenya Indli.Area, Banaglore-560058	Alpha Lene Mesh – Polypropylene- Sterile by EtO- 5-Export	CLAA/MD/KAR/02/2012-DC	<p>Mfg Licno No.KTK/28/401/2012 Validity from 18.05.2012 to 17.05.2017</p> <p>(30-05-2017) (Approval for additional Product)</p>

62.	M/s Sysko Asia Healthcare Pvt.ITd., Plot No.176,Sector-59,Industrial Area, Huda,Faridabad-121004	<ol style="list-style-type: none"> 1.Anaesthesia-medical Grade PVC- Sterile – 4- E &D 2.Central Cenous Access Cathaters- Medical Grade PVC- Sterile – 4- E&D 3.Gastroenterology- Medical Grade PVC- Sterile – 4- E&D 4.Urology- Sysko Foley Ballon Catheter – Silicon- Sterile – 4 – E&D 5. Urology- Sysko TUR Set– Medical Grade PVC- Sterile – 4- E&D 6.Surgery and Wound Dranage- Medical Grade PVC- Sterile – 4- E&D 7. Surgery and Wound Dranage – Silicon Manufal Resuscitators-Silicon- Sterile – 4- E&D 	CLAA/MD/HR/04/2016-DC	<p>Mfg Licno No.714-B (H) Validity from 27.05.2013 to 26.05.2018</p> <p>(06-06-2017) (Approval for additional Product)</p>
63.	M/s Meril Life Sciences Pvt.ltdl., Survery No.135/2a & 135/2/B,174/2 (M1-M2 Meril Park,4 th floor), Near G.M.Bilakhia Stadium, Muktanand Marg,Chala, Vapi,Valsad-369191,Gujarat	<ol style="list-style-type: none"> 1.Atrial Septal Defect Occluder (LotusTM Atrial Septal Defect (ASD) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 2. Atrial Septal Defect Occluder (FloretTM Atrial Septal Defect (ASD) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 3. Ventricular Septal defect (VSD) Occulder (Ventricular septal Defect (VSD) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 4. Muscular Ventricular Septal defect (VSD) Occulder (FloretTM Muscular Ventricular septal Defect (VSD) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 5. Patent Ductus Arteriosus (PDA)Occulder (LtousTM Patent Ductus Arteriosus (PDA) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 6. Patent Ductus Arteriosus (PDA)Occulder (FloretTM Patent Ductus Arteriosus (PDA) Occluder) – Nitinol with Polyethylene terephthalate Fabric – Sterile – 2 – E & D 	CLAA/MD/GUJ /01/2006-DC (Pt.IV)	<p>Mfg Licno No.G/28/1273 Validity from 20.06.2013 to 19.06.2018</p> <p>(13-06-2017) (Approval for additional Product)</p>
64.	M/s Sigma Surgicals Pvt.Ltd., Plot No.13,Shree Ram Industrial Estate, b/h CMC Sohni`'s Chawl,Cross Road, Odhav, Ahmedabad-382415,Gujarat	<ol style="list-style-type: none"> 1.Bone compression plate,small (2.7MM) – Titanium/SS316I/SS316LVM- Non sterile – E & D 2..Bone compression plate,small (3.5MM) – Titanium/SS316I/SS316LVM- Non sterile – E & D 3..Bone compression plate,small (4.5MM) – Titanium/SS316I/SS316LVM- Non sterile – E & D 4.Minimum contact bone compression plate,small (2.7mm) – Titanitum/ SS316I/SS316LVM- Non sterile – E & D 5. Minimum contact bone compression plate,small (3.5mm) – Titanitum/ SS316I/SS316LVM- Non sterile – E & D 6. Minimum contact bone compression plate,narrow –Titanitum/ SS316I/SS316LVM- Non sterile – E & D 7. Minimum contact bone compression plate,broad–Titanitum/ SS316I/SS316LVM- Non sterile – E & D 	CLAA/MD/GUJ /11/2010-DC	<p>Mfg Licno No.G/28/1352 Validity from 27.09.2015 to 26.09.2020</p> <p>(13-06-2017) (Approval for Renewal Product)</p>

- 8.Femur head compression plate (125),round hole,(38mmbareel)-Titanium/ SS316/SS316LVM- Non sterile – E & D
9. Distal Femoral compression plate,b.c.hole -Titanium/ SS316/SS316LVM- Non sterile – E & D
- 10.Condylar buttress plate(left n right)- Titanium/ SS316/SS316LVM- Non sterile – E & D
- 11.Lateral tibial plate(left n right)- Titanium/ SS316/SS316LVM- Non sterile – E & D
- 12.Semi Tubular plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 13.One third tubular plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 14."L"Buttress plate (left n righte)- Titanium/ SS316/SS316LVM- Non sterile – E & D
- 15."T"Simple plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 16."T"Buttress Plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 17."T"small plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
18. "T"Oblique plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 19.Cloverleaf plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 20.Reconstruction plate,round hole - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 21."Y"Reconstruction plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 22.Angle blade plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
23. Condylar plateTitanium/ SS316/SS316LVM- Non sterile – E & D
- 24.Clavicle hook plate (left and right)-Titanium/ SS316/SS316LVM- Non sterile – E & D
- 25.Hook plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 26.Spoon plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 27.Cobra Head Plate- Titanium/ SS316/SS316LVM- Non sterile – E & D
- 28.Ellis plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
- 29.Twin lock –BCP Clavicle hook plate (left and right)- Titanium/ SS316/SS316LVM- Non sterile – E & D
30. Twin lock –BCP Femur Head plate- Titanium/

SS316/SS316LVM- Non sterile – E & D
31.Twin Lock –BCP Distal femoral plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
32.Twin lock bone compression (small, broad and narrow) - Titanium/ SS316/SS316LVM- Non sterile – E & D
33.Twin lock BCP “L”Butress plate left and right - Titanium/ SS316/SS316LVM- Non sterile – E & D
34. Twin lock BCP “T” (Butress,simple,small, oblique and distal femoral ,lateral (left and right) plate- Titanium/ SS316/SS316LVM- Non sterile – E & D
35.Twin lock –BCP Lateral tibialplate (left and right) - Titanium/ SS316/SS316LVM- Non sterile – E & D
36. Twin lock –BCP Recon Plate(R/hole, B.C.Hole)- Titanium/ SS316/SS316LVM- Non sterile – E & D
37.Twin lock –BCP Proximal humrus plate-Titanium/ SS316/SS316LVM- Non sterile – E & D
38. Twin lock –BCP (Metaphyseal pate simple and distal medial tibia- Titanium/ SS316/SS316LVM- Non sterile – E & D
39. Twin lock –BCP Distal humrus plate Xyston-Left and right with lateral support - Titanium/ SS316/SS316LVM- Non sterile – E & D
40.Twin lock –BCP Distal Kopsis palte (left and right)- Titanium/ SS316/SS316LVM- Non sterile – E & D
41.Twin lock –BCP (Suraline distal kopsis plate left and right) - Titanium/ SS316/SS316LVM- Non sterile – E & D
42.Twin lock BCP Eight plate - Titanium/ SS316/SS316LVM- Non sterile – E & D
43.Cannulated canellous bone screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
44. canellous bone screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
45. Cortical bone screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
46.Interlocking Bone Screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
47.Twin lock and BCP Screw (self taping and self drilling)- Titanium/ SS316/SS316LVM- Non sterile – E & D
48.Fixator bone screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
49. Harp screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
50.Mono and Poly Axial Spinal Screw-Titanium/ SS316/SS316LVM- Non sterile – E & D
51.Wire Mount -Titanium/ SS316/SS316LVM- Non sterile – E & D

52.Spacer- Titanium/ SS316/SS316LVM- Non sterile – E & D
53.Washer- Titanium/ SS316/SS316LVM- Non sterile – E & D
54.Classic interlocking nail for femural bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
55.Recon interlocking nail for femural bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
56 Antrgrade interlocking nail for femural bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
57. Proximal interlocking nail for femural bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
58. Distal vinterlocking nail for femural bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
59. Classic interlocking nail for humral bone -Titanium/ SS316/SS316LVM- Non sterile – E & D
60.Enders Nail -Titanium/ SS316/SS316LVM- Non sterile – E & D
61.Rush Nail-Titanium/ SS316/SS316LVM- Non sterile – E & D
62.Squire Nail -Titanium/ SS316/SS316LVM- Non sterile – E & D
63.Kuntcher Nail -Titanium/ SS316/SS316LVM- Non sterile – E & D
64.K.Wire-Titanium/ SS316/SS316LVM- Non sterile – E & D
65.Bipolar Prosthesis-Titanium/ SS316/SS316LVM- Non sterile – E & D
66.Shoulder Prosthesis -Titanium/ SS316/SS316LVM- Non sterile – E & D
67.Tibial Base plate round and D Hole-Titanium/ SS316/SS316LVM- Non sterile – E & D
68.Vaya Cage- Titanium/ SS316/SS316LVM- Non sterile – E & D
69. .Vaya Expandable Cage- Titanium/ SS316/SS316LVM- Non sterile – E & D
70.Cagedisk - Titanium/ SS316/SS316LVM- Non sterile – E & D
71. Triton cage - Titanium/ SS316/SS316LVM- Non sterile – E & D
72.Spinal Road-Titanium/ SS316/SS316LVM- Non sterile – E & D
73.Inner Screw- Titanium/ SS316/SS316LVM- Non sterile – E & D
74.Spiked Washer- Titanium/ SS316/SS316LVM- Non sterile – E & D
75.Staple washer for pedicular screws-Titanium/ SS316/SS316LVM- Non sterile – E & D
76.Pedical Hook- Titanium/ SS316/SS316LVM- Non sterile – E & D
77.Lamina Hook (Transverse, offset)-Titanium/ SS316/SS316LVM- Non sterile – E & D
78.TPLO Plate (left and right) - Titanium/ SS316/SS316LVM- Non

		sterile – E & D 79.Spiked ligament staple for knee- Titanium/ SS316I/SS316LVM- Non sterile – E & D 80.Siwik Button- Titanium/ SS316I/SS316LVM- Non sterile – E & D 81.Puudu Plate - Titanium/ SS316I/SS316LVM- Non sterile – E & D		
65.	M/s Biotech Vision Care Pvt.Ltd., Plot No.555,556,557,Khatraj-Vadsar Road, At & Post Khatraj, Tal.Kalol,Distt.Gandhi Nagar	1.Hydrophobic Arcylic Foldable Intraocular Lens in Preloaded Delivery System (PAHFY600F) – Natural Yellow Hydrophobic Acrylic Foldabale Buttons – Sterile – Domestic 2. Hydrophobic Arcylic Foldable Intraocular Lens in Preloaded Delivery System (PLHF2) –Hydrophobic Acrylic Foldabale Buttons – Sterile – Domestic 3. Hydrophobic Arcylic Foldable Intraocular Lens in Preloaded Delivery System EYECRYL SERT (PLHF2) – Hydrophobic Acrylic Foldabale Buttons – Sterile – Export 4. Hydrophobic Arcylic Foldable Intraocular Lens in Preloaded Delivery System OPTIFLEX (MO/HF-011) –Hydrophobic Acrylic Foldabale Buttons – Sterile – Export 5. Hydrophobic Arcylic Foldable Intraocular Lens in Preloaded Delivery System OPTIFLEX (MO/HF-012) – Natural Yellow Hydrophobic Acrylic Foldabale Buttons – Sterile – Export	CLAA/MD/GUJ /03/2009-DC	Mfg Licno No.G/28/1298 Validity from (13-06-2017) (Approval for additional Product)
66.	M/s Harsoria Helathcare Pvt.Ltd.,-110-111,Phase-IV, Udyog Vihar,Gurgaon,Haryana	1.Long-Term Elite Harsoria Hemodialysis Catheter Kit Dual Lumen – Polyurethane Nitinol & SS PP SS Needle ABS Non Woven Fabric – Sterile – 3 – E & D	CLAA/MD/HR/ 04/2006-DC (Pt-2)	Mfg Licno No.DML No.637-B (H) Validity from 19.09.2016 to 18.09.2021 (13-06-2017)

				(Approval for additional Product)
67.	M/s Lars Medicare Pvt.Ltd., Kila No. 16 & 17,Village Sultanpur, G.T.Road, Bahalgarh CHowk, Distt.Sonepat-131021,Haryana	1.Safety IV Cannula (Ubi-one) – Radi-opaque strip/clear Teflon,PP,POM, Silicon,HDPE,LDPE & Stainless Steel – Sterile – 5- Export	CLAA/MD/HR/ 01/2008-DC	Mfg Licno No.660-B (H) Validity from 04.02.2013 to 18.09.2018 (13-06-2017) (Approval for additional Product)
68.	M/s Newtech Medical Devices, 14/5,1 st floor,Mathura Road,Faridabad- 121003,Haryana	1.Dilation Catheter – Medical Pebax/Medical Nylon – Sterile – 3 – E &D 2.Guiding Catheter – Medical Pebax – Sterile – 3 – E & D 3.Temporary Pacing Lead – Medical Polyurethane – Sterile – 3 – E& D 4.Central Venous Catheter Kit – Catheter (medical Polyurethane TPU), Dilator (medical PE/PP),Guidewire (SS304),Blade (Ss 304), Syringe (Medical PP) – Sterile – 3 – E& D 5.Dialysis Catheter Kit - Catheter (medical Polyurethane TPU), Dilator (medical PE/PP),Guidewire (SS304),Blade (Ss 304), Syringe (Medical PP) – Sterile – 3 – E& D 6.PD Catheter – Medical Silicone - Sterile – 3 – E& D 7.Introducer Sheath – Medical PTFE - Sterile – 3 – E& D 8.Diagnostic Catheter – Medical PEBAX- Sterile – 3 – E& D 9.Aspiration Catheter - Medical PEBAX- Sterile – 3 – E& D 10.Vascular Snare – Medical Nitinol + Medical Pebax - Sterile – 3 – E& D 11. Angio/PTCA Kit – Manifold (Medical poly carbonate), Y Connector (medical poly carbonate), PM Line (Medical Polyvinyl Chloride) HPT (Medical Polyvinyl Chloride),syringe (medical polypropylene),control syringe (medical polycarbonate) - Sterile – 3 – E& D 12.Guidewire – SS 304- Sterile – 3 – E& D	CLAA/MD/HAR /01/2017	Mfg Licno No. Validity from (29-06-2017) (Approval for additional Product)
69.	M/s Sahajanand Medical Technologies Pvt.Ltd., 2 nd & 3 rd floor, Plot No.33,34,35,52,53 & 54 ,Surat Special Economic Zone,Sachin,Surat- 394230,Gujarat	1.The Millennium Coronary Stent System – Stent system SS316 LVM Metal mounted on a PTCA balloon catheter – Sterile – 36 months – E &D 2.Sirolimus Eluting Coronary Stent System – Stent –Cobalt chromium metal alloy with blended coating of biodegradable polymers and sirolimus drug (Active ingredient),mounted on a PTCA balloon catheter – Sterile – 24 months – E &D 3.Percutaneous Transluminal angioplasty catheter – The catheter is made up to biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – Sterile –	CLAA/MD/GUJ /06/2011-DC	Mfg Licno No.G/28/1385 Validity from 29/08/2016 to 28/08/2021 (29-06-2017) (Approval for Renewal Product)

36months – E & D

4.Aspiration Catheter – The catheter is made up to biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – Sterile – 36months – E & D

5.Cobalt chromium coronary stent system – Stent made up to cobalt chromium alloy mounted on a balloon catheter – Sterile – 36 months – E &D

6.Cobalt Chromium coronary stent system – Stent made up of cobalt chromium alloy mounted on a balloon catheter – Sterile – 36months – E & D

7.Balloon Expandable cobalt –chromnum reneal and biliary stent system – stent made up tp cobalt chromnium alloy mounted on a balloon catheter- Sterile – 36months – E & D

8.Cobalt chromioum sirolimus eluting coronary stent system – Stent- cobalt chromium metal alloy with blended coating of biodegradeable polymers and sirolimus drug (Actuve ingredient),mounted on a balloon catheter – Sterile – 30months – E & D

9.Sirolimus eluting cobalt chromium coronary stent system – Stent-cobalt chromium metal alloy wth blended coating of biodegradeable polymers and sirolimus drug (active ingredient), mounted on a balloon catheter – Sterile – 30months – E & D

10. Sirolimus eluting cobalt chromium coronary stent system – Stent-cobalt chromium metal alloy with blended coating of biodegradeable polymers and sirolimus drug (active ingredient), mounted on a PTCA balloon catheter – Sterile – 3 months – E & D

11.Everolimus Eluting Cobalt chromium coronary stent system (EES-1) – Stent- cobalt chromium metal alloy with blended coating of biodegradeable polymers and everolimus drug (active ingredient), mounted on a balloon catheter – Sterile –24 months – E & D

12. Everolimus Eluting Cobalt chromium coronary stent system (EES-2) – Stent- cobalt chromium metal alloy wth blended coating of biodegradeable polymers and everolimus drug (active ingredient), mounted on a balloon catheter – Sterile –12 months – E & D

13.Cobalt Chromium coronary stent system on tetrinium platform – Stent-Cobalt chromium metal mounted on ballooon catheter – Sterile – 36months – E& D

14.Everolimus eluting cobalt chromium coronary stent system on tetrinium platform – Stent- Cobalt chromium metal alloy with blended coating of biodegradeable polymers and everolimus drug (active ingredient),mounted on a PTCA balloon catheter – Sterile – 24 months – E& D

15. Sirolimus eluting cobalt chromium coronary stent system on tetrinium platform – Stent-cobalt chromium metal alloy with blended coating of biodegradable polymers and sirolimus drug (active ingredient), mounted on a PTCA balloon catheter – Sterile – 30months – E & D
- 16.PTCA Balloon Catheter – The balloon catheter is made of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – NonSterile – 36 months – E & D
17. PTCA Balloon Dilatation Catheter (SC-Semi Compliant) – The balloon catheter is made of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – NonSterile – 36 months – E & D
18. PTCA Balloon Dilatation Catheter (NC-Non Compliant) – The balloon catheter is made of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – NonSterile – 36 months – E & D
- 19.PTCA Balloon Dilatation Catheter (CTO-Chronic Total Occlusion) – The balloon catheter of the stenotic portion of a coronary artery or bypass fraft stenosis – Sterile – 36 months – E & D
- 20.PTCA OTW Balloon Catheter – The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Sterile – 3 years – E & D
- 21.PTA Balloon dilatation cathter- The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Sterile – 36 months– E & D
- 22.PTA Balloon catheter- The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Sterile – 36 months– E & D
- 23.PTA Balloon dilatation catheter (NC-Non compliant) - The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Sterile – 36 months – E & D
24. PTA Balloon dilatation cathter- The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Sterile – 3 years – E & D
25. PTA Balloon dilatation cathter- The balloon catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shraft – Non- Sterile – 3 years – E & D
- 26.Percutaneous transluminal angioplasty catheter (RX)-Peripheral

		<p>– The catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – Sterile – 36 months – E & D</p> <p>27. Percutaneous transluminal angioplasty catheter (OTW)- Peripheral – The catheter is made up of biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft – Sterile – 36 months – E & D</p> <p>28. Guide Catheter – Extruded plastic components – Sterile – 42 months – E & D</p> <p>29. Micro Catheter - Extruded plastic tubing – Sterile – 36 months – E & D</p> <p>30. Introducer Sheath – Plastic components – Sterile - 36 months – E & D</p> <p>31. Guidewire – Metal core wire and plastic components – Sterile - 36 months – E & D</p> <p>32. Manifold – Plastic components - Sterile - 36 months – E & D</p> <p>33. Push Pull Y Connector - Plastic components - Sterile - 36 months – E & D</p> <p>34. Guide Catheter – Extruded plastic components – Sterile – 42 months – E & D</p> <p>35. Everolimus eluting cobalt chromium coronary stent system – Stent- Cobalt chromium metal alloy with blended coating of biodegradable polymers and everolimus drug (active ingredient), mounted on a PTCA balloon catheter – Sterile – 24 months – E & D</p> <p>36. Cobalt Chromium coronary stent system – Stent-Cobalt chromium metal mounted on balloon catheter – Sterile – 36 months – E & D</p> <p>37. Sirolimus eluting cobalt chromium coronary stent system - Stent-cobalt chromium metal alloy with blended coating of biodegradable polymers and sirolimus drug (active ingredient), mounted on a PTCA balloon catheter – Sterile – 30 months – E & D</p> <p>38.</p>		
70.	M/s Disposafe Health and life care Ltd., Plot No.1 & 2, Sector-	1. Intravenous Cannula – Polypropylene polyethylene polyacetal stainless steel (AISI304) – Sterile – 5 – E & D	CLAA/MD/HR/01/2013-DC	Mfg Lic No.714-B (H) Validity from 27.05.2013 to

	59, Phase-II, Ballabgarh-121004, Faridabad, Haryana-121004			26.05.2018 (11-07-2017) (Approval for Additional Product)
71.	M/s Vishal Surgitech Pvt.Ltd., Nr.Old Tall Plaza, Plot No.10, N.H.8B, Kothariya, Gondal Road, Rajkot	1.Cortex Screw – SS 316 L- Non Sterile – 5 – E& D 2.Cancellous Screw - SS 316 L- Non Sterile – 5 – E& D 3.Cannulated Cancellous Screw - SS 316 L- Non Sterile – 5 – E& D	CLAA/MD/GUJ/02/2017-DC	Mfg Lic No.G/28/1604 Validity from 18/03/2017 to 17/03/2022 (11-07-2017) (Approval for Additional Product)
72.	M/s Bombay Ortho Industries, Plot No.G-551, Lodhika G.I.D.C, Matoda, Dist.Rajkot-360021 (Gujarat)	1.Bone Screw – Stainless Steel/Titanium – Non Sterile – E & D 2. Bone Plates - Stainless Steel/Titanium – Non Sterile – E & D 3. Interloc King Nail - Stainless Steel/Titanium – Non Sterile – E & D 4. Pins & Wires System - Stainless Steel/Titanium – Non Sterile – E & D 5. Joint Replacement Prosthesis - Stainless Steel/Titanium – Non Sterile – E & D 6. Joint Replacement Prosthesis - SS/TI/UHMWPE/Cr-Co – Non Sterile – E & D 7.Spinal Implants - Stainless Steel/Titanium – Non Sterile – E & D	CLAA/MD/GUJ/03/2010	Mfg Lic No.G/28/11322 Validity from 04-01-2015 to 03.01.2020 (11-07-2017) (Approval for Additional Product)
73.	M/s Global Medikit Ltd., Khasra No.323 (MI), Camp Road, Selaqui, Dehradun, Uttarakhand	1.Latex Floeys Catheter (Avan) – Catheter, HARD Valve, Drainage, Funnel – Sterile – 5 – Export 2.IV Catheter with injection valve and wings – body, needle, catheter, silicon tube, port, cap, holder, threaded stopper, needle hub and cover - Sterile – 5 – Export 3.IV Catheter with wings for neonates – Body, needle, catheter, holder, threaded stopper, needle hub and cover - Sterile – 5 – Export 4.Silicon Foleys catheter – Catheter, hard valve, drainage, funnel - Sterile – 5 – Export 5.Suction Catheter – Tube and connector - Sterile – 5 – Export 6.Nelation catheter - Tube and connector - Sterile – 5 – Export 7. Safety IV Catheter with injection valve and wings – body, needle, N.Cover, Catheter, holder, threaded stopper, N Hub and Cover, Safety device - Sterile – 5 – Export 8.Tracheostomy cuffed – Tube, tube guide, tie ribbon, pilot balloon/valve assemble - Sterile – 5 – Export 9.Central venous catheter and IV Cannula – CVC Catheter, introducer cannula - Sterile – 5 – Export 10.Central venous catheter – single lumen – syringe, single lumen catheter, clamp, dilator, needle assembly, blade, three way stop	CLAA/MD/UA/01/2009-DC	Mfg DML No:-18/UA/SC/P-2009 Validity upto 23.03.2019 (11-07-2017) (Approval for Additional Product)

		cock,advancer, steeper - Sterile – 5 – Export 11. Central venous catheter – Double lumen – syringe, single lumen catheter, clamp, dilator, needle assembly, blade, three way stop cock,advancer, steeper - Sterile – 5 – Export 12. Central venous catheter – triple lumen – syringe, single lumen catheter, clamp, dilator, needle assembly, blade, three way stop cock,advancer, steeper - Sterile – 5 – Export		
74.	M/s Meril Life Sciences Pvt.ltdl., Survery No.135/2a & 135/2/B,174/2 (M1-M2 Meril Park,4 th floor), Near G.M.Bilakhia Stadium, Muktanand Marg,Chala, Vapi,Valsad-369191,Gujarat	1.Sirolimus Eluting Rx PTCA Balloon Dilatation Catheter – Sterile – 2 – E & D	CLAA/MD/GUJ /01/2008-DC (Pt.IV)	Mfg Lic No.G/28/1273 Validity from 20.06.2013 to 19.06.2018 (11-07-2017) (Approval for Additional Product)
75.	M/s Becton Dickinson India Pvt.Ltd., at Plot No.#1, Sector-3,IMT, Bawal,Distt.rewari-123501	1.IV Cannula (Venflonm) – Polypropylene, Methyl vinyl polysiloxane, Polytetrafluoroethylene Polycarbonate white,Polydimethylsiloxane – Sterile – 5 – E & D	CLAA/MD/HR/01/2009-DC	Mfg Lic No.671-B (H) Validity from 24.03.2014 to 23.03.2019 (25-07-2017) (Approval for Additional Product)
76.	M/s Kosdrug Pvt.Ltd., 105,MIE, Bahadugarh,Haryana	1.Aortic perfusion cannula – PVC Medial grade material- Sterile – 5 – Domestic 2.Venous perfusion cannula- PVC Medical grade material- Sterile – 5 – Domestic 3.Endotrachal tube – PVC Medical grade material, universal connector – Sterile – 5 - Domestic	CLAA/MD/HR/02/2013-DC	Mfg Lic No.742B (H) Validity from (25-07-2017) (Approval for Additional Product)
77.	M/s Lars Medicare Pvt.Ltd., Kila No.16 &17, Village sultanpur ,G.T.Road, bahalgarh chowk, sonapat-131021	1.IV Cannula (Fixo Flon) – radio opaque strip/clear teflon, PP,POM, Silicon, HDPE, LDPE and stainless steel – Sterile -5 – E & D	CLAA/MD/HR/01/2008-DC	Mfg Lic No.660-B (H) Validity from 04.02.2013 to 03.02.2018 (25-07-2017) (Approval for Additional Product)
78.	M/s Disposafe Health and life care ltd.,Plot No.1 & 2, Sector-	1.Sterile Scalp vein set – Polypropylene,Polyethylene, Polyacetal, Stainless Steel (AISI304),Poly Vinyl Chloride – Sterile – 5- E & D	CLAA/MD/HR/01/2013-DC	Mfg Lic No714-B (H) Validity from 27.05.2013 to

	59, Phase-II, Ballabgarh-121004, Faridabad, Haryana-121004	<p>2. Sterile Central Venous Catheter – Polyurethane, polycarbonate, Polypropylene – Sterile – 5 – E & D</p> <p>3. Sterile Foley Balloon Catheter- Silicon, Latex, Rubber, Polypropylene - Sterile – 5 – E & D</p> <p>4. Sterile Suction Catheter- Poly Vinyl Chloride polypropylene - Sterile – 5 – E & D</p> <p>5. Sterile Closed Suction Set - Polypropylene, Nylon, Polycarbonate, Polyoxymethylene, Poly Vinyl Chloride, Stainless Steel, Low density Polyethylene – Sterile – 5- E & D</p> <p>6. Sterile Oxygen Catheter- Poly Vinyl Chloride, Polypropylene - Polypropylene– Sterile – 5- E & D</p> <p>7. Sterile Umbilical Catheter – Polyurethane, Polypropylene- Sterile – 5 – E & D</p> <p>8. Sterile haemodialysis Catheter - – Polyurethane, Polypropylene- Sterile – 5 – E & D</p> <p>9. Sterile Peritoneal Dialysis Catheter Set – Polyamide, Stainless Steel, Polypropylene, ABS - Sterile – 5 – E & D</p> <p>10. Sterile Thoracic Catheter Set – Polyurethane, Polyvinylchloride, Polycarbonate, Stainless Steel - Sterile – 5 – E & D</p>		<p>26.05.2018</p> <p>(25-07-2017) (Approval for Additional Product)</p>
79.	M/s. Mediplus (India) Limited, 1261-62, M.I.E, Bahadurgarh, (Haryana).	<p>1. IV Cannula with injection valve with wings – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>2. IV Cannula without injection valve with wings – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>3. IV Cannula with injection valve – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>4. IV Cannula without injection valve and without wings – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>5. IV Cannula with port – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>6. IV Catheter – Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>7. CANNULA WITH INJECTION PORT WITH WINGS - Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>8. Cannula without injection valve with small wings - Polypropylene, polyethylene, Polyacetal, Stainless Steel (AISI304) – Sterile – 5 – E & D</p> <p>9. Safety IV Cannula with injection port and wings -</p>	CLAA/MD/HR/06/2006-DC	<p>Mfg Lic No:-646-B (H) Validity to 01.06.2017 to 31.05.2022</p> <p>(27-07-2017) (Approval for Renewal Product)</p>

		Polypropylene,polyethylene, Polyacetal, Stainless Steel(AISI304) – Sterile – 5 – E& D 10. Safety IV Cannula without injection port and wings - Polypropylene,polyethylene, Polyacetal, Stainless Steel(AISI304) – Sterile – 5 – E& D		
--	--	---	--	--

S. No.	Name of the Manufacturer	Product Licensed by CLAA	File No.	Manufacturing Licence No. & Validity
80.	M/s Bombay Ortho Industries ,Kishan Gate, Plot No.G-551,Lodhika,GIDC Metoda, Rajkot-360021,Gujarat,India	1. Arthroscopy Implants (Bombay Ortho) – Stainless Steel/Titanium + Non-Sterile – E & D	CLAA/MD/GU J/13/2011-DC	Mfg Licno No.G/28/1322 Validity from 04-01-2015 to 03-01-2020 (17-08-2017) (Approval for additional Product)
81.	M/s Disposafe Health and life care ltd.,Plot No.1 & 2, Sector-59, Phase-II,Ballabgarh-121004,Faridabad,Haryana-121004	1. Sterile Intravenous Cannula (Injexcath and Injexcath Safety) – Polypropylene, Polyethylene, Polyacetal, Stainless Steel (AISI304) - Sterile – 5- E & D	CLAA/MD/HR /01/2013-DC	Mfg Licno No.714-B (H) Validity from 27.05.2013 to 26.05.2018 (17-08-2017) (Approval for additional Product)
82.	M/s Envision Scitenific Pvt.Ltd., Plot No.C1B-108,Ichchhapore GIDC, Industrial Estate, Bhatpore, Sachin-34510,Gujarat,India	1.Sirolimus eluting coronary stent system (Abluminus)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D 2.Mitigator- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D 3.Sirolimus eluting coronary stent system (Nostrum)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D 4.Sirolimus eluting coronary stent system (AmazonuaSir)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO	CLAA/MD/GU J/13/2011-DC	Mfg Licno No.G/28/1396 Validity from 19/12/2011 to 18/12/2021 (22-08-2017) (Approval for Renewal Product)

		<p>Sterile – 2 – E& D</p> <p>5.Sirolimus eluting coronary stent system (ihtDEStinyBD)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>6. Sirolimus eluting coronary stent system (Agile) - Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – E& D</p> <p>7. Sirolimus eluting coronary stent system (Minoris)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>8. Sirolimus eluting coronary stent system (Serpentis)– Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>9.Sirolimus Coated Balloon Catheter (Magic Touch)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>10. Sirolimus Coated Balloon Catheter (Devoir)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>11. Sirolimus Coated Balloon Catheter (Vanquisher)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>12.Sirolimus Based Nano carrier eluting stent system (Focus np)- Sirolimus Coated Balloon Catheter (Magic Touch)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>13. Sirolimus Based Nano carrier eluting stent system (AmazoniaLex)-</p>		
--	--	---	--	--

		<p>Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>14. Sirolimus Based Nano carrier eluting stent system (Cynosure)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>15.Paclitaxel based Nano carrier eluting balloon catheter (Nano-pac)- Stent, Balloon,Hub,Strain relief,Proximal Tube, Inner distal tube, intermediate tube,marker brand, soft tip, sleeve, coating on stent and balloon,ETO Sterile – 2 – E& D</p> <p>16.Sirolimus Based Bifurcation stent system (Nile SIR) - Stent, Balloon,Hub,Strain relief,Proximal Tube, AUTO RELEASE SHEATH, Inner distal tube,marker brand, soft tip, stent and balloon,ETO Sterile – 2 – E& D</p> <p>17.Cobalt Chromium coronary stent system (Amazer)-Cobalt Chromium L605) – ETO Sterile -2.5 – E & D</p> <p>18.PTCA Dilatation Catheter (Mastery) –Polyamide Balloon- ETO Sterile – 2.5 – E &D</p> <p>19.PTCA Dilatation Catheter (Swimmer) – Polyamide Balloon - ETO Sterile – 2.5 – E &D</p> <p>20.PTCA Dilataion Cathter (Non-Compliant) (Swimmer NC) - Polyamide Balloon - ETO Sterile – 2.5 – E &D</p> <p>21.PTCA Dilatation Catheter (Chronic total occlusion)(Swimmer CTO) - Polyamide Balloon - ETO Sterile – 2.5 – E &D</p> <p>22.Aspiration catheter(TranspotAC) –Catheter, Stylet,ExtensionLine, Cell Strainer, Syringe-ETO Sterile -2.5 – E& D</p> <p>23.Sirolimus Coated PTA Balloon Catheter (OTW/RX) (Xtreme touch neo) – PTA Balloon,Polyether block amide ring maker, distal tubing, soft tip tubing, proximal shaft, pebax tubing,polycarbonate balloon coating – ETO Sterile – 2 – E & D</p>		
--	--	--	--	--

		&d 9.IV Cannula – PTFE/FEPPUR/PE,PP,SS,PA - Sterile for single use only – 5- E &d		
85.	M/s Sahajanand Medical Technologies Pvt.Ltd., 2 nd & 3 rd floor, Plot NO.33,34,35,52,53 & 54 ,Surat Speical Economic Zone, Sachin-394230,Surat,Gujarat	1.PTCA balloon dialation catheter (SC-Semi Complaint) – The balloon catheter is made up to biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft.The polymeric distal shaft is coated with biocompatible hydrophilic coating – Sterile -2 E & D 2. PTCA balloon dialation catheter (NC-Non Compliant) – The balloon catheter is made up to biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft.The polymeric distal shaft is coated with biocompatible hydrophilic coating – Sterile -2 E & D 3. PTCA balloon catheter– The balloon catheter is made up to biocompatible raw materials which mainly include extruded tubing made from polymers and metal shaft.The polymeric distal shaft is coated with biocompatible hydrophilic coating – Sterile -2 E & D	CLAA/MD/GU J/06/2011-DC	Mfg Licno No.G/28/1385 Validity upto 28-02-2021 (29-08-2017) (Approval for additional Product)
86.	M/s Newtech Medical Devices, 14/5,1 st floor Mathura Road, Faridabad-121003,Haryana	1.Temporary Pacing Catheter- Polyurethane – Sterile – 3 – E &D 2.Central Venous Catheter Kit – Catheter (medical Polyurethane,Dilator Medical,Guidewire Blade, (SS304),Needle – Sterile – 3- E&D 3. Dialysis Catheter Kit - Catheter (medical Polyurethane,Dilator Medical(PE/PP),Guidewire Blade, (SS304),Needle,Syringe – Sterile – 3- E&D 4.Daignostic Catheter –Medical PEBAX- Sterile – 3- E & D	CLAA/MD/HA R/01/2017-DC	Mfg Licno No.786-B (H) Validity from 11-07-2017 to 10-07-2022 (29-08-2017) (Approval for additional Product)
87.	M/s Biotech Vision Care Pvt.Ltd., Plot No.555,556,557,Khatraj-Vadsar Road, At & Post Khatraj, Tal:-Kalol,Distt-	1.Hydrophobic Acrylic foldable intraocular lens with Delivery System – (ASHFY6002)-Natural Yellow/Clear Hydrophobic acrylic foldable buttons-Sterile –Domestic 2. Hydrophilic Acrylic foldable intraocular lens with Hydrophobic Surface and Delivery System –(HSAS600)- Hydrophilic acrylic foldable	CLAA/MD/GU J/03/2009-DC	Mfg Licno No.G/28/1298 Validity from

	Gandhi Nagar	<p>buttons-Sterile –Domestic</p> <p>3. Hydrophilic Acrylic foldable intraocular lens with Delivery System – (S600MZ)- Hydrophilic acrylic foldable buttons-Sterile –Domestic</p> <p>4.Hydrophilic Acrylic Foldable Intraocular lens with delivery system (TP613)- Hydrophilic acrylic foldable buttons- Sterile –Domestic</p> <p>5.Hydrophobic Acrylic Foldable Intraocular lens with Delivery System (HFY-05) – Natural yellow hydrophilic acrylic foldable buttons-Sterile – Domestic</p> <p>6. Hydrophobic Acrylic Foldable Intraocular lens with Delivery System (HFYD-05) – Natural yellow hydrophilic acrylic foldable buttons-Sterile –Domestic</p> <p>7.Phakic Aspheric Hydrophobic Acrylic foldable intraocular lens with delivery system (PKC135NH)</p>		(29-08-2017) (Approval for additional Product)
88.	M/s Harsoria Healthcare Pvt.Ltd., 110-111,Phase-IV, Udyog Vihar, Gurgaon-122015,Haryana	<p>1.IV Cannula with catheter (Shaheen Flon) – Polypropylene HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 Silicone Tube)- Sterile -5- E& D</p> <p>2. IV Cannula with catheter (Medsorce) – Polypropylene HDPE,ABS,Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 Silicone Tube)- Sterile -5- E& D</p> <p>3.Harsoria (Central Venous Catheter Kit Quad Lumen) – Polyurethane nithol and SS,PP,SS Needle, ABS, Non woven Fabric – Sterile – 3- E & D</p> <p>4. Cadon (Central Venous Catheter Kit Quad Lumen) – Polyurethane nithol and SS,PP,SS Needle, ABS, Non woven Fabric – Sterile – 3- E & D</p>	CLAA/MD/HR /04/2006-DC (Pt-2)	Mfg Licno No.637-B (H) Validity from 19-09-2016 to 18-09-2021 (29-08-2017) (Approval for additional Product)
89.	M/s Prymas Healthcare LLP, Plot No.53/17,Industrial Area,NIT Faridabad-121001,Haryana	<p>1.Hemodialysis Catheter-Sterile-3- E& D</p> <p>2.Central Venous Catheter- Sterile-- E& D</p>	CLAA/MD/HR /01/2015-DC	Mfg Licno No.754-B (H),757-B (H) Validity from

				(29-08-2017) (Approval for additional Product)
92.	M/s Poly bond India Pvt.Ltd., Gate No.1088,Sanawadi,Tal:- Shirur,Pune-412210,Maharashtra	<p>1.IV Cannula with wings and port/IV Cannulae- Polypropylene,PTFE/FEP/PUR/ETFE,Stainless steel, Polyacetal,PP,Silcion Rubber,PP,/PC,HDPE- Sterile – 5- E& D</p> <p>2. IV Cannula with wings and without port/IV Cannulae- Polypropylene,PTFE/FEP/PUR/ETFE,Stainless steel, Polyacetal,PP,Silcion Rubber,PP,/PC,HDPE- Sterile – 5- E& D</p> <p>3. IV Cannula without wings and without port/IV Cannulae- Polypropylene,PTFE/FEP/PUR/ETFE,Stainless steel, Polyacetal,PP,Silcion Rubber,PP,/PC,HDPE- Sterile – 5- E& D</p> <p>4. IV Cannula with small wings and without port/IV Cannulae- Polypropylene,PTFE/FEP/PUR/ETFE,Stainless steel, Polyacetal,PP,Silcion Rubber,PP,/PC,HDPE- Sterile – 5- E& D</p> <p>5.Suction catheter/catheters-PVC/ABS-Sterile -5- E& D</p> <p>6.Needle Catheter/Catheters-PVC-Sterile -5- E&D</p> <p>7.Infant Feeding tube/catheters- PVC-Sterile -5- E&D</p> <p>8.Ryle tube/catheters-PVC-Sterile -5- E&D</p> <p>9.Foleys balloon catheter/cathers-Natural rubber altex/cilicon-Sterile -5- E& D</p> <p>10.Thoracic catheter(straight,angles,trocar)/catheters- PVC/HDPE/ABS/SS-Sterile -5- E& D</p> <p>11.Oxygen Catheter/catheters-PVC-Steril e-5- E& D</p> <p>12.Yankauer Suction set/catheters-Polyvinyl chloride,ABS Transparent/Rigid PVC,PVC-Sterile -5- E& D</p> <p>13.Close Wound Suction Unit/Catheters- LDPE,PP,PVC,ABS,Stainless Steel-Sterile 5- E &D</p> <p>14.Central Venous Catheter/Cathers-Polyethylene/PC/Nylon 6/PP/HDPE/SS)-Sterile -5- E&D</p> <p>15. Central Venous Catheter/Cathers-Polyurethane/PC/PE/PP/SS-Sterile -5- E&D</p>	CLAA/MD/MS A/03/2012-DC	<p>Mfg Licno No.PD/127 Validity from 01.06.2017 to 31.05.2022</p> <p>(30-08-2017) (Approval for additional Product)</p>

93.	M/s Devon Innovations Pvt.Ltd., Plot No.27A,Near State Bank of Inida,Electronic City, Phase-I,Hosur Road,Banagalore-560100	1.Catheters-Stainless Steel-Sterile by ETO -3 2.Biliary Duct Drainage Catheter-Medical Grade Polyurethane-Sterile by ETO-5 3.Percutaneous Pigtail Nephrostomy Catheter/Set-Medical grade Polyurethane-Sterile by ETO -5 4.External Biliary Drainage Catheter- Medical grade Polyurethane-Sterile by ETO -5 5. Percutaneous Malecot Nephrostomy Catheter/Set-- Medical grade Polyurethane,Polyvinyl Chloride,Ethylene-Vinyl Acetate-Sterile by ETO -5 6.Naso VBiliary Drainage Catheter- - Medical grade Polyurethane-Sterile by ETO -5 7.Nephrostomy Catheter- Medical grade Polyvin yl Chloride-Sterile by ETO -5 8.Percutaneous Biliary Drainage Catheter- - Medical grade Polyurethane-Sterile by ETO -5 9.Pleural Drainage Catheter/Set-- Medical grade Polyurethane-Sterile by ETO -5 10.Suprapubic Cystostomy Catheter/Set-Medical Grade Ethylene-vinyl acetate- Sterile ETO-3 11.Tieman Catheter-Medical Grade Polyvinyl chloride-Sterile by ETO-2 12.Percutaneous Pigtail Suprapubic Catheter/Set-Medical grade Polyuretahane-Sterile by ETO-5 13.Ureteral Catheter- - Medical grade Polyurethane-Sterile by ETO -5 16.Urinary Catheter- Medical grade Polyvinyl Chloride-Sterile by ETO -3 17.Dilator/Set-Medical grade Polyurethane,Ethylene-vinyl acetate,Polypropylene+Ethylene-vinylacetate,polyurethane-Sterile by ETO-3 18.Ureteral Indwelling double pigtail stenet/set (short/long term) – Medical grade polyurethane carbothane-Sterile by ETO-4 19.Long Abcess Drainage Catheter-- Medical grade Polyurethane-Sterile by ETO -5 20.Urethral stent/Set - - Medical grade Polyurethane-Sterile by ETO -5 21.Endoureterotomy stent/Set-- Medical grade Polyurethane-Sterile by	CLAA/MD/KAR/01/2010-DC	Mfg Licno No.KTK/28/395/2009 Validity from 16-06-2016 to 15-06-2021 (30-08-2017) (Approval for additional Product)
-----	--	--	------------------------	---

		<p>ETO -5</p> <p>22.Amplatz Renal Introducer/Set- Poly tetra Fluoro Ethylene-Sterile by ETO-4</p> <p>23.Biliary Stent/Set-- Medical grade Polyurethane-Sterile by ETO -5</p> <p>24.Internal external drainage catheter/Set-- Medical grade Polyurethane-Sterile by ETO -4</p> <p>24.Biliary Pancreatic stent/Set-- Medical grade Polyurethane-Sterile by ETO -4</p> <p>25.Cutting loop (single/double stent)-stainless steel-Sterile by ETO-4</p> <p>26.Multi-band logator barrel-Latex Free Protein-Sterile by ETO-4</p> <p>27.Foriegn Body retriever dormia basket-Stainless Steel-Sterile by ETO - 4</p> <p>28.Guidewire-Stainless Steel-Sterile by ETO-3</p> <p>29.Endoscopic Jejunal Feeding tube-- Medical grade Polyurethane,Medical grade polyvinyl chloride-Sterile by ETO -3</p> <p>30.Hemodialysis Catheter Set- Medical grade polyvinyl chloride-Sterile by ETO -5</p> <p>31.Vario Drain- Medical grade polyvinyl chloride-Sterile by ETO -3</p> <p>32.Orvac- Medical grade polyvinyl chloride-Sterile by ETO -5</p> <p>33.Suction Catheter- Medical grade polyvinyl chloride-Sterile by ETO - 5</p> <p>34.Klersite Cannula - Medical grade polyvinyl chloride-Sterile by ETO - 5</p> <p>35.Femoral Component Muller type STD- Stainless Steel 316L-Sterile by Gamma-5</p> <p>36.Femoral Component charmley type -STD- Stainless Steel 316L-Sterile by Gamma- 5</p> <p>37.Modular Bipolar Assemble with Head- Stainless Steel 316 and UHMWPE-Sterile by Gamma-5</p> <p>38.Endopro-Austin Moore- Stainless Steel 316L-Sterile by Gamma-5</p> <p>39.Acetular cup muller TypeLPW/STD-UHMWPE-Sterile by Gamma-5</p> <p>40.Integral bipolar for charnley- Stainless Steel316L and UHMWPE-Sterile by Gamma-5</p> <p>41.Acetular Cup Charnley type-LPW/STD- UHMWPE-Sterile by Gamma-5</p>		
--	--	--	--	--

		<p>42.Modular bipolar for muller – Stainless Steel 316L and UHMWPE- Sterile by Gamma-5</p> <p>43.Fixed Bipolar prosthesis Type-Fenestrated/Non fenestrated- Stainless Steel 316L and UHMWPE-Sterile by Gamma-5</p> <p>44.Collarless Stem with centralizer- Stainless Steel 316L- Sterile by Gamma -5</p> <p>45.Modular Head-Muller type STD- Stainless Steel 316L-Sterile by Gamma-5</p> <p>46.Enderails-Stainless Steel 316L- Unsterile</p> <p>47.Kirschner Wire/Threaded kirschner wire-Stainless Steel 316L-Sterile by Gamma -5</p>		
94.	M/s Bio-Med Healthcare Product Pvt.Ltd., Plot No.62-65,Sector-6A,DISCUL,Ranipur-249403,Haridwar,Uttarakh and	<p>1.Ritemed-Polypropylene,stainless Steel-Sterile –Export</p> <p>2. Manaflo- Polypropylene,stainless Steel-Sterile –Export</p>	CLAA/MD/UA/02/2007-DC	<p>Mfg Licno No.123/UA/SC/P-2007</p> <p>Validity from 02/11/2012 to 01/11/2017</p> <p>(12-09-2017)</p> <p>(Approval for additional Product)</p>
95.	M/s Bio-Med Healthcare Product Pvt.Ltd., Plot No.62-65,Sector-6A,DISCUL,Ranipur-249403,Haridwar,Uttarakh and	<p>1.Kings- Polypropylene,stainless Steel-Sterile –Export</p> <p>2.Health- Polypropylene,stainless Steel-Sterile –Export</p> <p>3.Top Care- Polypropylene,stainless Steel-Sterile –Export</p>	CLAA/MD/UA/02/2007-DC	<p>Mfg Licno No.123/UA/SC/P-2007</p> <p>Validity from 02/11/2012 to 01/11/2017</p> <p>(12-09-2017)</p> <p>(Approval for additional Product)</p>

96.	M/s Disposafe Health and Lifer Care Ltd., Plot No.1 &2, Sector-59,Ballabgarh-121004	1.Periflow(with/without wings and with/without injection port)- Polypropylene ,POlythylene,Polyacetal,Stainless Steel(AISI304)-Sterile -5- E& D 2..Periflow(with integrated three way stopcock)- Polypropylene ,POlythylene,Polyacetal,Stainless Steel(AISI304)-Sterile -5- E& D 3..Periflow(with safety feature)- Polypropylene ,POlythylene,Polyacetal,Stainless Steel(AISI304)-Sterile -5- E& D 4..Periflow(with small wings and with injection port)- Polypropylene ,POlythylene,Polyacetal,Stainless Steel(AISI304)-Sterile -5- E& D	CLAA/MD/HR /01/2013-DC	Mfg Licno No.714-B (H) Validity from 27-05-2013 to 26-05-2018 (12-09-2017) (Approval for additional Product)
97.	M/s SPM Medicare Pvt.Ltd., B-40,phase-II,Noida,Gautam Budh Nagar-201305,UP	1.IV Cannula with/without injection port and with/without wings- Teflon,PP,LDPE,Silicon,Tube,Needle-Sterile-5- E& D 2.IV Cannula without injection port and with small wings- eflon,PP,LDPE,Needle-Sterile-5-E & D	CLAA/MD/UP 02/2017-DC	Mfg Licno No. Validity from (14-09-2017) (Approval for additional Product)
98.	M/s Jay Harsh Ortho Tech, 965/K-1,GIDC, Makarpura, Vododara-390010,Gujarat,India	1.Bone Screws(Corical Screw, Cancellous Screw,Malleolar Screw,Cannulated Cancellous Screw, Hip Screw, Locking Screw, Locking Bolt, Hip Pin,Hip Bolt,Interference Screw)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D 2.Bone Plates (Dynamic Compression Hole, Low Contect Dynamic Compression Holes,Tubular Plates,Reconstruction Plates,Locking Plate,Hip Plates,Endo Buttons)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D 3.Nails(Femur Nails,Tibia Nails,Humerus Nails,R/U Nails,Supracondylar Nails,SSTN Nails,Trochantric Nails,Ankle Nails,Fibula Nails- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D 4.Pins and Wires (K-Wires,Threaded Wires,Tapered Threaded pins, Schanz Pins,Other Misc.Wires,S.T.Pins)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D 5.Spinal Impants (Spinal Plates, Spinal AScrews, Spinal Hooks, Pedicle Screw, Spinal Rod, Spinal Cages, Spinal Mesh,connectors)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D	CLAA/MD/GU J/04/2016-DC	Mfg Licno No. G/28/1611 Validity from 02-05-2017 to 01-05-2022 (22-09-2017) (Approval for fresh License)

99.	M/s Innovative Otho Surgicals Pvt.Ltd., Plot No.41,Road No.5,Sardar Estate,Vododara-390019,Gujarat	<p>1.Bone Nail (Femur Nails,Tibia Nails, Humerus Nails, Radius Ulna Nails,Supracon dylar Nails,All fracture femur nails,all fracture tibia nail,trochantric nails,ankle nails, fibula nails - S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D</p> <p>2.Bone Plate (Dynamic Compression Plate,Round Hole Plate,Tubular Plate,Reconstruction Plate, Angle stable plate, Hip plate, endo buttons)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D</p> <p>3.Bone Screws (Cortical Screw,Cancellous Screw, Malleolar Screw, Angle stable cancellous screw, Cannulated Cancellous Screw, Angle stable Screw, Locking Bolt, Hip Pin- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D</p> <p>4.Pins and Wires (K-Wires, Threaded Wires, Tapered threaded pins, Schanz Pins, Others Misc.Wires and S T Pins)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D</p> <p>5.Spinal Implants (Spinal Plate,spinal Screw, pedicle screw,spinal cages, spinal mesh,connectors)- S.S,316L,S.S.,316LVM,Titanium/Titanium Alloy)- Non-Sterile – E& D</p>	CLAA/MD/GU J/05/2011-DC	<p>Mfg Licno No. G/28/1373 Validity from 27-05-2016 to 26-05-2021</p> <p>(22-09-2017) (Approval for fresh License)</p>
100.	M/s Sharma Surgical and Engg.Pvt.ltd.,654 & 655,GIDC Estate, Waghodia-391760,Gujarat,India	<p>1.Locking Plates (DHS Plates,DCS 95Plate, Proximal Femoral Plate – left and right ,Femoral Condylar Buuttres plates (left and right LISS Plate),Dynamic COMression Plate(small,narrow and broad),Latreal tibial head buttress plates 4.5mm (left and right),Medial proximal tibial head buttress paltes (T type left and right)(4.0mm and 5.0mm), Distal Medial tibial plate (metaphyseal distal tibia plate,anatomical distal medial tibia plate (admt.)left and right), proximal humerus plate,low profile humerus plate,Philosh 4.0 & 5.0mm proximal humerus plate (small and larg), distal radius plates,T-Oblique plate left and right), reconstruction plates 3.5mm and 4.5mm, L-Plate (left and right)4.5mm,L-Buttress plate (left and right 4.5mm), T-Plate (4.5mm), T-Buttress plate(4.5mm), cloverleaf palte, calcanium locked palte,1/3rd tubular plates,semi tubular plares, distal humerus locking pplate, anatomical distal humerus locked plate (left and right), medial distal humerus loked plate (left and right),distal humerus lateral posterior (left and right),4mm DCLC Distal lateral tibial head buttress paltes (left and</p>	CLAA/MD/GU J/14/2009-DC	<p>Mfg Licno No. G/28/1373 Validity from 27-05-2016 to 26-05-2021</p> <p>(26-09-2017) (Approval for Renewal Products)</p>

		<p>right antrolateral dstal tibia),DCLC Low band distal medial tibia head buttress plates (left and right),Extra articular distal humerus locking plate (left and right), DCLC Proximal ulna plates (left and right), DCLC 20,23,26mm angle distal radius plate (left and right), DCLC 23,28mm T-Oblique plate(left and right), DCLC 24,31mm T Plate, Clavicle hook plate (left and right)(small,medium,large)Clavicle lateral medial anterior plate,distal, fibula plate(left and right),postermedial proximal tibia locking plate(left and righ), fuddu tibial ostemotomy plate,fuddu femor distal osteotomy locking plate,distal pilon plate,DHS Trochanteric support plate (with or without flange), recon curved round and DC Hole locking plate (3.5mm) (single and dual lock),Ankle Arthodesis Nail,GKT Distal tibia Nail (cannulated and solid)(A.P.&M.L),Locking Bolt, Hlical Blade, Flexible Nail,Internal Head (cemented and uncemented total hip replacement), Acetabular cup-10-20life cemented total hip replacement,Acetabulum acge (Cliff and without sliff)(cemented) total hip replacement,Glenoid (cemented and uncemented shoulder prosthesis system, eccenter (cemented or uncemented shoulder prostehsis system,Tibial inserts (fixed bearing PS, cemented and uncemented,total hip replacement,Tibial Wedges Hinged knee (5.0MM(Thickness cemented and uncemented total knee replacement),Proximal Femur Prosthesis,Total Femur Mega Prosthesis,Proximal Humerus Mega Prothesis,Distal Humerus Mega Prosthesis,Total Humerus Mega Prothesis, Radial Hea (Radial Head Rweplacement)-S.S,316L,S.S, 316LVM,Titantium Alloy,CoCr Alloy-Non-Sterile –E & D</p>		
101.	M/s Ribbel International Limited,20 th Mile, Jatheri Road,P.O.Rai,Sonipat,Haryana	1.Eetacathfoley catheter- Natural Latex-Sterile and Non Sterile -5-E&D	CLAA/MD/HR /03/2014-DC	Mfg Licno No. Validity from (26-09-2017) (Approval for Additional Products)

102.	M/s Meril Healthcare Pvt.Ltd., Survey No.135/2/B AND 174/2,1 st floor, H1-H3, Meril Park, Muktanand Marg, Chala,Vapi-396191	1.Hip Replacement System-Ti-6Al ELI &Hydroxyapatite(coating material) , High N2 SS, Commerically pure titanium, UHMWPE,HXLPE,E-CiMa-Non Sterile (to be sterilized by Gamma Radiation) & Sterile (ETO Gas)-5- E& D	CLAA/MD/GU J/05/2013-DC	Mfg Licno No. G/28/1465 Validity from 01.04.2014 to 31.01.2019 (26-09-2017) (Approval for Additional Products)
103.	M/s Nebula Surgical Pvt.Ltd., GIDC ,Plot NO.G/739-740,b/h Telephone tower, Metoda Dist.Rajkot	1.Bone Screws (Cortical Screw 2.5mm, 2.7mm,3.5mm,4.5mm , DHS Lag Screw, Mini Screw, Proximal Bolt for 6.4mm,8.0mm,Locking screw 2.7mm,3.5mm,4.9mm, Locking cancellous screw 3.5mm,4.9mm,6.5mm,4.9mm,7.3mm,7.3 C.C. Non licking screw mm)-SS 316L & Titanium- Non Sterile – E &D 2.Bone Plates (VAP Screw, VAP Cortical Screw, VAP cannulated cancellous screw (VAP C.C Screw), VAP Locking screw 4.9mm fixed head, VAP cancellous Screw, VAP Spacer, VAP locking cap, PFNA2 Blade, Small D.C Plate, Narrow D.C Plate, Broad D.C Plate, Small L.C.D.C Plate, Broad L.C.D.C Plate,Olique T Plate, DHS Plate-Long/short Barrel, DHS Plate-Long/short barrel with trochentric support,95 DCS Plate, DCS Condylar balde plate, small locking plate,pelvic curved plate, pelvic straight plate, Plevic J-Shape plate, Pelvic Symphysis plate, Pelvic straight locking plate,Proximal femur locking plate,Anterolateral distal tibial locking plate,Posterior Medical Proximal Tibial Locking plate, Distal tibial medial locking plate, superior anterior clavical locking plate,clavical hook locking plate, distal medial humerus metaphyseal locking plate,Distal humerus Dorsolateral locking plate without support, Distal radius volar locking plate,Olecronan Locking plate, U Distal Radius Dorsal Locking plate, Straight distal radius dorsal locking plae, Oblique distal radius dorsal locking plate, T Distal Radius Dorsal locking plate, VAP distal femur locking plate, Extra articular distal humerus locking plate. Fibula Locking plate, Proximal tibial medial locking plate,Tomofix Locking plate,Distal Radius Locking plate,Optional Buttress locking plate,Extra Articular distal radius locking plate,SS 316L & Titanium-Non-Sterile –E & D	CLAA/MD/GU J/15/2010-DC	Mfg Licno No. G/28/1355 Validity from 29/10/2015 to 28/10/2020 (28-09-2017) (Approval for Renewal Products)

		<p>3.Endo Prosthesis- ,SS 316L & Titanium-Non-Sterile –E & D</p> <p>4.IM Nails (Tibia Nail-Gundrill,Distal Jig, Multifunction Tibia Nail,Cyrus Tibia Nail,Femur Nail-13mm. P.F Nail-25cm (short,Left and right),Prximal Humerus Nail)- SS 316L & Titanium-Non-Sterile –E & D</p> <p>5.Spine Implants (Mono Axial Screw, Poly Axial Screw, TLIF Cage-Peek, PLIF Cage-peek, Facet Hook,Supra Laminar Hook,Infra Laminar Hook, Laminar Hook Narrow Balde, Laminar Hook Wide Balde, Spike Washer) – SS 316L & Titanium,Titanium/Peek) – Non Sterile – E &D</p> <p>6.Cervical Implants -,SS 316L & Titanium-Non-Sterile –E & D</p>		
104.	M/s Surgimade ,Plot No.459/15/B,GIDC, Por-Ramangamdi,At-Por, Dist.Vadodara	<p>1.Tibia Nail – SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>2.Femur Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>3.Retrograde Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>4.Gam Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>5.Expert Gam Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>6.Recon Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>7.P.F.Nail- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>8.Supracondy Iar Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>9.Humerus Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>10.Radius/Ulna Nail - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>11.Enders Nail- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>12.Ankle Nail- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>13.Locking Bolt- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p>	CLAA/MD/GU J/01/2017-DC	<p>Mfg Licno No. G/28/1591</p> <p>Validity from 04/01/2017 to 03/01/2022</p> <p>(28-09-2017)</p> <p>(Approval for License Products)</p>

	<p>14.DHS Plate - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>15.DCS Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>16.Proximal Femur Plate - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>17.CBP Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>18.Proximal humerus locking plate - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>19.Low Profile Proximal Humerus locking plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>20.5.0mm philosh plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>21.T-Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>22.L-Plate - SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>23.L-Buttess Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>24.Cloverleaf plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>25.1/3rd Tubular plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>26.Volar Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>27.Olccranon Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>28.Clavical Hook Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>29.Recon Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>Locking Screw- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>30.Trochnetric Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p>		
--	---	--	--

		<p>31.Hebert Screw- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>32.Washer- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>33.Mono Reduction Screws- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>34.San Cage- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>35.Stapple- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>36.Cervical Spinal Plate- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>37.K-Wires- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>38.S T Pins- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>39.Dehnam Pin- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>40.Knowels Pin- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>41.Rush Nail- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p> <p>42.Wire Reels- SS 316L,S.S. 316 L.V.M ,Titanium/Titanium Alloy – Non Sterile – E&D</p>		
105.	M/s Siri Surgical Parts Pvt.Ltd., Plot No.21/5,1 st Floor,Phae-III,IDA Jeedimetla ,Hyderabad-500055	<p>1.Bone Plates (DCP Small Plates 3.5mm, DCP Narrow Plates 4.5mm, DCP Broad Plates 4.5mm, LCDCP Small Plates, LCDCP Narrow Plates, LCDCP Broad Plates, Reconstruction Plates, Small T Plates 3X3 Hole, Small T Plate 4X4 Hole,Cloverleaf plates, Ellis Buttress plate for radius, Smaill T plate right angles 3.5mm, Small T Plate Oblique Angles 3.5mm,Caclcanium Plate,T-Plate, T-buttress Plate Left, Y Reconsturction plate, Lateral tibial plates leftD.H.S Barrel Plates SB/RH120 ,125,130,135,140,145Degree, Condylar Femoral Buttress plate-left and right, Cobra head plates,,)- SS 316L/Titanium-Non-Sterile/Sterile - 5- E&D</p> <p>2. Bone Screws (2.7mm cortical Screws, 3.5mm cortical screws-14 TPI,</p>	CLAA/MD/AP/01/2017-DC	<p>Mfg Licno No. 16/RR/TS/2017/S/G Validity from 25/04/2017 to 24/04/2022</p> <p>(04-10-2017) (Approval for License Products)</p>

		<p>3.5mm Cortical screws-20 TPI, 4.5mm Cortical Screws-20 TPI, 4.5mm Cortical screws,Cancellous screws, malleolar screws, cannulated cancellous screws,CAP Screws, ACL Screws, Herbert Screws, Interlocking Screws,DHS Screws)- SS 316L/Titanium – Non Sterile /Sterile -5- E & D</p> <p>3.Bone Nails (interlocking tibia nail, solid interlocking for tibia, syrus type tibia nail,Femur nails, humerus nail,proximal humerus nails, elastic nails,square nails, enders nails, rush nails, PFN Long Nails,PFN Short Nails,Recon Nails, K Nails,Schanz PINS,Tapered threaded pins, Steinman pins) SS 316L/Titanium – Non Sterile /Sterile -5- E & D</p> <p>4.Locking Plates (LC DCP Small 3.5mm, 4.5mm, LCDCP Broad 4.5mm,LC Humrus Plates, LC Tibia Plates, LC Femoral Plates, LCP cortical screws, LCP Cancellous srews, ,LC DHS Plates, CLC Olecrenon plates, LC _DISTAL DORSAL HUMURAL PLATES, lc Calneal Plates) SS 316L/Titanium – Non Sterile /Sterile -5- E & D</p>		
106.	M/s Harsoria Healthcare Pvt.Ltd., 110-111,Phae-IV, Udyog Vihar,Gurgaon-122015,Haryana	<p>1.IV Cannula Safety (Safety IV Cannula with catheter and injection valve and with wings) – Polypropylenen HDPE ABS Poly Acteal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 lokcing strip in AISI302) – Sterile -5- E &D</p> <p>2. IV Cannula Safety (Safety IV Cannula with catheter and without injection valve and with small wings) – Polypropylenen HDPE ABS Poly Acteal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 lokcing strip in AISI302) – Sterile -5- E &D</p> <p>3. IV Cannula Safety (Safety IV Cannula with catheter and without injection valve and without small wings) – Polypropylenen HDPE ABS Poly Acteal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 lokcing strip in AISI302) – Sterile -5- E &D</p> <p>4.Suresite IV Auto (Sterile Disposable safety IV Cannula with catheter for single use only) Polypropylenen HDPE ABS Poly Acteal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 lokcing strip in AISI302) – Sterile -5- E &D</p> <p>5.Suresite IV Slide (Safety IV Cannula with catheter and without injection valve and without small wings)- Polypropylenen HDPE ABS Poly Acteal Catheter (PTFE/PUR/FEP/SS Needle in AISI304 lokcing strip in AISI302), Porous Plug – Sterile -5- E &D</p>	CLAA/MD/HR /04/2006-DC (Pt-2)	<p>Mfg Licno No. 16/RR/TS/2017/S/G</p> <p>Validity from 19/09/2016 to 18/09/2021</p> <p>(04-10-2017)</p> <p>(Approval for License Products)</p>

107.	M/s MerilLife sciences Pvt.Ltd., Survey No:-135/a and 135/B, 174/2 (M1-M2,Meril park, 4 th floor, near G.M.Bilakhia stadium, muktanand marg, chala, vapi,valsad-396191,Gujarat	1.Barre Metal aortic stent and stent system (little angels)- L 605 Cobalt chromium stent- Sterile – 5- E& D 2.Covered Bare Metal Aortic Stent and Stent System (Little Angels) – L 605 Chromium Stent covered with Eptfe (Poly Tetra Fluoro Ethylene) Material –Sterile -3 – E& D	CLAA/MD/GU J/01/2008-DC (Pvt.4)	Mfg Licno No. G/28/1273 Validity from 20/06/2013 to 19/06/2018 (04-10-2017) (Approval for License Products)
108.	M/s Sterimed MD (P) Ltd., 38 KM. Stone Delhi – Rohtak Highway ,Bahadurgarh (Haryana)	1.Foley Ballon Catheter 100 % Silicone – Medical grade silicone – Sterile -5- E& D	CLAA/MD/HR /02/2008-DC	Mfg Licno No. 661-B (H) Validity from 23/06/2013 to 22/06/2018 (10-10-2017) (Approval for License Products)
109.	M/s Disposafe Helath and Life Care Ltd., Plot No.1 &2, Sector -59,Pahse – II,Ballabgarh-121004,Faridabad,Haryana	1.Medeco (Cannula with/without wings and with/without injection Port) – Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) – Sterile- 5- E&D 2.Ipd Lifecare (Cannula with/without wings and with/without injection Port,Sterile IV Cannula with safety Feature) – Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D 3.Pen Cath (Sterile intravenous cannula without wings and without injection port)– Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D 4.All Pro (Sterile IV Cannula with wings and injection port, Sterile IV Cannula with safety feature) - Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D 5.Shilohcath (Sterile IV Cannula with wings and injection port, , Sterile IV Cannula with safety feature) - Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D 6.Mecyflon,Medisafe, Nirvana,IDL,Precision,Golden+,Q-Flon,Lifecare,IPD Lifecare(Sterile Intravenous cannula with/without small wings and with/without port)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D 7. Mecyflon,Medisafe, Nirvana,IDL,Precision,Golden+,Q-	CLAA/MD/HR /01/2013-DC	Mfg Licno No. 714-B (H) Validity from 27/05/2013 to 26/05/2018 (10-10-2017) (Approval for Additional Products)

		<p>Flon,Lifecare,IPD Lifecare (Sterile Intravenous Cannula with safety feature, SSterile Intravenous Cannuls with integrated three way stopcock)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D</p> <p>8.Hexacare (Sterile IV Cannula with/without wings and with/without injection port, Sterile IV Cannula with safety Feature) - Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D</p> <p>9.TOPDESC (Sterile IV Cannula without wings and without injection port) Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D</p> <p>10.FAB (Sterile IV Cannula with wings and with injection port)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) – Sterile- 5- E&D</p> <p>11.FAB Cath (Sterile IV Cannula integrated three way stopcock)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) – Sterile- 5- E&D</p> <p>12.NET-ON (Sterile IV Cannula without wings and without injection port)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D</p> <p>13.Venfic (Sterile IV Cannula with/without wings and with/without injection port)- Polypropylene ,Polyethylene,Polyacetal, Stainless Steel (AISI304) –Sterile- 5- E&D</p>		
110.	M/s AL Optics Ltd.,84/1,Aavin Dairy Firm Road, SIDCO Industrial Estate, Ambattur,Chennai	Intra Ocular Lens-Hydroxy Ethyl Methacrylate (HEMA)-Yellow- Sterile-3- E	CLAA/MD/TN /02/2016-DC	<p>Mfg Licno No. D.DIS No.2017/D1/2/2017 Validity from</p> <p>(16-10-2017) (Approval for Additional Products)</p>

111.	M/s Neo Surgical,Plot No.G-306,GIDC Almighty Gate, matoda-360021,Rajkot	<p>1.Expert Tibia Nail (Tibia Nail)- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>2.EXPERT Femur NI (Femoral Nil)- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>3.Humorous Nail - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>4.PFNA Nail - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>5.PFNA Nail left and Right - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>6.Denhan Pin - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>7.Schanz Pin - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>8.Lateral Tibia Locking plate left and right- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>9.DyNAMIC Reconsturction plate - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>10.Dynamic Reconstruction locking plate - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>11.Anatomical distal medial tibia locking plate – left and right- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>12.Locking head screw - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>13. Cancellous locking screw- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>14.Cancellous Locking Screw-32 MM Thread- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>15.Cacellous locking screw full thread - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>16.Cannulated cacellous screw – full thread- SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>17.PFNA Blade - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p> <p>Anchor Screw - SS316L,SS316 LVM, Titanium- Non Sterile – 5- E &D</p>	CLAA/MD/GU J/01/2016-DC	<p>Mfg Licno No. G/28/1585 Validity from 22-09-2016 to 21-09-2021</p> <p>(16-10-2017) (Approval for Additional Products)</p>
------	---	--	-------------------------	---

112.	M/s Poly Medicure Ltd., Plot No:-115- 116,HSIDC,SECTOR- 59,Ballabgarh,Faridabad- 121004,Haryan	1.Gatrosopic flushing catheter – PTFE/FEP/PUR,PE,PP,SS,PA – Sterile for single use only- 5- E & D 2.Fionavet - PTFE/FEP/PUR,PE,PP,SS,PA – Sterile for single use only- 5- Export 3.Exelint (Sterile Haemodialysis catheter)- PTFE/FEP/PUR,PE,PP,SS,PA – Sterile for single use only- 5- E & D 4. Exelint (SterileEndotracheal tube)- PTFE/FEP/PUR,PE,PP,SS,PA – Sterile for single use only- 5- Export 5.Exevent (Double Lumen Central Venous Catheter) –PVC, ABS, PC,SS,PP-Sterile for single use only -5-Export 6. .Exevent (triple Lumen Central Venous Catheter) –PVC, ABS, PC,SS,PP-Sterile for single use only -5-Export 7.Medmax (IV Cannula)-PTFE/FEP/PUR,PE,PP,SS,PA- Sterile for single use only -5-Export 8.Petit FLO(IV Cannula) - PTFE/FEP/PUR,PE,PP,SS,PA- Sterile for single use only -5-Export	CLAA/MD/HR /01/2006-DC	Mfg Licno No.627-B (H) Validity from 06.03.2016 to 05.03.2021 (16-10-2017) (Approval for Additional Products)
113.	M/s Disposafe Health and Life Care Ltd., Plot No.1 &2, Sector-59,Phase-II,Ballabgarh-121004	1.Intravenous Cannula-Polypropylene Polyethylene polyacetal,PC Stainless Steel (AISI304)-Sterile -5- E &D 2.BPI - Polypropylene Polyethylene polyacetal,PC Stainless Steel (AISI304)-Sterile -5- E &D 3.Periflow- Polypropylene Polyethylene polyacetal,PC Stainless Steel (AISI304)-Sterile -5- E &D 4. Periflow (Sterile closed wound suction set)- Polypropylene,Nylon,Polycarbonate, Polyxymethylene,Poly Vinyl Chloride,Stainless Steel Low density polyethylene-Sterile -5- E &D	CLAA/MD/HR /01/2013-DC	Mfg Licno No714-B (H) Validity from 27.05.2013 to 26.05.2018 (16-10-2017) (Approval for Additional Products)
114.	M/s Dalvi and Sons, G-10.Additional MIDC, Kudal, Dist:-Singhudurg, Maharshtra	1.Artho Button- TI.grade -5,UHMWP,Biocompatible for loop- Non Sterile -Domestic	CLAA/MD/MS /01/2015-DC	Mfg Licno No. MH/101512 Validity from 23/10/2015 (16-10-2017) (Approval for Additional Products)

115.	M/s Lars Medicare Pvt.Ltd., Kila No.16 & 17,village sultanpur, G T Road, Bahalgarh chowk-131021	1.Lifemed – Radio-opqae strip/clear Teflon, PP,POM, Silicon,HDPE,LDPE and Stainless Steel- Sterile -5- Export	CLAA/MD/HR /01/2008-DC	Mfg Licno No. D660-B (H) Validity from 04.02.2013 to 03.02.2018 (16-10-2017) (Approval for Additional Products)
116.	M/s Intermedics Devices Pravte Limited, #74,opp. Srinidhi College, CK Palya Road,Hommadevanahalli, Bannerghatta Road, Bangalore-560083	1.PCN Catheter- Polyurethane,PVC and ABS, PP- ETO-3- E &D 2.Malecot Catheter- Polyurethane,PVC and ABS, PP- ETO-3- E &D 3.Utheral Catheter- Polyurethane, ABS, Stainless Steel- ETO-3- E &D 4.NBDC Catheter- Polyurethane, ABS - ETO-3- E &D 5.PTBD Catheter- Polyurethane,PVC and ABS, PP- ETO-3- E &D 6.Internal External Drainage Catheter- Polyurethane,PVC and ABS, PP- ETO-3- E &D 7.PCN Catheter with Needle- Polyurethane,PVC and ABS, Stainless Steel, PP- ETO-3- E &D 8.Suprapubic Catheter- Polyurethane,PVC and ABS, Stainless Steel, PP- ETO-3- E &D 9.RE-Entry Malecot Catheter- Polyurethane,PVC and ABS, PP- ETO-3- E &D 10.Suprapublic Drainage Catheter Set-Polyurtahne,Urine Bag Connector, Guide Wire,IP Needle, Dilator,Schapel- Polyurethane,PVC and ABS, Stainless Steel, PP- ETO-3- E &D 11.Nephrostomy Drainage Cathter Set- Polyurethane,PVC,PTFE,SS,PU,SS - ETO-3- E &D	CLAA/MD/KR /02/2017-DC	Mfg Licno No. KTK/28/442/2017 Validity from 17/04/2017 to 16/04/2022 (21-10-2017) (Approval for Products)

117.	M/s Vertias Bioventions Pvt.Ltd., Plot No.5305,Phase-IV,GIDC, Vapi-396195	<p>1.Endorser-Stainless Steel tube-SS316L (Hard Drawn)-Sterile-24 months-Domestic</p> <p>2.Patron- Stainless Steel tube-SS316L (Hard Drawn)-Sterile-24 months-Domestic</p> <p>3.Hotron- Stainless Steel tube-SS316L (Hard Drawn)-Sterile-24 months-Domestic</p> <p>4.Flextron- Stainless Steel tube-SS316L (Hard Drawn)-Sterile-24 months-Domestic</p> <p>5. Endorser CrCo-Cobalt Chromium Tube-L605 (Hard Drawn)-Sterile-24 months-Domestic</p> <p>6. Patron Patron- Stainless Steel tube-SS316L (Hard Drawn)-Sterile-24 months-Domestic</p> <p>7. Hotron CrCo-Cobalt Chromium Tube-L605 (Hard Drawn)-Sterile-24 months-Domestic</p> <p>8. Flextron CrCo-Cobalt Chromium Tube-L605 (Hard Drawn)-Sterile-24 months-Domestic</p> <p>9.SiNew- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide)-Sterile-18 months-Domestic</p> <p>10.Fluence- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide)-Sterile-18 months-Domestic</p> <p>11.Freenix- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide)-Sterile-18 months-Domestic</p> <p>12.Flexicore- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide)-Sterile-18 months-Domestic</p> <p>13.Elimin8-R- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide)-Sterile-18 months-Domestic</p> <p>14.Forel- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>15.Fleam- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus</p>	CLAA/MD/GU J/08/2011-DC	<p>Mfg Licno No. G/28/1379 Validity from 22/07/2016 to 21/07/2021</p> <p>(06-11-2017) (Approval for Renewal Products)</p>
------	---	---	-------------------------	--

		<p>Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>16.Fistula- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>17.Fictile- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>18.Enovo- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>19.Ensconce- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>20.Espalier- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>21.Expediter- Cobalt Chromium Tube-L605 (Hard Drawn), Sirolimus Drug, 50/50 Poly (DL Lactide-co-Glycolide),poly (L—Lactide)-18 months- Domestic</p> <p>22.Inspire- Cobalt Chromium Tube-L605 (Hard Drawn), Paclitaxel Drug, 50/50 Poly (DL Lactide-co-Glycolide,Sterile-18 months- Domestic</p> <p>23.Zealous- Cobalt Chromium Tube-L605 (Hard Drawn), Paclitaxel Drug, 50/50 Poly (DL Lactide-co-Glycolide,Sterile-18 months- Domestic</p> <p>24.Zynoa- Cobalt Chromium Tube-L605 (Hard Drawn), Paclitaxel Drug, 50/50 Poly (DL Lactide-co-Glycolide,Sterile-18 months- Domestic</p> <p>25.Zontrax- Cobalt Chromium Tube-L605 (Hard Drawn), Paclitaxel Drug, 50/50 Poly (DL Lactide-co-Glycolide,Sterile-18 months- Domestic</p>		
--	--	--	--	--

118.	M/s Advanced Life Sciences Pvt.Ltd.,D-22,Okhla Industrial Area,Phase-1,New Delhi-110020	<p>1.Hemodialysis Catheters-Single/double/triple lumen catheter-PU Extensio Line-silicone /PVC Hub-PP,Clip-PP,G Wire-SS-304V,Dilator-PE,Hub-PP,Intiducer needle-SS,HUB-PC,Injection site cap-body-PC,Gasket natural rubber-Sterile-5-Export as well domestic market</p> <p>2.P.D Catheters/sets-trocar-SS-304,Catheter-PU,Scalpel blade-SS-304,Handle-PP-Sterile-5- Export as well domestic market</p> <p>3.CVC Catheters-It is used for giving intravenous fluids,blood, medication or nutrition-Sterile-5- Export as well domestic market</p> <p>4.PTCA Balloon Catheters-Shaft-SS 304V, Balloon Shaft-PA/PE, Proximal Shaft coating-PTFE, Balloon-PA/Nylon 44,X Ray marker-750 gold ring Tip-PEBAX 5533,Guidewire-PEBAX-7033-Sterile-5- Export as well domestic market</p> <p>5.Guiding Catheters-Hub-PC/PET Shaft-PEBAX,Catheter Tip-PEBAX 5533/Tungsten carbide-Sterile-5-Export as well domestic product</p> <p>6.Angiogrphery/R Adiology Catheters-Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>7.Thromus Aspiration Catheters- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>8.PTMC Catheters- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>9.PTA Baloon Catheters- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>10.Tenckhof Catheters for dialysis- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>11.Forgarty/Embo Lectomy Catheters- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>12.Wedge pressure/angio graphic balloon catheter/berman catheters-Sterile-5-Export as well domestic market</p> <p>13.Balloon dialation catheter- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>14.PTBD Catheters- Nylon/Teflon/Pebax/PU/PE-Sterile-5-Export as well domestic market</p> <p>15.Introducer sheath-Sheath-PU,Tube-PVC,Handle-PE,GWire- core wire-SS 304 SPRING WIRE-ss304 COATING-ptfe Dispenser-HD/PP,HUB-Hub Tip cover,pusher and grip-PP,Needle-SS,HUB-</p>	CLAA/MD/DL /01/2011-DC	Mfg Licno No.S(0024)/12/MD Validity from 18/05/2017 to 17/05/2022 (06-11-2017) (Approval for Renewal Products)
------	---	--	------------------------	---

		<p>PC,Protector-HDPE-Sterile-5- Export as well domestic market</p> <p>16.Guidewire/S Steel/PTFE/Hydrophilic- GWire- core wire-SS 304 SPRING WIRE-ss304 COATING-ptfe Dispenser-HD/PP,HUB-Hub Tip cover,pusher and grip-PP, Hydrowire-Core wire nitinol coating-hydrophilic dispenser-HD/PP,HUB-Hub, Tip cover, usher and grip-PP-Sterile-5- Export as well domestic market</p> <p>17.Occlusion devices ASD,PDA,VSD,PLUGS-Wire mesh-Nitinol polyester febric membrane-teraline SS Bushing –stainless steel 316L-Sterile-5-Export as well domestic market</p> <p>18.Occlusion devices delivery system/cable- Sterile-5-Export as well domestic market</p> <p>19.Coronary stent (S Steel) systems-Stainless Steel-Sterile-5- Export as well domestic market</p> <p>20.Coronary sten(cobalt chromium) systems-Phynox (Co-Cr Alloy)-Sterile -5- Export as well domestic market</p> <p>21.Peripheral/renal stent balloon mounted/self expanding-Stainless Steel/Nitinol Co-Cr alloy- Sterile -5- Export as well domestic market</p> <p>22.Double J/Uthral stent, open end/closed multioop end multiloop both end multiloop-Nylon/Teflon/PEBX/PU/PE- Sterile -5- Export as well domestic market</p> <p>23.Biliary stent standard,amsterdam,straight,single pigtail, double pigtail,tenaben pancreatic,pancreatic with single pigtail-Nylon/Teflon/PEBX/PU/PE, Stainless Steel/nitinol Co-Cr alloy-Sterile-5- Export as well domestic market</p> <p>24.Urology/gastro/esophageal/biliary stents-Sterile-5-Export was well as domestic market</p> <p>25Arterial/venous/cadioplegia cannula-Stainless Steel-Sterile-5-Export was well as domestic market</p> <p>26.PTBD Set/Drainage Set-Nylon/Teflon/PEBXPU/PE-Sterile-5-Export as well as domestic market</p> <p>27.IV Cannula/IV Port-Catheter-PU/PTFE/FPE Needle-SS 304-Sterile -5- Export as well as domestic market</p>		
--	--	---	--	--

119.	M/s Degania Medical Devices Pvt.Ltd., Plot No:- 275,Sector-6,IMT Manesar, Gurugram-121004,Haryana	1.Bare Catheter(Maya)/Vega- Stainless STEEL-Non Sterile –Export 2.PEGASE Aspiration catheter-Thermoplastic Polymer-Non Sterile – Export 3.Y-Conenctor-Polypropylene-Non Sterile –Export	CLAA/MD/HR /04/2013-DC	Mfg Licno No. 730 B (H) Validity from 07/07/2014 to 06/07/2019 (08-11-2017) (Approval for Additional Products)
120.	M/s Unisur Lifecare Pvt.Ltd.,No.15/1,2,3,Andra halli Main Road, Acharya Indl.,Complex, Vishwanedam Post,Bangalore-560091	1.M611-Polypropylene-Sterile by ETO-5-E&D 2.M 715- Polypropylene-Sterile by ETO-5-E&D 3. M 115- Polypropylene-Sterile by ETO-5-E&D 4.M 151- Polypropylene-Sterile by ETO-5-E&D 5.M 303- Polypropylene-Sterile by ETO-5-E&D	CLAA/MD/KA R/01/2017-DC	Mfg Licno No. KTK/28/441/2016 Validity from 17/03/2017 to 18/03/2022 (09-11-2017) (Approval for Additional Products)
121.	M/s Harsoira Healthcare Pvt.Ltd.,110-111,Phase-IV,Udyog Vihar,Gurgaon-122015,Haryana	1.Biociv(IV Cannula with Catheter)-Polypropylene HDPE,ABS,Polyacetal catheter(PTFE/PUR/FEP/SS Needle in AISI304)-Sterile -5-E&D 2. Biociv(Safety IV Cannula with catheter)-Polypropylene HDPE,ABS,Polyacetal catheter(PTFE/PUR/FEP/SS Needle in AISI304,locking strip AISI302 Porous Plug)-Sterile -5-E&D 3.Bullpup(Safety IV Cannula with catheter)-Polypropylene HDPE,ABS,Polyacetal catheter(PTFE/PUR/FEP/SS Needle in AISI304,locking strip AISI302-Sterile- 5- E &D 4.Harsoria (Central Venous Catheter kit single lumen)-Polyurethane Nitinol & SS PP SS Needle ABS Non woven fabric PU Film Polycarbonate-Sterile-3- E&D 5.Harsoria (Central Venous Catheter kit double lumen)-Polyurethane Nitinol & SS PP SS Needle ABS Non woven fabric PU Film Polycarbonate-Sterile-3- E&D 6. Harsoria (Central Venous Catheter kit triple lumen)-Polyurethane	CLAA/MD/HR /04/2006-DC (Pt-2)	Mfg Licno No. 637-B (H) Validity from 19/09/2016 to 18/09/2021 (10-11-2017) (Approval for Additional Products)

		Nitinol & SS PP SS Needle ABS Non woven fabric PU Film Polycarbonate-Sterile-3- E&D		
122.	M/s iLife MD PVT.Ltd.,#86,Sector- 56,Phase- IV,I.E.,Kundli,Sonepat- 131028	1.DEN-V-Catheter, Polypropylene,HDPE,LDPE and Stainless Steel(AISI304)-Sterile-5- E &D 2.Becton- Catheter, Polypropylene,HDPE,LDPE and Stainless Steel(AISI304)-Sterile-5- E &D	CLAA/MD/HR /03/2009-DC	Mfg Licno No. 673-B (H) Validity from 13/08/2014 to 12/08/2019 (10-11-2017) (Approval for Additional Products)
123.	M/s Orthotech,Plot No:- 06,N.H.No.8,Nandavala,K ewada,Valsad-396001	1.ECO(Uncemented HA Coated cup/Modular shoulder prosthesi/Elbow Prosthesis/Acetabular Cup-Shell and liner,cup,hinge,rod)- Shell:COBALT Chromim/stainless steel, coated with hydroxyapatite (HA), Liner-UHMWPE/Cross link polymers-Sterile-5- E&D 2.Orthosis(Uncemented HA Coated cup,porous coated cup,HA Coated Elbow prosthesis,HA Coated modular shoulder prosthesis (Cup,hinge,rod,shell and liner)-Shell-Titanium,stainless steel,coated with hydroxyapatite (HA),Liner:-UHMWPE/Cross link polymers- Sterile-5- E&D	CLAA/MD/GU J/01/2013-DC	Mfg Licno No. 673-B (H) Validity from 13/08/2014 to 12/08/2019 (10-11-2017) (Approval for Additional Products)
124.	M/s Appasamy Ocular Devices Pvt.Ltd., R.S.No.9/1,2 & 3,4 &5, NH-45-A,Villupuram Main Road,Vadamangalam- 605102,Puducherry	1.Intra Ocular Lens(PMMA Single piece anterior or posterior chamber intraocular lens-PMMA (Poly Methyl Methacrylate)-The lens is an Ophthalmic replacement of human lens- Sterile 2. Hydrophilic Foldable single pieceintra ocular lens-HEMA (Hydroxy ethyl methacrylate)-The lens is an Ophthalmic replacement of human lens- Sterile 3. Hydrophilic Foldable single pieceintra ocular lens-HEMA (Hydroxy ethyl methacrylate) -Yellow-The lens is an Ophthalmic replacement of human lens- Sterile 4. Hydrophilic Foldable single pieceintra ocular lens- Hydrophobic Acrylic -The lens is an Ophthalmic replacement of human lens- Sterile 5. Hydrophilic Foldable single pieceintra ocular lens- Hydrophobic Acrylic-Yellow-The lens is an Ophthalmic replacement of human lens- Sterile	CLAA/MD/PU D/01/2009-DC	Mfg Licno No. 07 22 2062 Validity from 10/09/2007 to 09/09/2022 (20-11-2017) (Approval for Additional Products)

125.	M/s Surgimedik Healthcare, Surgimedik House, Plot No.223,230,Surat Special Economic Zone, Sachin, Surat-394230	<ol style="list-style-type: none"> 1.Suprapubic Strainage catheter and sets-LatexSilicone elastomer— Sterile-5-E &D 2.Balloon Dilatation Catheter-Rilsan Benso, Polyamide,PTFE- Sterile-5-E &D 3.Instrauterine Insemination catheter-Polyethylene-Sterile -5- E&D 4.Selective Injection Catheter-PEBAX- Sterile -5- E&D 5.Endometrial Biopsy Catheter-Polypropylene- Sterile -5- E&D 6.Dual Lumen Dialysis Catheter and Sets-Polyurethene- Sterile -5- E&D 7.Femoral Catheter- Polyurethene- Sterile -5- E&D 8.Naso Biliary Drainage Catheter-PEBAX and Polyurethene- Sterile -5-E&D 9.Biliary Dilation catheter-PEBAX- Sterile -5- E&D 10.Naso Pancreatic Drainage Catheter- PEBAX- Sterile -5- E&D 11.Injector Needle catheter-PTFE- Sterile -5- E&D 12.Small Bowel Enema Catheters and sets-Polyurethene- Sterile -5-E&D 13.Percutaneous Transp hepatic Cholangiography drainage catheters and sets Polyurethene- Sterile -5- E&D 14.Home Catheter-Polyvinyl chloride- Sterile -5- E&D 15.Spilt-o-cath-Silicone &SS -Sterile -5- E&D 16.Centesis catheter needle-Polyethylene-Sterile-5-E &D 17. Liver Abscess drainage catheter- Polyurethane-Sterile -5- E&D 18.Triple Lumen Catheter- Polyurethane-Sterile -5- E&D 19.Embryo Transfer catheter- Polyurethane & HDPE-Sterile -3- E&D 20.Hystero Salpingography catheter-Polyethlene and polypropylene- Sterile -2- E &D 21. Sonohysterography uterine catheter-Polyethylene and tylon-Sterile - 3- E&D 	CLAA/MD/GU J/09/2009-DC	Mfg Licno No. G/28/1384 Validity from 17/08/2016 to 16/08/2021 (21-11-2017) (Approval for Additional Products)
------	--	--	-------------------------	--

126.	M/s Sahajanand Medical Technologies Pvt.Ltd., 2 &3 Floor, Plot No:- 33,34,35,52,53 &54, Special economic zone, Surat-394230	1.Turquoise +-Dialator (Plastic components),Needle (Stainless Steel and Angiography wire:-Metal core-wire and plastic components)-Sterile-3 – E &D	CLAA/MD/GU J/06/2011-DC (Pt-1)	Mfg Licno No. 07 22 2062 Validity till 28/08/2021 (21-11-2017) (Approval for Additional Products)
127.	M/s ACE Osteomedica, Plot No:-1999,MIE,Pasrt-B,Bahadurgarh,District Jhajjar (Haryana)	<ol style="list-style-type: none"> 1.Austeofix (mini cortex screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 2. Austeofix (Cortex screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 3. Austeofix (Cancellous screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 4. Austeofix (Malleolar screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 5. Austeofix (Washer for cancellous screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 6. Austeofix (Ashnis screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 7. Austeofix (Cannulated Cancellous screw thread)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 8. Austeofix (Large Cannulated Cancellous screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 9. Austeofix (Small Cannulated screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 10. Austeofix (Cannulated screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 11. Austeofix (Shft Screw for LC-DCP/DC)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 12. Austeofix (Steffee Screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 13. Austeofix (Hip Screw with compression screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 14. Austeofix (DHS/DCS Compression screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D 15. Austeofix (Locking bolt)- SS 316L VM,SS,316L,Titanium-Non 	CLAA/MD/HR /08/2010-DC	Mfg Licno No. 707-B (H) Validity upto 17/12/2017 (21-11-2017) (Approval for Additional Products)

Sterile-5- E&D

16. Austeofix (Plug Screw for Interlocking I.M.L.Nails)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

17. Austeofix (End Cap for undreamed humeral nail)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

18. Austeofix (Cancellous traction screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

19. Austeofix (Recon bolt)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

20. Austeofix (Herbert screw all sizes)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

21. Austeofix (ACL Screw All sizes)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

22. Austeofix (Lag Screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

23. Austeofix (Locking head screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

24. Austeofix (Cortexscrew)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

25. Austeofix (Slef Tapping Screws)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

26. Austeofix (Self drilling screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

27. Austeofix (1.5MM/2.0MM,2.5MM,1.7MM, Self drilling screw)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

28. Austeofix (Mini L/T/Straight Plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

29. Austeofix (Dynamic Compression plates)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

30. Austeofix (Straight reconstruction plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

31. Austeofix (Onr third tubular plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

32. Austeofix (Semi tubular plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D

		<p>33. Austeofix (Calcaneal plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>34. Austeofix (Angle blade-plate 130/95 degree with dynamic holes)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>35. Austeofix (Lateral tibial head buttress polate for left,rigt leg)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>36. Austeofix (Spoon PLate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>37. Austeofix (Lengthening plate 4.5MM)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>38. Austeofix (Jewtt nail plate 135 degree/MC Laughlin plate)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>39. Austeofix (Locking DHS Plate 130/135/140/145/150 All sizes/Locking DCS Plate 95 degree all sizes)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>40. Austeofix (Locking anatomical distal tibia plate right,left)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>41. Austeofix (Spacer)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>42. Austeofix (Locking anatomical distal tibia plate right,left)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>43. Austeofix (Universal humeral/femoral nail)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>44. Austeofix (Smith Patterson nail with bolt and washer)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>45. Austeofix (Spare Bolt and washer for S.P.Nail)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>46. Austeofix (Kuntsher cloverlead nail-for femur)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p>		
		<p>47. Austeofix (Smith patterson nail with bolt and washer)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>48. Austeofix (Spare bolt and washer for S.P.Nail)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>49. Austeofix (Cerclage Wire with Loop)- SS 316L</p>		

		<p>VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>50. Austeofix (Drummers Wire)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>51. Austeofix (Centrally threaded Steinman Pin)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p> <p>52. Austeofix (Harrington Rods)- SS 316L VM,SS,316L,Titanium-Non Sterile-5- E&D</p>		
128.	<p>M/s Romsons International,59 J © Noida,Special Economic Zone,Noida,Phase-II,Dadri Road, G.B.Nagar,UP</p>	<p>1.Scalp Vein Set- PVC,MB,PE,SS-Sterile-5- E&D</p> <p>2. Scalp Vein Set (IV Cannula)- FEP,SS,Silicone,PP,POM and MB-Sterile-5- E&D</p> <p>3. Scalp Vein Set(IV Cannula safety cannula with jjection/without injection with/without injection port)-- FEP,SS,Silicone,PP,POM and MB-Sterile-5- E&D</p> <p>4. Scalp Vein Set (Central Venous Catheter)-PU,ASB,PP,PC,SS,PE,MB-Sterile-5- E&D</p> <p>5. Umbilical Catheter-PVC-Sterile-5-E &D</p> <p>6.Suction Catheter(Endobroncheal Suction catheter)-PVC-Sterile-5-E&d</p> <p>7.Laryngeal mask airways PVC-Sterile-5- E&d</p> <p>8.Infant Feeding Tube-PVC,MB-Sterile-5-E &D</p> <p>9.Mucus Extactor-PP,PVC,PE,MB-Sterile-5- E&D</p> <p>10.Oxy Cath-PVC,MB-Sterile-5- E&D</p> <p>11.Foley Balloon Catheter-Latex-Sterile-5- E&D</p> <p>12.Male Cath- Latex-Streile-5- E&D</p> <p>13.Chest Drainage Catheter(GS-5010)-PVC,MB,PP-Sterile-5-E&D</p> <p>14.Turpexe-PP,PE,ABS,PVC,MB-Sterile-5-E&D</p> <p>15.Gude Wire-SS,PE-Sterile-5-E &D</p> <p>16.RomoVac Set-PP,PE,PVC,ABS,SS,MB-Sterile-5-E&D</p> <p>17.Flar Drain Set-PP,PE,PVC, ABS,SS,MB-Sterile-5-E&D</p> <p>18.Tracheostomy Tube-Plain-PP,PVC,PE,Polyamide MB-Sterile-5-E&D</p> <p>19.Gelato-Plian-PVC,MB-Sterile-5-E&D</p> <p>20.EPI Kit-PEBAX,ABS,PC,PE,SS,MB-Sterile-5-E&D</p> <p>21.Romo Seal-K-Resin,PVC,PE,PP,SS,MB,Adhesive Paper-Sterile-5-E&D</p>	<p>CLAA/MD/UP/01/2011-DC</p>	<p>Mfg Licno No. Validity from 22-12-2016 to 21-12-2021</p> <p>(27-11-2017)</p> <p>(Approval for Renewall Products)</p>

		22.Karman Cannula-PP,MB-Sterile-5- E&D		
130.	M/s Disposafe Health and Life Care Ltd.,Plot No:-1 &2,Sector-59,Phase-II,Ballabgarh-121004,Faridabad,Haryana	<p>1.Medeco(Sterile Thoracic Catheter Set)-Polyurethane,Polyvinylchloride,Polycarbonate,Stainless Steel-Sterile-5- E&D</p> <p>2. Medeco(Sterile Suction Catheter)-Polpropylene & Polyvinylchloride, -Sterile-5- E&D</p> <p>3. Medeco(Sterile Foley Balloon Catheter)-Silicon,Latex Rubber,Polypropylene-Sterile-5- E&D</p> <p>4. Medeco(Sterile Oxygen Catheter)- Polpropylene & Polyvinylchloride, -Sterile-5- E&D</p> <p>5.Medeco (Sterile IV Cannula with/without injection port with/without injection stopper)-Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304)-Sterile-5-E &D</p> <p>6. Venfic WP/SW/Leadcare (Sterile IV Cannula with/without injection port with/without injection stopper)-Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304)-Sterile-5-E &D</p> <p>7. E-Care (Sterile IV Cannula with/without injection port with/without injection stopper)-Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304)-Sterile-5-E &D</p> <p>8. Elvacath (Sterile IV Cannula with/without injection port with/without injection stopper)-Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304)-Sterile-5-E &D</p> <p>9. All Pro (Sterile IV Cannula with/without injection port with/without injection stopper)-Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304)-Sterile-5-E &D</p> <p>10.All Pro (Sterile Foley Calloon Catheter)-Silcion,Latex Rubber, Polypropylene-Sterile-5- E&D</p> <p>11. All Pro (Sterile Suction Catheter)- Polpropylene & Polyvinylchloride, -Sterile-5- E&D</p> <p>12.All Pro (Sterile Thoracic Catheter Set)-</p>	CLAA/MD/HR /01/2013-DC	<p>Mfg Licno No. 714-B (H) Validity from 22-12-2016 to 21-12-2021</p> <p>(27-11-2017) (Approval for Additional Products)</p>

		Polyurethane,Polyvinylchloride,Polycarbonate,Stainless Steel-Sterile-5-E&D		
131.	M/s Global Medikit Ltd.,Khasra No:-323 (MI),Camp Road,Selaqui,Dehradun,Uttarakhand	1.Neocan Pro-Body,needle,catheter,holder,threded stopper,needle hub and cover-Sterile-5- E&D 2.Centum Uno/acent-1/2/3,Syringe,Single lumen cathter,clamp,dilator,needle assembly, blade, three way stop cock, advancer,stopper - Sterile-5- E&D	CLAA/MD/UA /01/2009-DC	Mfg Licno No. DML No:- 18/UA/SC/P-2009 Validityupto 23-03-2019 (27-11-2017) (Approval for Additional Products)
132.	M/s Unsisys Ortho, 51,52,53, Bhagirath-4,Industrial Estate, Op..Jawahar Nagar, Nageshwar Estate,Amraiwadi,Ahemed abad	1.Bone Plates-SS/Titanium-Non-Sterile-E &D 2.Bone Screws- SS/Titanium-Non-Sterile-E &D 3.Spinal Implants- SS/Titanium-Non-Sterile-E &D 4.Intramedullary Nils- SS/Titanium-Non-Sterile-E &D 5.Arthroscopy Implants- SS/Titanium-Non-Sterile-E &D 6.Maxilofacial Implants- SS/Titanium-Non-Sterile-E &D	CLAA/MD/GU J/02/2010-DC	Mfg Licno No. G/28/1319 Validity from 18-11-2014 to 17-11-2019 (27-11-2017) (Approval for Renewal Products)
133.	M/s Polymedicure Ltd.,Plot No:-115-116,HSIDC, Sector-59,Ballabgarh,Faridabad-121004,Haryana	1.Single lumen Central Venous Catheter- PVC,ABS,PC,SS,PP-Sterile-5- E&D 2.Double Lumen Central Venous Catheter- PVC,ABS,PC,SS,PP-Sterile-5- E&D 3. Triple Lumen Central Venous Catheter- PVC,ABS,PC,SS,PP-Sterile-5- E&D 4. Quadruple Lumen Central Venous Catheter- PVC,ABS,PC,SS,PP-Sterile-5- E&D 5.Sterile Single lumen central venous catheter (CVC) kit for single use only/Polymed- PVC,ABS,PC,SS,PP-Sterile-5- E&D 6. Sterile Double lumen central venous catheter (CVC) kit for single use only/Polymed - PVC,ABS,PC,SS,PP-Sterile-5- E&D 7. Sterile Triple lumen central venous catheter (CVC) kit for single use only/Polymed - PVC,ABS,PC,SS,PP-Sterile-5- E&D 8. Sterile Quadruple lumen central venous catheter (CVC) kit for single use only/Polymed - PVC,ABS,PC,SS,PP-Sterile-5- E&D	CLAA/MD/HR /01/2006-DC	Mfg Licno No. 627-B (H) Validity from 06-03-2016 to 05-03-2021 (27-11-2017) (Approval for Renewal Products)

		<p>9.I.V.Cannula with/without injection port and wings/with small wings/three ways stop cock(Polycath/Polyflex/Polywin/Ultra/Polyon Adva./Tro-Venocath Sine,novo)-PTFE/FEP/PUR,PE,PP,SS,PA-Sterile-5-Export</p> <p>10.IV Cannula without injection port and with wings (Sentracan Safety)-PTFE/FEP/PUR,PE,PP,SS,PA-Sterile-5- Export</p> <p>11.Heuer (IV Cannula with/without wings with/withouy injection port)-PTFE/FEP/PUR,PE,PP,SS,PA-Sterile-5-Export</p> <p>12.Kilani CVC Kit Central Venous Access Cathert (Single/Double/Triple Lumen Central Venous Access Catheter Kit)-PVC,ABS,PC,SS,PP-Sterile-5-Export</p> <p>13.Foley Balloon Cathter “for single use only-Latex Rubber/Silicon Elastomer- PTFE/FEP/PUR,PE,PP,SS,PA -Sterile-5- E&D</p> <p>14.Sumitex (Catheter sin Puerto de inyeccio y sin alas)-PTFE/FEP/PUR,PE,PP,SS,PA -Sterile-5- E&D</p> <p>15.Idealcare IV Catheter with injection valve (quick flesh back)-PTFE/FEP/PUR,PE,PP,SS,PA -Sterile-5- E&D</p> <p>16.Neo Venopic 1- PTFE/FEP/PUR,PE,PP,SS,PA -Sterile-5- E&D</p>		
134.	M/s Harsoria Helatchare Pvt.Ltd., 110-111,Phase-IV, Udyog Vihar, Gurgaon-122015,Haryana,India	<p>1.JCM-MED (With injection valve and with wings)-Polyproylene HDPE,ABS, Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle) in AISI304 Locking strip in AISI302-Sterile-05-Export</p> <p>2.JCM-MED (Safety IV Cannula catheter and without injectionvalve and with small wings)- Polyproylene HDPE,ABS, Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle) in AISI304 Locking strip in AISI302-Sterile-05-Export</p> <p>3. JCM-MED (Safety IV Cannula with and without injection valave and without wings)-Polyproylene HDPE,ABS, Poly Acetal Catheter (PTFE/PUR/FEP/SS Needle) in AISI304 Locking strip in AISI302-Sterile-05-Export</p>	CLAA/MD/HR /04/2006-DC 9Pt-2)	<p>Mfg Licno No. 637-B (H) Validity from 19-09-2016 to 18-09-2021</p> <p>(19-12-2017) (Approval for Additional Products)</p>
135.	M/s Uma Surgicals,21/22/27,Gayatri Industrial Estate, Village-Navghar,Vasai (E),Maharashtra-401210	<p>1.Bone plates (638874)-Non Sterile- Domestic</p> <p>2. BoneScrews -Non Sterile- Domestic</p> <p>3.Endo Prosthesis Replacement (638872)- Non Sterile- Domestic</p> <p>4.Intra Medullary Nails (638875)- Non Sterile- Domestic</p> <p>5.Pins and Wires (638873)- Non Sterile- Domestic</p> <p>6.Spine Implants - Non Sterile- Domestic</p>	CLAA/MD/MS /01/2010-DC	<p>Mfg Licno No. KD/491 Validity from 19-06-2010 to 09-06-2015</p> <p>(19-12-2017) (Approval for Additional</p>

				Products)
136.	M/s Hollister India Pvt.Ltd.,Plot No:- 27,Sector-5,Phase-II,HSIDC Growth Center, Bawal,Rewari-123501,Haryana	1.Intermittent Catheter without blister pack-Medical grade PVC Tube-Medical grade, Gamma Resistant,Polyvinyl Chloride Tri Trimellitate (TEHTM Palsticizer) compound, Funnel-A Blend of 96% raumedic noDOP-Compopund in plasticizedpolyvinyl chloride (PVC)-Un-Sterile-Export	CLAA/MD/HR /04/2009-DC	Mfg Licno No. 676-B (H) Validity till 11-11-2019 (19-12-2017) (Approval for Additional Products)
137.	M/s Envision Scientific Pvt.Ltd.,Plot No:-C1B-108,Ichchhapore GIDC, Industrial estate, Bhatpore,Sachin-Hazira Road,Surat-394510,Gujarat,India	1.Sirolimus eluting coronary stnt system (Crocus ninth brand)-L605 Co-Cr Alloy ,Polyamide Rilsan-Sterile-2-Export 2.Sirolimus Based bifurcated stent system (DESA) (Bifur)- L605 Co-Cr Alloy ,Polyamide Rilsan,Slip coat Hydrophilic coating,Polyamide grilamid,Polyamide rilsan,blue pigment,Stainless steel 304L,Polyamide Rilsan,Trilayer tubing-Sterile-2-E &D	CLAA/MD/GU J/13/2011-DC	Mfg Licno No. G/28/1396 Validity from 19-12-2011 to 18-12-2021 (20-12-2017) (Approval for Additional Products)
138.	M/s Global Medikit Ltd., Kasra No:-323(MI), Camp Road, Selaqui,Dehradun,Uttarkhand	1.IV Catheter with injectionvalve and wings (Minut)-Body,needle, catheter, silicon, tube, port cap,holder,threded stopper,needle hub and cover-Sterile-5-Export 2.IV Catheter with wings for neonates- Body,needle, catheter, holder,threded stopper,needle hub and cover-Sterile-5-Export 3.IV Cannula with wings and injection port- Body,needle, catheter, silicon, tube, port cap,holder,threded stopper,needle hub and cover-Sterile-5-Export	CLAA/MD/UA /01/2009-DC	Mfg Licno No. 18/UA/SC/P-2009 Validity till 23-03-2019 (20-12-2017) (Approval for Additional Products)
139.	M/s Meril Life Sciences Pvt.Ltd.,Survey No:- 135/2a and 135/2/B,174/2(M1-M2,meril Park,4 th floor),Near G.M.Bilkhia Stadium, Muktanand MArg, CHala,Vapi-396191,Valsad,Gujarat	1.Pacilitaxerl eluting PTA Ballonn catheter(Paximus)-Sterile-2- E &D	CLAA/MD/GU J/01/2008-DC (Pt-IV)	Mfg Licno No. G/28/1273 Validity from 20-06-2013 to 19-06-2018 (20-12-2017) (Approval for Additional Products)

140.	M/s Nipro India Corporation Pvt.Ltd., Plot No:-E1,MIDC, Kesurdi Khaandala,Tal:-Khandala,Dist:-Satara,Maharashtra-412802	1.Nipro wing cath/Sterile IV Cannula (IVC) with injection port and wing-Polypropylene,Ethylene Tetra Fluoro Ethylene/Polyurethane,Stainless Steel,Silicone,Poly Butylene Terephthalate ISO (-)-Export 2. Nipro wing cath/Sterile IV Cannula (IVC) with injection port and wing-Polypropylene,Ethylene Tetra Fluoro Ethylene/Polyurethane,Stainless Steel,Silicone,Poly Butylene Terephthalate ISO (-)--Domestic	CLAA/MD/MS /04/2012-DC	Mfg Licno No.PD/128 Validity from 03-07-2017 to 02-07-2022 (20-12-2017) (Approval for Renewal Products)
141.	M/s Hindustan Syringes and MD, Ltd.,174,Sector-25,Ballabgarh,Faridabad-121004,Haryana	1.Safety Intravenous Cannula with/without Injection port (Clip wing ed and ported)-Polypropylene ,Satinless Steel(SUS 304 or equivalent)-Sterile and Non Sterile-5-Export 2.Safety Intravenous Cannula with/without injection port with SIP CLIP for sharp injury prevention(Cathy and Cathula)- Polypropylene ,Satinless Steel(SUS 304 or equivalent)-Sterile and Non Sterile-5-Export 3.Intravenous Cannula with and without injection port(Kitkath)- Polypropylene ,Satinless Steel(SUS 304 or equivalent)-Sterile and Non Sterile-5-Export 4.Intravenous Cannula with/without injection port(Cathy and Cathula)- Polypropylene ,Satinless Steel(SUS 304 or equivalent)-Sterile and Non Sterile-5-Export 5.Intravenous Cannula (with/without injection port and with and without safety clip)- Polypropylene ,Satinless Steel(SUS 304 or equivalent)-Sterile and Non Sterile-5-Export	CLAA/MD/HR /01/2014-DC	(20-12-2017) (Approval for Additional Products)
142.	M/s Prymax Healthcare LLP,Plot NO:-53/17,Industrial Area,NIT Faridabad-121001,Haryana	1.Catheter(Prymax dialysis hemodialysis catheter kit /Prymax ANgi Introducer sheath kit)-PU,Silicone,PC,SS304,PP,Nitinol/PTFE coated-Sterile-3-E &D	CLAA/MD/HR /01/2015-DC	Mfg Licno No.757-B (H) Validity from 01-02-2016 to 31-01-2021 (20-12-2017) (Approval for Additional Products)
143.	M/s Becton Dickinson India Pvt.Ltd.,Plot No:-1,Sector-3,IMT,Bawal,Rewari (Haryana)	1.Venflon (IV Cannula)-Polypropylene,Polyethylene, Methyl vinyl polysiloxane,polytetrafluoroethylene,polycarbonate white,stainless steel,polydimethylsiloxane,Tyvek Medicalgrade paper, PVC Blister film APET,Duplex boxboard, Corrugated box-Sterile-5- E&D	CLAA/MD/HR /01/2009-DC	Mfg Licno No.671-B (H) Validity from 24-03-2014 to 23-03-2019 (20-12-2017)

				(Approval for Amendment)
144.	M/s Mediplus India Ltd.,1261,1262,mie, Bahdurgarh,Dist:- Jhajjar,Haryana-124507	1.IV Cannula with injectionvalve with wings (Vitasta)- Polyproylene,Polyethylene,POLyacetal,Stainless Steel (AISI304)-Sterile-5- Export	CLAA/MD/HR /06/2006-DC	Mfg Licno No.646-B (H) Validity from 01-06-2017 to 31-05-2022 (20-12-2017) (Approval for Additional Products)
145.	M/s Orthocare and Cure 9I) Pvt.Ltd.,Plot No:- 22,S.No.327,B/1,Opp Vedh textile,At Post Mudshingi,Tal:- akrvir,Kolhapur-416119	1.Hip prosthesis-SS 316L-Non Sterile-2- E &D 2.Spinal Implant- SS 316L/Titanium-Non Sterile-2- E &D 3.	CLAA/MD/MS /02/2017-DC	Mfg Licno No.MH 102033 Validity from 09-01-2017 to 08-01-2022 (20-12-2017) (Approval for products)
146.	M/s Geemed Healthcare,A-4/3/28,South Side Industrial Area,NH-24,Ghaziabad-226001,Uttar Pradesh	1. E.T.Tube plain and cuffed (Geemed)-Manfctured from silcionised non toxic, non irritant,kink resistant PVC material-Sterile-5- E&D 2.Guedel Airways(Geemed)-Manufactured from EVA t oac hieve extra soft surface finish- Sterile-5- E&D 3.Clave Connectors(Geemed)-Clave connectors are made upto sicion material- Sterile-5- E&D 4. Connectors(Geemed)-Extracorporeal medical tubing connectors are maniufatcured from crystal clean polycarbonate USP class IV in plant tested with GE Latex 144 grade- Sterile-5- E&D 5.Chest Drainage cathter (Geemed)-Made from toxic medical grade PVC Meterial- sicion material- Sterile-5- E&D 6.Chest drainage cathter angeles(Geemed)- Made from toxic medical grade PVCMaterial- Sterile-5- E&D 7.Chest drainage catheter with trocar- Made from toxic medical grade PVCMaterial like box made up t poly carbonates- Sterile-5- E&D 8.Rigid Suction devices(Sucker) (Geemed)- The suction devices consist of special stainless steel (grade SS 404) tube with a ribbed mini trip and	CLAA/MD/UP/ 03/2016-DC	Mfg Licno No./SC/P of 2017 Validity from 09-01-2017 to 08-01-2022 (20-12-2017) (Approval for products)

		<p>ergonomically designed handle-Sterile-5- E &D</p> <p>9.Hemoconcentrator(Geemed)-Made up of medical grade PVC polysulfide membrane offers excellent biocompatibility-Sterile-5- E&D</p> <p>10.A.V.Tubing (Geemed)-made form medical PVC Material vented protective caps are also made of soft PVC-Sterile -5- E&D</p> <p>11.Disposable Chesr Driange system (C D Bottle) (Geemed)-made up of kink resistant non toxic PVC-Sterile-5- E&D</p> <p>12.MTP Cannula(Geemed)-manfuatcured from soft medical grade LDPT-Sterile-5- E&D</p> <p>13.Umbilical catheter(Geemed)-manfuatcured from soft medical grade PVC-Sterile-5-E &D</p> <p>14.Yankeur suction set(Geemed) – made from non toxic medical grade PVC material-sterile-5- E&D</p> <p>15.Abdominal drainage set (Geemed)-made from non toxic medical VC Material and medical grade PVC Sheet-Sterile-5-E&D</p> <p>16.Corrugated drain sheet(Geemed)-manfauctured from non toxic PVC sheet-Sterile-5- E&d</p> <p>17.Venous cannula straight(Geemed)-Made up of non toxic medical grade PVC Material connetcormade of PC material SS Wire-Sterile-5-E&D</p> <p>18.Foley Catheter 2way and 3way(Geemed)--manufactured with latex-Sterile-5-E&D</p> <p>19.Suction catheter plain and with thumb control(Geemed)--Xatheter are manfauctured from non toxic non irritants soft PVC Compatible with catheter-Sterile-5-E &d</p> <p>20.Ryle`S tube(Geemed)--manfuactured from non toxic PVC-Sterile-5-E&d</p> <p>21.Levin`s ube(Geemed)--manfuatcured from low friction medical grade PVC-Sterile-5-E&D</p> <p>22.Infant feeding tube(Geemed)-mafnauctured from non toxic PVC-Sterile-5- E&D</p> <p>23.Mucus extractor (Geemed)-manfauctured from non toxic medical grade PVC-Sterile-5- E&D</p> <p>24.Peritoneal dialysis set and catheter(Geemed)-manufactured of stainless spring-Sterile-5- E&D</p>		
--	--	---	--	--

		<p>25. Guide Wire (Geemed)-manufactured of stainless spring steel wires- Sterile-5-E&D</p> <p>26. Hemodialysis catheter kit (Geemed)--Hemolysis catheter is polyurethane catheter radio opaque with length markers, color coded safety mental clamps, lure connection plastic hubs and PP protective caps scapel, lurelock PP/PE Syringes, introducer needle, guide wire- Sterile-5-E&D</p> <p>27. Femoral catheter kit (Geemed)-Made up of soft medical PVC and PP- Sterile-5-E&D</p> <p>28. Hemodialysis blood tubing set (Geemed)--consists of two parts arterial and venous used during dialysis and attached with fistula and dialyzer- Sterile-5-E&D</p> <p>29. TUR Set (Geemed)-manufactured from non toxic medical grade PVC- Sterile-5-E&D</p> <p>30. Nelaton Catheter (Geemed)--Manufactured from non toxic non irritant, medical grade PVC- Sterile-5- E&D</p> <p>31. Urethral catheter (Geemed)--manufactured from non toxic, non irritant, medical grade- Sterile-5-E&D</p> <p>32. Male catheter (Geemed)-made up of latex- Sterile-5-E &D</p> <p>33. Suprubic balloon catheter (Geemed)-made from non toxic medical grade PVC- Sterile-5-E&D</p> <p>34. Central venous catheter (CVC) KIT (Geemed)-polyurethane catheter with fixation on device (single lumen, double lumen, triple lumen)- 1 item, introducer needle, J-Trip guide wire, vessel dilator, additional clamp for attachment on patient`s skin, connector with lure lock system (needleless connector) for needle less syringe attachments- Sterile-5-E &D</p> <p>35. IV Cannula with wings and with port, without wings and without port, with wings without port, without wings with post (Geemed)-manufactured of PP caps, PTFE/FEP/PU Tubes, SS Needles. Silcion rubber, LDPE port cap, PP hub, PP Wing body and POM Teflon holder- Sterile-5- E&D</p>		
--	--	--	--	--

147.	M/s Meril Healthcare Pvt.Ltd., Survey No:- 135/2/b AND 174/2,1 ST FLOOR, h1-h3, Meril Park, Muktanand Marg,Chala,Vapi-396191	1.Proximally coated uncemented stem-Ti-6Al-4V ELI, CPTI Coating,CPTI+HA Coating,HA Coating-To be sterilized by gamma radiation-5- E&D 2. Proximally coated uncemented stem-Ti-6Al-4V ELI, CPTI Coating,CPTI+HA Coating,HA Coating-Non Sterile-5- E&D	CLAA/MD/GU J/05/2013-DC	Mfg Licno No.G/28/1465 Validity from 01-02-2014 to 31-01-2019 (21-12-2017) (Approval for Additional products)
148.	M/s Unisys Ortho ,51,52,53,Bhagirath-4,industrial estate, Opp.JawaharnagarmNages hwar Estate,Amraiwadi,Ahmedabad-380026	1.Sternum locking plates-Titanium-Non Sterile-E &D 2.2.4mm Cross drive lock cancellous screw- Titanium-Non Sterile-E &D 3.2.7mm Cross drive lock cancellous screw- Titanium-Non Sterile-E &D 4.Cranial Bone plate- Titanium-Non Sterile-E &D 5.Cranial orbital plate- Titanium-Non Sterile-E &D 6.Cranial cross drive bone screw- Titanium-Non Sterile-E &D 7.Maxilifical bone plate- Titanium-Non Sterile-E &D 8.Reconstruction bone plate- Titanium-Non Sterile-E &D 9. Maxilifical cross drive bone screw- Titanium-Non Sterile-E &D	CLAA/MD/GU J/02/2010-DC	Mfg Licno No.G/28/1319 Validity till 17-11-2019 (22-12-2017) (Approval for Additional products)
149.	M/s Pragati Surgitech ,#49,50.Annapurna industrial estate, inside Zaveri Estate/H.C.L.High School, Singarva Katheada Rpad,Ahmedabad-382430	1.Cortex Screw-SS316L,3316LVM,Titanium-Non Sterile-E&D 2.Cancellous SCREW- SS316L,3316LVM,Titanium-Non Sterile-E&D 3.Cannulated cancellous screw- SS316L,3316LVM,Titanium-Non Sterile-E&D 4.Locking head Screw- SS316L,3316LVM,Titanium-Non Sterile-E&D 5.Locking Bolt- SS316L,3316LVM,Titanium-Non Sterile-E&D 6.Cortical Screw- SS316L,3316LVM,Titanium-Non Sterile-E&D 7.Cancellous locking screw- SS316L,3316LVM,Titanium-Non Sterile-E&D 8.Cannulated Cancellous Screw (16 th , 32 th)- SS316L,3316LVM,Titanium-Non Sterile-E&D	CLAA/MD/GU J/11/2017-DC	Mfg Licno No.G/28/1638 Validity from 10-10-2017 to 09-10-2022 (22-12-2017) (Approval for Products)

150.	M/s Parshwa Corporation ,Plot No:-2,Survey No:- 133,Village ODE, Taluka- Daskroi,Dist:-Ahmedabad- 382427	1.Cortex Screw- 3316LVM,Titanium-Non Sterile-E&D 2.Cancellous Screw- 3316LVM,Titanium-Non Sterile-E&D 3.Cannulated Cancellous Screw- SS316L,3316LVM,Titanium-Non Sterile-E&D	CLAA/MD/GU J/08/2017-DC	Mfg Licno No.G/28/1644 Validity from 13-11-2017 to 12- 11-2022 (22-12-2017) (Approval for Products)
151.	M/s Osseous Surgicals,Plot No:-3,amruti Industrial Estate, Near Uma cotton,hadmatiya road, Lajai,Morbi Rajkot Highway,Tal;-Tankara, Dist:-Morbi-363650	1.Cortex Screw- 3316LVM,Titanium-Non Sterile-E&D 2.Cancellous Screw-3316LVM,Titanium-Non Sterile-E&D 3. Cancellous Screw Cannulated Cancellous Screw- SS316L,3316LVM,Titanium-Non Sterile-E&D	CLAA/MD/GU J/13/2017-DC	Mfg Licno No.G/28/1648 Validity from 28-11-2017 to 27- 11-2022 (22-12-2017) (Approval for Products)
152.	M/s Freedom Ophthalmic Pvt.ltd.,Plot No:-31,phase- 1,SIPCOT Industrial Complex, SIPCOT,Hosur- Krishnagiri Dist:- TamiNadu-635126	1.PMMA Posterior chamber single piece intraocular lens-POlyMethyl Methacrylate with UV Blocker-Sterile-5-E &D 2. PMMA Posterior chamber single piece square edge intraocular lens- POlyMethyl Methacrylate with UV Blocker-Sterile-5-E &D 3.IRIS Claw intraocular lens-PMMA- PMMA Posterior chamber single piece intraocular lens-POlyMethyl Methacrylate with UV Blocker- Sterile-5-E &D 4.Hydrophilic acrylic posterior chamber single piece foldable intraocular lens-Copolymer of hydroxyethyl methacrylate and PMMA Posterior chamber single piece intraocular lens-Methyl Methacrylate with UV Blocker-Sterile-3-E &D 5.Hydrophilic acrylic posterior chamber single piece square edge aspheric foldable intraocular lens-Copolymer of hydroxyethyl methacrylate and PMMA Posterior chamber single piece intraocular lens-Methyl Methacrylate with UV Blocker-Sterile-3-E &D 6.Hydrophilic acrylic posterior chamber single piece square edge aspheric foldable intraocular lens-Copolymer of hydroxyethyl methacrylate and PMMA Posterior chamber single piece 360 degree square edge intraocular lens-Methyl Methacrylate with UV Blocker- Sterile-3-E &D 7.Yellow hydrophilic acrylic posterior chamber single piece square aspheric foldable intraocular lens-Yellow copolymer of hydroxythly	CLAA/MD/TN /02/2007-DC	Mfg Licno No.TN00002315 Validity from 13-04-2017 to 12- 04-2022 (22-12-2017) (Approval for Renewal Products)

		<p>methacrylate and methyl methacrylate and UV Blocker-Sterile-3-E &D</p> <p>8.Hydrophobic Intraocular lenses-Cross linked hydrophobic acrylic co-polymer-Sterile-5- E&D</p> <p>9.Yellow Hydrophobic acrylic aspheric foldable intraocular lenses- Yellow cross linked hydrophobic acrylic Co-Polymer-Sterile-5-E &D</p> <p>10.Yellow hydrophobic acrylic aspheric foldable 360 degree square intraocular lenses- Yellow cross linked hydrophobic acrylic co-polymer-Sterile-5-E &D</p> <p>11.Hydrophilic acrylic posterior chamber multifocal diffractive – refractive aspheric square edge single piece foldable intraocular lens- Copolymer of hydroxyethyl methacrylate and methyl methacrylate with UV Blocker-Sterile-3- E&D</p> <p>12.Yellow hydrophilic acrylic posterior chamber multifocal diffractive-refractive aspheric square edge single piece foldable intraocular lens- Yellow copolymer of hydroxyethyl methacrylate and methyl Methacrylate with UV Blocker-Sterile-3-E &D</p> <p>13.hydrophobic acrylic posterior chamber multifocal diffractive refractive aspheric square edge single piece foldable intraocular lens- Cross linked hydrophobic acrylic Co-Polymer-Sterile-5-E&D</p> <p>14.Yellow hydrophobic acrylic posterior chamber toric aspheric square edge piece foldable intraocular lens-yellow cross linked hydrophobic acrylic co-polymer-Sterile-5- E&D</p> <p>Hydrophobic acrylic posterior chamber multifocal toric aspheric square edge single piece foldable intraocular lens-Cross linked hydrophobic acrylic co-polymer-Sterile-5-E &D</p>		
153.	<p>M/s Shree Hari Surgical, Plot No:-3, Survey No:-66/4, Gurudev Industrial Area, Nr Leyras Backry, Ravki Road, Tal:- Lodhika, Dist:-Rajkot-360004, Gujarat</p>	<p>1.Cortex Screw- 3316LVM, Titanium-Non Sterile-E&D</p> <p>2.Locking Screw-316LVM, Titanium-Non Sterile-E &D</p> <p>3.Maleolar Screw-316LVM, Titanium-Non Sterile-E &D</p> <p>4.Cancellous Screw-3316LVM, Titanium-Non Sterile-E&D</p> <p>5.locking Bolt-3316LVM, Titanium-Non Sterile-E&D</p> <p>6.Dhs Screw-3316LVM, Titanium-Non Sterile-E&D</p> <p>7.p f n Screw-3316LVM, Titanium-Non Sterile-E&D</p> <p>8.herbert Screw-3316LVM, Titanium-Non Sterile-E&D</p>	<p>CLAA/MD/GU J/07/2016-DC</p>	<p>Mfg Licno No.G/28/1647</p> <p>Validity from 20-11-2017 to 19-11-2022</p> <p>(26-12-2017)</p> <p>(Approval for Products)</p>

154.	M/s Translumina Therapeutics, Plot No:- 9, Pharmacity, Selaqui, Dehradun-248197, Uttarkhand	1. Sirolimus Drug eluting coronary stent system (Vivo PE)-Stainless Steel (SS-316 LVM)-Sterile (by EO)-24 months –E &D 2. Sirolimus Drug eluting coronary stent system (Vivo ISAR)-L 605 Cobalt Chromium alloy- Sterile (by EO)--24 months –E &D	CLAA/MD/UA /01/2010-DC	Mfg Licno No.12/UA/SC/P-2016 Validity from 05-02-2016 to 04-02-2021 (26-12-2017) (Approval for Additional Products)
155.	M/s Madras Medical Devices SF, No:-150, Site NO:-063, Salasundaram Nagar, K.R.Puram, Coimbatore-641006, Tamil Nadu	1. Peritoneal dialysis catheter-Silastic (Silicone rubber, non coated), material radiopaque medical grade tubing – Sterile-2- E&D	CLAA/MD/TN /03/2017-DC	Mfg Licno No.TN00004758 Validity from (26-12-2017) (Approval for Products)
156.	M/s Rivarp Medical Pvt.Ltd., #34, Azeez Sait Industrial Town, 6 th Mile, Mysore Road, Nayandahalli, Bangalore-560039	1. Central Venous Catheterization Set (Single Lumen, Double Lumen, Triple Lumen, Quadra Lumen)-Medical grade polyurethane-Sterile (Ethylene oxide gas)-3 E&D 2. Hemodialysis catheterization Set (Straight, curved and M Shape)-Medical grade polyurethane-Sterile (Ethylene oxide gas)-3 E&D 3. Double J Shape Catheter/ureteral stent (open end, close end, multi loop, with sutures, set)- Medical grade polyurethane-Sterile (Ethylene oxide gas)-3 E&D 4. IUI Catheter Cannula (Straight with syringe, straight without syringes and curved with syringe, curved without syringe)-Medical grade Polyethylene-Sterile (Ethylene Oxide Gas)-3- E&D	CLAA/MD/KA R/03/2017-DC	Mfg Licno No.KTK/28/452/2017 Validity from 07-12-2017 to 06-12-2022 (26-12-2017) (Approval for Products)
157.	M/s Sharma Pharmaceuticals Pvt.Ltd., 444, GIDC Estate, Waghodia-391760, Gujarat, India	1. Locked plates- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&D 2. Locked Screws- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&d 3. Non Locked screws- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&d 4. Interlocking nails and accessories- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&d 5. Spinal Implants- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&d	CLAA/MD/GU J/03/2008-DC	Mfg Licno No.G/28/1302 Validity till 03-05-2019 (27-12-2017) (Approval for Additional Products)

		<p>6. Joints and Prosthesis- SS316L, SS316LVM, Titanium Alloy, CoCr Alloy-Non Sterile-E&d</p> <p>7. Tibial Inserts (fixed bearing)(PS)(Cemented and UN-cemented) total knee replacement (Joints and Prosthesis- UHMWPE, XLPE-Non Sterile-E&d)</p> <p>8. Tibial inserts (fixed bearing)(CR)(Cemented and uncemented) total knee replacement(Joints and Prosthesis)- UHMWPE, XLPE-Non Sterile-E&d)</p> <p>9. Patella Components (cemented and uncemented) total knee replacement-(Joints and Prosthesis)- UHMWPE, XLPE-Non Sterile-E&d)</p> <p>10. Tibial inserts mobile bearing (mobile bearing) total knee replacement-(Joints and Prosthesis)- UHMWPE, XLPE-Non Sterile-E&d)</p> <p>11. Tibial inserts hinged knee (cemented and uncemented) total knee replacement-(Joints and Prosthesis)- UHMWPE, XLPE-Non Sterile-E&d)</p>		
158.	M/s Uteshiya Medicare Pvt.Ltd., 7F1, F2 Dashrathbh ai Estate, B/h Part Plot, Amaraiwadi, ahmedabad-380026, Gujarat, India	<p>1. Button with loop-SS316L, SS316LVM, Titanium, Peek-Sterile-3- E&D</p> <p>2. All sutures Anchor-UHMWPE-Sterile-3-E &D</p> <p>3. Anchor-SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-Ha-Sterile-3-E &D</p> <p>4. Preloaded suture anchor-SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-Ha-Sterile-3-E &D</p> <p>5. Sutures Post-SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-A-HA, PLLA-BitaTCP- Sterile-3-E &D</p> <p>6. Sutures Disc- SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-A-HA, PLLA-BitaTCP- Sterile-3-E &D</p> <p>7. Spiked Washer- SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-A-HA, PLLA-BitaTCP- Sterile-3-E &D</p> <p>8. Spiked ligament staples-SS316L, SS316LVM, Titanium, Peek, PEKK, Plla, Plla-A-HA, PLLA-BitaTCP- Sterile-3-E &D</p> <p>9. Tibia Plate-SS316L, SS316LVM and titanium-Sterile-3- E&D</p> <p>10. ACL Screw-ss316l, ss316lvm, Titanium, PEEK, PEKK, Plla, plla-HA, Plla-Bita TCP-Sterile-3-E &D</p>	CLAA/MD/GU J/07/2013-DC	<p>Mfg Licno No.G/28/1464</p> <p>Validity till from 16-01-2014 to 15-01-2019</p> <p>(27-12-2017)</p> <p>(Approval for Additional Products)</p>

		<p>10. Austin Moore Prosthesis-SS316L,SS316LVM and Titanium-Sterile-3-E &D</p> <p>11. Bipolar- SS316L,SS316LVM, UHMWPE and Titanium-Sterile-3-E &D</p> <p>12. Modular Bipolar- SS316L,SS316LVM, UHMWPE and Titanium-Sterile-3-E &D</p> <p>13. Total Hip replacement prosthesis- SS316L,SS316LVM, UHMWPE and Titanium-Sterile-3-E &D</p> <p>14. Total knee replacement prosthesis- SS316L,SS316LVM, UHMWPE,Cobalt Chrome and Titanium-Sterile-3-E &D</p> <p>15. Shoulder Prosthesis- SS316L,SS316LVM, Titanium-Sterile-3-E &D</p> <p>16. Cervical Spacer-SS316,SS316LVM, Totanium, Peek,pEkk,Plla,Plla-HA,Plla-Bita TCP-Sterile-3-E&D</p> <p>17. TLIF Spacer- SS316,SS316LVM, Totanium, Peek,pEkk,Plla,Plla-HA,Plla-Bita TCP-Sterile-3-E&D</p> <p>18. PLIF Spacer- SS316,SS316LVM, Totanium, Peek,pEkk,Plla,Plla-HA,Plla-Bita TCP-Sterile-3-E&D</p> <p>19. Thomson Prosthesis-SS316L,SS316LVM, Titanium-Sterile-3-E&D</p>		
159.	M/s BioMed Healthcare Products Pvt.Ltd.,Plot No:-30,DLF Industrial Estate, Faridabad-121003, Haryana	<p>1. Bio-Flon-Polypropylene, Polyethylene, SS Needle, Catheter, Polyacetal, Silcion Tube-Sterile-5-E &D</p> <p>2. Bio-Can- Polypropylene, Polyethylene, SS Needle, Catheter, Polyacetal, -Sterile-5-E &D</p> <p>3. Bio-On- Polypropylene, Polyethylene, SS Needle, Catheter, POLYacetal -Sterile-5-E &D</p> <p>4. Neo-Kath-Polyporylene, Polyethylene, SS Needle, Catheter, PolyAcetal-Sterile-5- E&D</p> <p>5. Solidor- Polyporylene, Polyethylene, SS Needle, Catheter, PolyAcetal-Sterile-5- E&D</p> <p>6. Bel-Flon- Polypropylene, Polyethylene, SS Needle, Catheter, Polyacetal, Silcion Tube-Sterile-5-E &D</p> <p>7. Maxima- Polypropylene, Polyethylene, SS Needle, Catheter, Polyacetal, Silcion Tube-Sterile-5-E &D</p> <p>8. Euro-Flon- Polypropylene, Polyethylene, SS Needle, Catheter, Polyacetal, Silcion Tube-Sterile-5-Export</p> <p>9. Euro-Can-Polypropylene, Polyethylene, SS Needle, Catheter, Poly</p>	CLAA/MD/HR /01/2007-DC	<p>Mfg Licno No.651-B (H)</p> <p>Validity till from 25-07-2017 to 24-07-2022</p> <p>(27-12-2017)</p> <p>(Approval Products)</p>

		<p>Acetal-Sterile-5- Export</p> <p>10.Venax-2-Polyprylene,Polyethylene,SS Needle,CDatheter,PolyAcetal Silicon Tube-Sterile-5-Export</p> <p>11. .Venax-1Polyprylene,Polyethylene,SS Needle,CDatheter,PolyAcetal Silicon Tube-Sterile-5-Export</p> <p>12.Gren-Flon-*Polypropylene,Polyethylene,SS Needle,Catheter,Polyacetal,Silicon Tube-Sterile-5-Export</p> <p>13.Comcaset P-Polypropylene,Polyethylene,SS Needle, Catheter, Polyacetal-Sterile-5- E&D</p> <p>14.Vagon Plus-Polypropylene,Polyethylene,SS Needle, Catheter, Poly acetal, SiliconeTube-Sterile-5-Export</p> <p>15.Comcaset-Polypropylene,Polyethylene,SS Needle, Catheter, Poly acetal-Sterile-5-E &D</p> <p>16.NovaFlon- Polypropylene,Polyethylene,SS Needle, Catheter, Poly acetal, Silicon Tube-Sterile-5-Export</p> <p>17.Net-Can-Polypropylene,Polyethylene,SS Needle, Catheter,poly acetal-Sterile-5- Export</p> <p>18.Swan Medicare- Polypropylene,Polyethylene,SS Needle, Catheter,poly acetal-Sterile-5- Export</p> <p>19.Romed-Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Cilicon Tube-Sterile-5-Export</p> <p>20.Kay-Flon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal Silicon Tube-Sterile-5-Export</p> <p>21.Holden- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal, Silicon Tube-Sterile-5-Export</p> <p>22.Morningside- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal, Silicon Tube-Sterile-5-Export</p> <p>23.Vasocath- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal, Silicon Tube-Sterile-5-Export</p> <p>24.Demotek- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>25.Meeco- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>26.Farcocath- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p>		
--	--	--	--	--

		<p>27.HM Healthcare- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>28. HM Healthcare- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,-Sterile-5-Export</p> <p>29.Medicare- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>30.Irismed- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>31.Re-Flon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>32.Cannatech- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>33.Agary- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>34.Ambi-Flon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>35.Ambi-Can- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>36.Ambi-kath- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>37.Indo-Flon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>38.IndoPLas- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>39.Ivecath- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>40.Kkaremax- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>41.Medin- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>42.Bex- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>43.Helmflon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>44. HelmflonPolypropylene,Polyethylene,SS Needle, Catheter,Poly</p>		
--	--	--	--	--

		<p>Acetal -Sterile-5-Export</p> <p>45.Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>46. Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>47.iGLACATH- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>48.Iglacath- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>49.VM Flon- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>50.Medicann- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>51.Medicann- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>52.Sidra- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>53.Jrz- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>54.JRZ- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p> <p>55.Harmed- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>56.Asco- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal,Silicon Tube-Sterile-5-Export</p> <p>57.Asco- Helmsyte- Polypropylene,Polyethylene,SS Needle, Catheter,Poly Acetal -Sterile-5-Export</p>		
160.	M/s Disposafe Health and Life Care Ltd., Plot No:-1 & 2,Sector-59,Phase-II,Ballabgarh-121004,Faridabad	<p>1.Cathtec (Sterile scalp vein set)- Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304),Polyvinyl,Chloride ABS-Sterile-5-E &D</p> <p>2. Cathtec (Sterile foley balloon catheter)- Polypropylene,Polyethylene,Silicon,Latex Rubber- Sterile-5-E &D</p> <p>3. Cathtec (Sterile sunction catheter)-Polypropylene, Polyvinyl Chloride -Sterile-5-E &D</p>	CLAA/MD/HR /01/2013-DC	<p>Mfg Licno No.714-B (H) Validity till from 25-07-2013 to 26.05.2018</p> <p>(27-12-2017) (Approval Additional Products)</p>

161.		<p>4.Leadcare (Sterile Scalp vein set)- Polypropylene,Polyethylene,Polyacetal,Stainless Steel (AISI304),Polyvinyl,Chloride ABS-Sterile-5-E &D</p> <p>5.Leadcare (Sterile Foley Balloon catheter)- Polypropylene,Polyethylene,Silicon,Latex Rubber- Sterile-5-E &D</p> <p>6.Leadcare (Sterile Suction Catheter)- Polypropylene, Polyvinyl Chloride -Sterile-5-E &D</p> <p>7.Cathtec(IV Cannula with integrated three way stopcock)- Polypropylene,Polyethylene,Polyacetal,Stainless Steel(AISI304),Polycarbonate,Telfon,PUR,PTFE,FEP-Sterile-5-E &D</p>		
162.	M/s Shajanand Laser Technology Ltd.,A-8,GIDC Electronic Estate,Sector-25,GandhiNagar-382016	1.Cobalt Chromium Everolimus eluting coronary stent system –Cobalt chromium-Sterile and Non Sterile-24 months-E &D	CLAA/MD/GU J/08/2009-DC	<p>Mfg Licno No.G/28/1314 Validity till from 27-08-2014 to 26-08-2019</p> <p>(27-12-2017) (Approval Additional Products)</p>
163.	M/s NRV Orthotech Pvt.Ltd.,Plot No:-18,Aneri Industrial Park,b/h,Super Gas,Bavla Bagodra Road,Ranesar,Ahmedabad-382240,Gujarat	<p>1.N-Fix Cortex Screw-SS316L,SS316LVM and titanium-Sterile and Non Sterile-3-E &D</p> <p>2.N-Fix Cancellous Screw- SS316L,SS316LVM and titanium-Sterile and Non Sterile-3-E &D</p> <p>3. N-Fix Cannulated Cancellous Screw-SS316L,SS316LVM and titanium-Sterile and Non Sterile-3 & 5-E &D</p> <p>4.N-Link Interference Screw- Titanium-Sterile and Non Sterile-3 &5-E &D</p> <p>5.N-Knot Button--Sterile -3 &5-E &D</p>	CLAA/MD/GU J/12/2017-DC	<p>Mfg Licno No.G/28/1642 Validity till from 16-10-2017 to 15-10-2022</p> <p>(27-12-2017) (Approval Products)</p>
164.	M/s Perfect Medicare Products,Plot No:-26,27,28,Mahadev Estate,Part-4/A,B/H,Ranol Police Station,Amraiwadi,Ahmedabad-380026	<p>1.Bone Screw (Cortex Screw/LockingHead Screw/Cancellous Locking head screw/Interlocking Belt/PFN Screw/Angle locking head screw (Variable)/C.C.Screw/C.C.Screw(16 threads/32Threads)/DHS Leg Screw-SS316L,Titanium-5-Non Sterile-E &D</p> <p>2.Bone Plating Locking (Tlocking plate,T oblique locking plate/small locking plate/reconstruction combi locking plate/reconstruction locking plate/broad locking compression plate/narrow locking compression plate/T-Buttress locking plate/L-Buttress locking plate/distal femur</p>	CLAA/MD/GU J/14/2017-DC	<p>Mfg Licno No.G/28/1645 Validity till from 15-11-2017 to 14-11-2022</p> <p>(27-12-2017) (Approval Products)</p>

		locking plate/lateral tibial locking plate/distal medial locking plate/distal humeral dorsoleteral locking plate with and without lateral /Olecranon locking plate/Proximal tibia locking plate/particular proximal humerus locking plate/Lep Distal radius straight,Angles and oblique locking plate/Volar buttress locking plate/Volar distal radius locking plate/Anterilateral distal tibia locking plate/Extra-Articular distal humerus locking plate/Clavicle hook locking plate/superior anterior clavicle locking plate and extension plate/Variable angle lateral anerior clavicle locking plate/Variable Angle Dossolateral distal humeral locking plate with support and without support/Variable angle proximal ulna locking plate extra/variable angles dorsal distal radius L Oblique locking plate/Variable angle volar rim distal radius locking plate/Variable angle volar extra-articular distal radius locking plate/Variable angle lcp two-column volar distal radius locking plate/Variable angle condylar locking plate-SS316L and titanium-Non Sterile-E &D 3.Nail (Femoral Nail/Tibial Nail/PFN Nail/PFN Nail Short/TFN Nail/Expert Lateral Femoral Nail/PFNA Nail/Humres Nail/expert Humers Nail/Recon Nail/Supra Nail/Expert tibia Nail/Enders Nail/Square Nail/Square Radius Nail/Rush Nail/Eladstic Nail/Dynamic compression plate/Dyna Hip and Condylar./Cover leaf compression plate/Locking compression plate-SS316L,titanium-Non Sterile-E &D		
165.	M/s Sagar Surgical ,Plot No:-116,Amrut Udhoyognagar,Opp.Thakor dwar farm,vill:- Laskana,Tal:- Kamrej,Diist:- Surat,Gujarat	1.Interl ocking nail-SS316/Titanium/LVM-Non Sterile-E &D 2.Intter locking bolt- SS316/Titanium/LVM-Non Sterile-E &D 3.Locking Plates- SS316/Titanium/LVM-Non Sterile-E &D 4.LockinG Screw- SS316/Titanium/LVM-Non Sterile-E &D 5.Bone Plate- SS316/Titanium/LVM-Non Sterile-E &D 6.Bone Screw- SS316/Titanium/LVM-Non Sterile-E &D 7.Spine- SS316/Titanium/LVM-Non Sterile-E &D 8.Arthroscopy implants- SS316/Titanium/LVM-Non Sterile-E &D 9.Pins and Wires- SS316/Titanium/LVM-Non Sterile-E &D 10.Maxillofacial Implants- SS316/Titanium/LVM-Non Sterile-E &D	CLAA/MD/GU J/06/2017-DC	Mfg Licno No.G/28/1625 Validity till from 19-07-2017 to 18-07-2022 (27-12-2017) (Approval Products)

166.	<p>M/s Suhradam Ortho,Survey No:- 183,B/H, Satyam Mall,Nr.Cadila Pahrmaceuticals Ltd., 1st floor, Shop No:-1 to 11,Krishna corporation compound, sakhej-bavla highway,changodar- 3882213</p>	<ol style="list-style-type: none"> 1.Cortex Screws- SS316/Titanium/SS316LVM-Non Sterile-E &D 2.Cancellous Screws- SS316/Titanium/SS316LVM -Non Sterile-E &D 3.Cannulated Cancellous Screws- SS316/Titanium/SS316LVM -Non Sterile-E &D 4.ACL Screws- SS316/Titanium/SS316LVM -Non Sterile-E &D 5.Malleolar Screws- SS316/Titanium/SS316LVM -Non Sterile-E &D 6.Herbet Screws- SS316/Titanium/SS316LVM -Non Sterile-E &D 7.Acetabular screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D 8. Lokcing Screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D 9.PFA Cannulated blades- SS316/Titanium/ SS316LVM-Non Sterile-E &D 10.Schanz Screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D 11.Threaded tapper pins- SS316/Titanium/ SS316LVM-Non Sterile-E &D 12.Steinmann Pins-SS316/Titanium/ SS316LVM-Non Sterile-E &D 13.Denham Pins-SS316/Titanium/ SS316LVM-Non Sterile-E &D 14.Kirschnewr Wires-SS316/Titanium/ SS316LVM-Non Sterile-E &D 15.Washers-SS316/Titanium/ SS316LVM-Non Sterile-E &D 16.End Caps-SS316/Titanium/ SS316LVM-Non Sterile-E &D 17.Spacer and Wire mount-SS316/Titanium/ SS316LVM-Non Sterile-E &D 18.Intramedullary nail K Nails femur and tibia-SS316/Titanium/ SS316LVM-Non Sterile-E &D 19.V Nails femur and tibia-SS316/Titanium/ SS316LVM-Non Sterile-E &D 20.Enders Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 21.Square Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 22.Rush Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 23.Universal Tibia and femur nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 24.Dynamic Tibia Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 25.Upper bend tibia nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 26.Solid tibia nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D 	CLAA/MD/GU J/03/2017-DC	<p>Mfg Licno No.G/28/1606 Validity till from 01-04-2017 to 31-03-2022</p> <p>(27-12-2017) (Approval Products)</p>
------	---	--	----------------------------	---

27.Solid humeral nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
28.Femoral Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
29.Supra Condylar Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
30.Proximal Femoral Nails (PFN and APEN)- SS316/Titanium/ SS316LVM-Non Sterile-E &D
31.Trochantric femoral nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
32.Antegrade femoral nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
33.Sirua Femoral Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
34.Gamma Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
35.Elastic Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
36.TTC Nails-SS316/Titanium/ SS316LVM-Non Sterile-E &D
37.Bone Plates-SS316/Titanium/ SS316LVM-Non Sterile-E &D
38.One third plates-SS316/Titanium/ SS316LVM-Non Sterile-E &D
39.Semi Tubular plates-SS316/Titanium/ SS316LVM-Non Sterile-E &D
40.Reconstruction plates-SS316/Titanium/ SS316LVM-Non Sterile-E &D
41.T Plates-SS316/Titanium/ SS316LVM-Non Sterile-E &D
42.L Buttress Plate- SS316/Titanium/ SS316LVM-Non Sterile-E &D
43.Condylar Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
44.Angles Balde Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
45.DHS Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
46.DCS Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
47.Hook Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
48.Pediatric Hip Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
49.DSistal Radius Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
50.Distal fibula plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
51.Proximal Tibia plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D
52.Distal tibia plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D

		<p>53.Humerus Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>54.Clavicle Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>55.Pilon Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>56.Spoon Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>57.Olecranon Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>58.Cloverleaf plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>59.proximal femur plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>60.Distal femur plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>61.Stabilization plate- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>62. Metaphseal plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>63.Calcaneal plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>64.Cobra Plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>65.Spinal implants (PHI Systems) Monoaxial Screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>66.Polyaxial Screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>67.Reduction polyaxial screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>68.Reduction mono screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>69.Inner screws- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>70.Hooks- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>71.Connectors- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>72.Cervical screws and plates- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>73.Polyaxials screw occipital plate- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p> <p>74.Ankle Implants- SS316/Titanium/ SS316LVM-Non Sterile-E &D</p>		
167.	M/s Ortho Care,A-31,Sector-A-4,Tronica City Loni, Ghaziabad,Uttar Pradesh	<p>1.LCP Clavicular plates,3.5mm Right-Titanium/SS316L/SS316LV</p> <p>2.LCP Clavicular plates,3.5mm Right-Titanium/SS316L/SS316LV</p> <p>3.Locking plate superior anterior clavicle with lateral extention 2.7/3.5mm-right- Titanium/SS316L/SS316LV</p> <p>4. Locking plate superior anterior clavicle with lateral extention 2.7/3.5mm-Left- Titanium/SS316L/SS316LV</p> <p>5.LCP Superior anterior clavicle plate 3.5mm-right-</p>	CLAA/MD/UP/03/2010-DC	Mfg Licno No.05/SC/P of 2011 Validity till from 07-09-2016 to 06-09-2021

	<p>Titanium/SS316L/SS316LV</p> <p>6. LCP Superior anterior clavicle plate 3.5mm-left- Titanium/SS316L/SS316LV</p> <p>7.LCP Clavicle S-type plate-right- Titanium/SS316L/SS316LV</p> <p>8. LCP Clavicle S-type plate-left- Titanium/SS316L/SS316LV</p> <p>9.3.5mm metaphyscial lacking compression plate- Titanium/SS316L/SS316LV</p> <p>10.Proximal humeral locking plate,short- Titanium/SS316L/SS316LV</p> <p>11. Proximal humeral locking plate,long -Titanium/SS316L/SS316LV</p> <p>12.3.5mm locking compression plate- Titanium/SS316L/SS316LV</p> <p>13.LCP Hook Plate 3.5mm 3 Holes length 62mm- Titanium/SS316L/SS316LV</p> <p>14.LCP Humeral distal/dorsal plate 3.5mm-right- Titanium/SS316L/SS316LV</p> <p>15. LCP Humeral distal/dorsal plate 3.5mm-left- Titanium/SS316L/SS316LV</p> <p>16.LCP Extra articular distal humerus right- Titanium/SS316L/SS316LV</p> <p>17.LCP Extra Articular distal humerus left- Titanium/SS316L/SS316LV</p> <p>18.LCP Plecranon plate 3.5mm type-1 Right- Titanium/SS316L/SS316LV</p> <p>19.LCP Olecranon plate plate 3.5mm type-1 Right- Titanium/SS316L/SS316LV</p> <p>20.Locking t Plate right angle for 3.5mm screws- Titanium/SS316L/SS316LV</p> <p>21.Locking T Plate right angle for 3.5mm screws- Titanium/SS316L/SS316LV</p> <p>22.LCP T Plate oblique angle 3.5mm right- Titanium/SS316L/SS316LV</p> <p>23. LCP T Plate oblique angle 3.5mm left- Titanium/SS316L/SS316LV</p> <p>24.LCP Volar column distal radius 2.7mm,8 holes right- Titanium/SS316L/SS316LV</p> <p>25. LCP Volar column distal radius 2.7mm,8 holes left- Titanium/SS316L/SS316LV</p>		<p>(27-12-2017)</p> <p>(Approval for Renewal Products)</p>
--	--	--	--

26. LCP Volar column distal radius 2.7mm right-
Titanium/SS316L/SS316LV

27. LCP Volar column distal radius 2.7mm left-
Titanium/SS316L/SS316LV

28.LCP Extra Articular volar distal radius plate 2.7mm right type-1-
Titanium/SS316L/SS316LV

29. LCP Extra Articular volar distal radius plate 2.7mm left type-1-
Titanium/SS316L/SS316LV

30. LCP Extra Articular volar distal radius plate 2.7mm left type-1I-
Titanium/SS316L/SS316LV

31. LCP Extra Articular volar distal radius plate 2.7mm right type-1I-
Titanium/SS316L/SS316LV

32.LCP wrist fusion plate standard/short/straight bend-
Titanium/SS316L/SS316LV

33.3.5mm reconstruction locking plate- Titanium/SS316L/SS316LV

34.4.5mm reconstruction locking plate- Titanium/SS316L/SS316LV

35.LCP Pediatric hip plate 3.5, 100 degree, 110 degree, 120 degree, 150
degree- Titanium/SS316L/SS316LV

36.Locking Dcs plate 95 degree- Titanium/SS316L/SS316LV

37.Locking condylar plate 4.5/5.0mm right- Titanium/SS316L/SS316LV

38.Narrow locking plate for 4.5mm/5.0mm screws-
Titanium/SS316L/SS316LV

39.Broad locking plate for 4.5mm/5.0mm screws-
Titanium/SS316L/SS316LV

40.Locking distal femur plate right/left-4.5mm/5.0mm screws-
Titanium/SS316L/SS316LV

41.Lateral distal femur LCP for osteotomy -Titanium/SS316L/SS316LV

42.Lateral high tibia LCP for osteotomy 34-
Titanium/SS316L/SS316LV

43.LCP targeting proximal tibia plate 4.5/5.0 left-
Titanium/SS316L/SS316LV

44.3.5/4.5mm metaphyseal locking compression-
Titanium/SS316L/SS316L

45.LCP Anterolateral distal tibi plate right/left-
Titanium/SS316L/SS316LV

46.Fibular distal lateral locking plate- Titanium/SS316L/SS316LV
47.Locking anatomical distal tibia plate for fibulae-right/left- Titanium/SS316L/SS316LV
48.LCP Calcaneal plates 3.5mm left/right- Titanium/SS316L/SS316LV
49.Lag screws 12/60 to 120mm- Titanium/SS316L/SS316LV
50.Compression screws 4X36- Titanium/SS316L/SS316LV
51.Mini cortex screw dia 2.4mm/self tapping- Titanium/SS316L/SS316LV
52. Cortex screw dia 4.5mm self tapping- Titanium/SS316L/SS316LV
53.5.0mm locking head screw self tapping-Titanium/SS316L/SS316LV
54.5.0mm locking head cancellous screws thread-16mm/32mm self tapping- Titanium/SS316L/SS316LV
55.5.0/3.5mmmm locking head cancellous screws full /short threaded self tapping- Titanium/SS316L/SS316LV
56.3.5mm/5.0mm locking head screws,self tapping,self drilling- Titanium/SS316L/SS316LV
57.Sapcer,3.5mm/5.0mm 2mm length- Titanium/SS316L/SS316LV
58. Mini cortex screw dia-2.0/2.7mm- Titanium/SS316L/SS316LV
59. Cortex screw dia-3.5/4.5mm- Titanium/SS316L/SS316LV
60.Cortex screw dia-3.5/4.5mm self tapping- Titanium/SS316L/SS316LV
61. Cancellous screw dia-2.7/3.5mm/4.0/4.5/6.5 mm full/short threaded- Titanium/SS316L/SS316LV
62.Malleolar screw dia 4.5mm- Titanium/SS316L/SS316LV
63.Washer for cancellous screws- Titanium/SS316L/SS316LV
64.Ashnis screw- Titanium/SS316L/SS316LV
65.Large cannulated cancellous screw dia-7.0mm 16/32 mm thread- Titanium/SS316L/SS316LV
66. Small cannulated cancellous screw dia-7.0mm 16/32 mm short thread- Titanium/SS316L/SS316LV
67.Shaft screws for LC-DCP dia 3.5/4.5mm- Titanium/SS316L/SS316LV
68.Steffee screws dia-4.5/5.55mm- Titanium/SS316L/SS316LV
69.Hip screw with compression screw- Titanium/SS316L/SS316LV
70.DHS/DCS Compression screw- Titanium/SS316L/SS316LV

71.Locking bolt dia 3.9/4.9mm- Titanium/SS316L/SS316LV
72.Locking bolt dia 6.0mm length-20mm to 120 mm for dfn- Titanium/SS316L/SS316LV
72.Plug screw for interlocking I.M.L.Nails- Titanium/SS316L/SS316LV
73.End Cap for undreamed humeral nail- Titanium/SS316L/SS316LV
74.Cacellous traction screws- Titanium/SS316L/SS316LV
75.Mini straight plate for 2.0mms crews- Titanium/SS316L/SS316LV
76.Mini L-Plates fpr 2.00mm screws-left and right- Titanium/SS316L/SS316LV
77.Dynamic compression plates for 2.0/2.7mms crews- Titanium/SS316L/SS316LV
78.Quater tubular plates for 2.7mms crews- Titanium/SS316L/SS316LV
79.1/3 tubular plate-3.5mm /screws with collar- Titanium/SS316L/SS316LV
80.Clover leaf plate- Titanium/SS316L/SS316LV
81.Y-Reconstruction plate for humeral condyle- Titanium/SS316L/SS316LV
81.Calcaneal plate-3.5mm- Titanium/SS316L/SS316LV
82.Pelvic reconstruction plate for 4.5mms crews- Titanium/SS316L/SS316LV
83.Narrow dynamic compression plate for 4.5mms crews- Titanium/SS316L/SS316LV
84.Angle blade-plate 130degree with dynamic holes- Titanium/SS316L/SS316LV
85.Spoon Plate- Titanium/SS316L/SS316LV
86.Hook Plate -4.5mm- Titanium/SS316L/SS316LV
87.Proximal tibia plate-right/left- Titanium/SS316L/SS316LV
88.Lengtheing p[late 4.5mm-narrow holes/Broad- Titanium/SS316L/SS316LV
89.H-Plate for cervical spine- Titanium/SS316L/SS316LV
90.Steffee plate- Titanium/SS316L/SS316LV
91.D.H.S Barrel plate 120,125,130,135,140,145,150 degree double bend with dynamic holes- Titanium/SS316L/SS316LV
92. D.H.S Barrel plate,130,135,140,145,150 degree with dynamic hole short barrel- Titanium/SS316L/SS316LV

		<p>93.D.C.SBarrel plate 95 degree with dynamic holes- Titanium/SS316L/SS316LV</p> <p>94.Mclaugin plate-- Titanium/SS316L/SS316LV</p> <p>95.Jewett nail plate 135 degree-- Titanium/SS316L/SS316LV</p> <p>96. Wain right osteotomy plate-- Titanium/SS316L/SS316LV</p> <p>97.Thompson Femoral Head standard/Narrow stem-Non/gamma sterile-- Titanium/SS316L/SS316LV</p> <p>98.Austin moore femoral head standard/narrow/broad stem-Non/gamma sterile-- Titanium/SS316L/SS316LV</p> <p>99.Bipolar prosthesis-talwar type-hard top-- Titanium/SS316L/SS316LV</p> <p>100.Bipolar prosthesis-talwakar type-hard top fenestrated gamme sterile-- Titanium/SS316L/SS316LV</p> <p>101.Universal tibial nail solid non cannulated- Titanium/SS316L/SS316LV</p> <p>102.Universal femoral nail-- Titanium/SS316L/SS316LV</p> <p>103.Unreamed humeral nail-- Titanium/SS316L/SS316LV</p> <p>104.Intramedullary supra condylar nail-5 holes-- Titanium/SS316L/SS316LV</p> <p>105.INtramodullary supra condylar nail-multiple holes-- Titanium/SS316L/SS316LV</p> <p>106.Spare bolt and washer for S.P.Nail-- Titanium/SS316L/SS316LV</p> <p>107.Smaith Patterson nail with bolt and washer-- Titanium/SS316L/SS316LV</p> <p>108.Bend K-Nail Dia-0.6mm to 13mm-length 20cm to 46cm-with 1 cm difference--SS316L</p> <p>109.V Nail dia-0.3m mto 11mm, length-15 cm to 40cm-with 1cm difference- SS316L</p> <p>110.Ulna DIa-2.0mm to 4.00mm length-16cm to 30cm-with 1 cm difference- SS316L</p> <p>111.Radius DIa-2.0mm to 4.00mm length-16cm to 30cm-with 1 cm difference- SS316L</p> <p>112.Rush nail DIa-2.0mm to 4.00mm length-16cm to 30cm-with 1 cm difference- SS316L</p> <p>113.Stainless Steel wire 16swg to 30 swg- DIa-2.0mm to 4.00mm</p>		
--	--	--	--	--

length-16cm to 30cm-with 1 cm difference- SS316L
114.Kriscner wire both side trocar - SS316L
115.Kriscner wire on side trocare-round end-SS316L
116.Threaded kriscner wire dia 1.5mm to 2.5mm length-150mm to 225mm –SS316L
117.Steinman pin dia 3.0mm to 5.0mm length 150mm to 230mm-SS316L
118.Centrally threaded Steinman pin dia 4.0mm to 6.0 length 225mm-SS316L
119.Knowles pin dia 3.0/4.0mm-SS316L
120.Cerclege witre with loop 175swg to 21swg-SS316L
Hart Shell-SS316L
121.Schanz Screw-Trocare Points-SS316L
122.Clip washer-SS316L
123.Harrington rods-SS316L
124.Sacral rod with six nuts-SS316L
125.Illizarov wire-SS316L
126.Illizarov wire bayonet point-SS316L
127.Olive Wire-SS316L
128.Canceelous Schanz Screw-Trocare points-SS316L
129.Steinman pin Dia 3.0mm to 5.0mm –SS316L
130.Opening Wedge osteotomy plate-SS316L
131.Cross link bar 40mm to 80mm-Titanium
132.ROD-Titanium
133.Pedicle screw inner nut-Titanium-Cross link hooks-Titanium
134.Mono reduction screws-Titanium
135.POLyaxial screws-Titanium
136.Anterior certical locking plate-Titanium
137.Mesh-Titantium
138.Cortex Screws 0 1.5mm length-6mm to 14mm –Tiatnium
139.Staight plate 0.85mm thick W/1.5mm screws-Titanium
140.L-Plate 4h,0.85mm thick W/1.5MM screws-Titanium
141.Angles L plate 4h 0.85mm thick W/1.5mm screws-titanium
142.Crtex Screws 0 0.2mm-Titanium
143.Short T plate 4h 1mm thicl W/2.0mm screws-Titanium

		<p>144.Hex cortex ,partially threased 2.7mm –Titanium Straight Plate,2mm thick W/2.7MM screws-Titanium</p> <p>145.Angled L-Plate, 5h,1.2mm thick W/2.7mm screw left/right-Titanium</p> <p>146.Staight L Plate 5h,1.2mm thick W/2.7mm screw left/right-Titanium</p> <p>147.Mini T Plate,5h,1.2mm thick W/2.7mm-Titanium</p> <p>148.W-plate,7h,1.2mm thick W/2.7mm screws-Titanium</p>		
168.	M/s Angioplast Pvt.Ltd.,Plot No:-4083and 4804 Phase-IV, GIDC Vatva, Ahmedabad-382445.Gujarat,India	<p>1.Haemodialysis blood tubing set-Discharge hook-PP,Rotating luer lock-ABS, C Clamp-HDPE,Injection port-ABS, Stopper-Latex Free,Pump connector-PVC, Pump tube-PVC, Female luer lock-ABS,Cap-PP, Chamber Cap-PVC, Chamber-PVC, Dialysis Connector-ABS,Cap-PP,Tubing-PVC-Sterile-5 from the date of manufacturing-E &D</p> <p>2.IV Cannula-Catheter-PTFE,Protective Cap-PP,Injection Port with Cap-PE,Non return valve-EVA,Cannula-SS- Sterile-5 from the date of manufacturing-E &D</p> <p>3.Scalp vein set-Wings/Tubing-PVC,Cannula-SS,Guard-LDPE- Sterile-5 from the date of manufacturing-E &D</p> <p>4.3 Way stop cock-Body-Polycarbonate,Handle-HDPE,ProtectiveCap-PP/HDPE- Sterile-5 from the date of manufacturing-E &D</p> <p>5. Peritoneal dialysis set-Spike guard-pp,Spike,Drip chamber, WYE Connector, Tubing,funnel open-PVC,Male fitting adaptor-ABS,Protective Cap-HDPE- Sterile-5 from the date of manufacturing-E &D</p> <p>6.Infant feeding tube-Catheter-PVC/DEHP FREE,Female connector-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>7.Ryle`s tube-Catheter,funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>8.Deodenal ryle`s tube- Catheter,funnel,Feedingc onnector-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>9.Stomach Tube-Catheter,funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>10.Levins Tube- Catheter,funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>11.Corrugated drainage sheet-Corrugated drainage sheet-PVC-- Sterile-5</p>	CLAA/MD/GU J/17/2010-DC	<p>Mfg Licno No.G/28/1357 Validity till from 25-11-2010 to 24-11-2020</p> <p>(27-12-2017) (Approval for Renewal Products)</p>

		<p>from the date of manufacturing-E &D</p> <p>12.Male External Catheter size various-Male catheter-natural rubber latex-- Sterile-5 from the date of manufacturing-E &D</p> <p>13.Infant Mucus extractor-Mucus cylinder-PPCP,Cylinder cap,suction catheter,suction tube,funnel-PVC-- Sterile-5 from the date of manufacturing-E &D</p> <p>14. Infant Mucus extractor with filter-Mucus cylinder-PPCP,Cylinder cap,suction catheter,suction tube,funnel-PVC,filter-PP-- Sterile-5 from the date of manufacturing-E &D</p> <p>15.Nasal oxygen cannula-Nastril, Star tubing,two way connector,Machine connector-PVC-Sterile-5 from the date of manufacturing-E &D</p> <p>16.Yankauer suction set-Yankauer handle-rigid PVC, Connector,tubing-PVC--Sterile-5 from the date of manufacturing-E &D</p> <p>17.TUR Set-Spike guard-PP,Spike and extra spike-ABS, Tubing, WYE Connector-PVC,Tubing –PVC with cotton onside, injection site-Latex-- Sterile-5 from the date of manufacturing-E &D</p> <p>18. Nelaton Catheter-Catheter &Funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>19. Suction Catheter- Catheter &Funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>20.Umbical Catheter- Catheter &Funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>21.Thoracic catheter with trocar-thoracic catheter-PVC, Trocar-Aluminium,Connector for thoracic catheter and cap for trocar-ABS-Sterile-5 from the date of manufacturing-E &D</p> <p>22.Rectal Tube- Catheter &Funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>23.Karman cannula with and without adaptor size:-FG4,5,6,7,8,10-Adaptor and cannula-LDPE- Sterile-5 from the date of manufacturing-E &D</p> <p>24.Closed wound cution unit- Bellow container –LDPE,Perforated catheter-PVC,Curved Needle-SS,Connecting Tube-PVC, C Clamp-HDPE, Honger-HDPE,Bellow container cap-HDPE,Inner Cap-PVC-Sterile-5 from the date of manufacturing-E &D</p>		
--	--	--	--	--

		<p>25. Closed wound suction unit- Bellow container –LDPE, Perforated catheter-PVC, Curved Needle-SS, Connecting Tube-PVC, C Clamp-HDPE, Honger-HDPE, Bellow container cap-HDPE, Inner Cap-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>26. Mini Vac Set- Bellow container –LDPE, Perforated catheter-PVC, Curved Needle-SS, Connecting Tube-PVC, C Clamp-HDPE, Honger-HDPE, Bellow container cap-HDPE, Inner Cap-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>27. Foley Balloon Catheter-Foley Catheter-Latex with elastomer coated-- Sterile-5 from the date of manufacturing-E &D</p> <p>28. Silicon Floey Balloon Catheter-Foley Catheter-100% Silicon-- Sterile-5 from the date of manufacturing-E &D</p> <p>29. Endotracheal tube size:-FG2,5,3,3,5,4,4,5,5,5,5,6,6,5,7,7,5,8,8,5,9- Tube –PVC Connector-HDPE-- Sterile-5 from the date of manufacturing-E &D</p> <p>30. Endotracheal tube size:-FG2,5,3,3,5,4,4,5,5,5,5,6,6,5,7,7,5,8,8,5,9- Reinforced tube-PVC Connector-HDPE-- Sterile-5 from the date of manufacturing-E &D</p> <p>31. Tracheostomy tube size:-FG2,5,3,3,5,4,4,5,5,5,5,6,6,5,7,7,5,8,8,5,9- Tube-PVC, Connector-HDPE, Neck Plate-Cloth Tie, Obturator-LDPE--- Sterile-5 from the date of manufacturing-E &D</p> <p>32. Thoracic catheter-Catheter-PVC, Zipped connector-PVC, Thoracic connector-ABS- Sterile-5 from the date of manufacturing-E &D</p> <p>33. Abdominal drain kit-Bag-Soft PVC Cloth, Catheter-PVC, Tubing-PVC, Connector-PVC, Hanger-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>34. Urethral catheter-Catheter and funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>35. Nasal oxygen catheter- Catheter and funnel-PVC- Sterile-5 from the date of manufacturing-E &D</p> <p>36. A.V. Fistula needle-Catheter-PTFE, Protective cap-PP, Injection port with cap-PE, Non return valve-EVA, Cannula-SS- Sterile-5 from the date of manufacturing-E &D</p> <p>37. Catheter Mount (Double Swivel)-Female Connector-K esin, Double Swivel Mount-LDPE, Corrugated Tube connector-K Resin, Corrugated</p>		
--	--	---	--	--

		<p>tube-EVA, cap with closure Cap-PVC, Female connector-K Resin-Sterile-5 from the date of manufacturing-E &D</p> <p>38. Peritoneal dialysis catheter Kit-Trocar-SS,Indicator ball-HDPE, Trocar Handle-HDPE, Catheter-PVC, Duodenal stopper-PVC, Bulb latex-Natural rubber latex,Tubing-PVC, Junction Unit-LDPE, Flow Regulator-ABS- Sterile-5 from the date of manufacturing-E &D</p> <p>39.Pigtail Catheter urethral tube/set-Drainage catheter-PVC/PU,Pusher-LDPE/PP,Dilator Set-Flexible PU,Guide Wire-Metal coated,Umbilical cord clamp-PP-Sterile-5 from the date of manufacturing-E &D</p> <p>40.Silicon Wound suction unit-Bellow container-Silicon,Trocar-SS,Catheter-PVC-Sterile-5 from the date of manufacturing-E &D</p> <p>41.Nasopharyngeal Airway- Nasopharyngeal aiway-PVC-5 from the date of manufacturing-E &D</p> <p>42.Laryngeal mask airway-Tube-Silicom/PVC,Machine End connector-Tr.ABS, Pilot Balloon with inflation line-silicon/PVC, Non return valve,Ring-Latex- Sterile-5 from the date of manufacturing-E &D</p> <p>43.Safety Cannula-Catheter-PTFE,Protection Cap-PP,Injection port with cap-PE, Non return valve-EVA,Cannula-SS- Sterile-5 from the date of manufacturing-E &D</p> <p>44.Haemodialysis blood tubing set-Discharge hook-PP, Rotating laser lock-ABS, C Clamp-HDPE, Injection port-ABS, Stopper-LLatex Free, Pump connector-PVC, Pump tube-PVC, Female luer lock-ABS, Cap-PP, Chambe Cap-PVC, Chamber-PVC, Dialyser connector-ABS, Cap-PP, Tubing-PVC Sterile-5 from the date of manufacturing-E &D</p> <p>45.IV Cannula-Catheter –PTFE,Protection cap-PP, Injection port with ca-PE, Non return valve-EVA,Cannula-SS- Sterile-5 from the date of manufacturing-E &D</p> <p>46.Scalp vein set-wings-PVC,Tubing-PVC,Cannula-SS,Guard-LDPE-Sterile-5 from the date of manufacturing-E &D</p> <p>47.Karman cannula with and without adapter size:-FG4,5,6,7,8,10-Adaptor and Canula-LDPE- Sterile-5 from the date of manufacturing-E &D</p> <p>48.A.V.Fistula needle-Catheter-PTFE,Protective Cap-PP, Injection port with cap-PE, Non rteun Vlave-EVA,Cannula-SS- Sterile-5 from the date of manufacturing-E &D</p>		
--	--	---	--	--

		49.Percutaenous pigtail nephrostomy catheter set-Drainage catheter-LDPE,Dilator Set-Flexible PU,Needlw ith style-SS,Guide wire-metal coated,Scalpel-SS,1 way stop cock (body)-Polycarbonate,1 way stop cock(Handle)-Handle ,Cap-LDPE- Sterile-5 from the date of manufacturing-E &D 50.Silicon Wound suction unit-Bellow container –Silicon,Trocare-SS,Catheter-PVC- Sterile-5 from the date of manufacturing-E &D 51.		
169.	M/s HLL Life Care Ltd.,Sreekariyam akkulam Road,Thiruvananthapuram	1.Surgical mesh/hernia mesh-Polypropylenen yarn-Sterile-5- E&D	CLAA/MD/KE R/01/2012-DC	Mfg Licno No.27/28/92 Validity till 31-12-2021 (28-12-2017) (Approval for Renewal Products)
170.	M/s Shajanand Laser technology Ltd.,A-8,GIDC Electronic Esttae, Sector-25,Gandhi Nagar-382016	1.Cobalt chromium coronary stent(with delivery system)-Cobalt Chromium-Sterile and non sterile-24 months-E &D 2.Cobalt chromium rapamycin eluting- Cobalt chromium -Sterile and non sterile-18months-E &D 3. Cobalt chromium rapamycin eluting coronary system (with delivery system)- Cobalt chromium -Sterile and non sterile-18months-E &D 4.PTCA Balloon catheter-Nylon- Sterile and non sterile-18months-E &D	CLAA/MD/GU J/08/2009-DC	Mfg Licno No.G/28/1314 Validity from 27/08/2014 to 26/08/2019 (29-12-2017) (Approval for Additonal Products)
171.	M/s Meril Life Sciences Pvt.Ltd., Survey No:- 135/2a and 135/2/B,174/2 (M1-M2,Meril Park,4 th floor),Near G.M.Bilakhia stadium,Muktanand Marg,Chala,Vapi-	1.Rapid exchange percutaneous transluminal coronary angioplasty balloon dilatation catheter-Nylon, Pebax, PTFE/Uncoated stainless steel-Sterile-3- E &D 2.Rapid exhcnge Percutaneous translumlnial coronary angioplasty balloon dilatation catheter-Nylon,Pebax,PTFE ro uncoated stainless steel-4- E&D	CLAA/MD/GU J/01/2008-DC (Pt.V)	Mfg Licno No.G/28/1273 Validity from 02/06/2013 to 19/06/2018

	396191, Valsad, Gujarat			(28-12-2017) (Approval for Additional Products)
172.	M/s Meril Life Sciences Pvt.Ltd., Survey No:- 135/2a and 135/2/B,174/2 (M1-M2,Meril Park,4 th floor),Near G.M.Bilakhia stadium,Muktanand Marg,Chala,Vapi-396191, Valsad,Gujarat	1.Biodegradable based sirolimus eluting coronary stent system (drug eluting coronary stent)-L605 Cobalt chromium alloy-Sterile-2- E &D	CLAA/MD/GU J/01/2008-DC (Pt.V)	Mfg Licno No.G/28/1273 Validity from 02/06/2013 to 19/06/2018 (28-12-2017) (Approval for Additional Products)
173.	M/s Care Group Sight Solution LLP,Block No:- 310/C &E, Village:- Dhabasa,Tal:-Padra,Dist :- Vadodara-391440	1.Foldabale Hydrophilic acrylic intraocular lens (sterile) GENERIC-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D 2. Foldabale Hydrophilic acrylic intraocular lens (sterile) Idiff pLUS-Hydrophilic acrylic co polymer with UV blocker and or yellow dye-Sterile- E &D 3. Foldabale Hydrophilic acrylic intraocular lens (sterile) IPCL-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D 4. Foldabale Hydrophilic acrylic intraocular lens (sterile) OCUFLEX-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D 5. Foldabale Hydrophilic acrylic intraocular lens (sterile) TRiDIFF-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye-Sterile- E &D 6. Foldabale Hydrophilic acrylic intraocular lens (sterile) TRiDIFF TORIC- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye -Sterile- E &D 7. Foldabale Hydrophilic acrylic intraocular lens (sterile) ULTIMA	CLAA/MD/GU J/15/2017-DC	Mfg Licno No.G/28/1321 Validity from 12/05/2017 to 11/05/2022 (28-12-2017) (Approval for Products)

		<p>PLUS SMART TORIC- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye -Sterile- Export</p> <p>8. Foldabale Hydrophilic acrylic intraocular lens (sterile) ACRIOL- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye- Sterile- E &D</p> <p>9. Foldabale Hydrophilic acrylic intraocular lens (sterile) ACRIOL TD- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D</p> <p>10. Foldabale Hydrophilic acrylic intraocular lens (sterile) AcriDiff - Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye- Sterile- E &D</p> <p>11. Foldabale Hydrophilic acrylic intraocular lens (sterile) ACRIOL EC- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye)-Sterile- E &D</p> <p>12. Foldabale Hydrophilic acrylic intraocular lens (sterile) ACRIOL EC TORIC- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye)-Sterile- E &D</p> <p>13. Foldabale Hydrophilic acrylic intraocular lens (sterile) ATM- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D</p> <p>14. Single piece PMMA Intraocular lens(sterile)Eye-O-Care-Poly Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- E &D</p> <p>15. Foldabale Hydrophilic acrylic intraocular lens (sterile) EYEDIFF- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- E &D</p> <p>16. Foldabale Hydrophilic acrylic intraocular lens (sterile) EYEDIFF plus- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>17.Foldabale Hydrophilic acrylic intraocular lens (sterile) EYEOL- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>18. Foldabale Hydrophilic acrylic intraocular lens (sterile) EYEPCL- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye - Sterile- Export</p>		
--	--	---	--	--

		<p>19. Foldabale Hydrophilic acrylic intraocular lens (sterile) ULTIMA-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>20. Foldabale Hydrophilic acrylic intraocular lens (sterile)GOLD-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>21. Foldabale Hydrophilic acrylic intraocular lens (sterile) HYDIFF-T-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>22. Foldabale Hydrophilic acrylic intraocular lens (sterile) HYFLEX-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>23. Foldabale Hydrophilic acrylic intraocular lens (sterile) HYDIFF TORIC- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>24. Foldabale Hydrophilic acrylic intraocular lens (sterile) HYFLEX Y-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>25.Single Piece PMMA Intraocular lens (Sterile) EYEOL– Polymethacrylate with UV locker and or yellow dye-Sterile-Export</p> <p>26. .Single Piece PMMA Intraocular lens (Sterile) PMMA GOLD – Polymethacrylate with UV locker and or yellow dye-Sterile-Export</p> <p>27. Foldable Intraocular lens hydrophilic UV (MEDICAL)-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>28. Foldable Intraocular lens hydrophilic UV (PRESBYSMART PLUS)-Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>27.PMMA Intraocular lens-Poly(methyl methacrylate) with UV balocker and or yellow dye-Sterile-Export</p> <p>28.Foldable Hydrophilic acrylic intracolar lens(Sterile)- Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>29. Foldable Hydrophilic acrylic intracolar lens(SMARTIOL CRYSTAL) Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>30. Foldable Hydrophilic acrylic intracolar lens(Bluevision) Hydrophilic</p>		
--	--	--	--	--

		<p>acrylic Co-Polymer with UV blocker and or yellow dye -Sterile- Export</p> <p>31. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA ASPPERHIC CLEAR)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>32. Foldable Hydrophilic acrylic intracolar lens(Sterile) (vista SPHERE)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>33. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA YOUTH HYDROPHILIC/HYDROPHILIC AD)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>34. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA ASPPERHIC GOLD)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>35. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA ACRYPHOBIC CLEAR)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>36. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA TORIC)Hydrophilic acrylic Co-Polymer with UV blocker and or yellow dye Sterile- Export</p> <p>37. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA AC IOL)Poly (methyl methacrylate with UV blocker and or yellow dye Sterile- E &D</p> <p>38. Foldable Hydrophilic acrylic intracolar lens(Sterile) (VISTA PMMA)-Poly (methyl methacrylate with UV blocker and or yellow dye Sterile- E & D</p>		
174.	M/s GPC Medical Ltd., C-3, Sector B-1, Trans Delhi Sinature City, Loni, Ghaziabad-201102, Uttar Pradesh	<ol style="list-style-type: none"> 1. Bone Screws 2. Knee Surgery Implants 3. D.H.S/Lag Screw 4. D.H.S Compression screw/blocket screw 5. Paediatric D.H.S.Screw 6. Locking Head Screws 7. Bone Plates 8. Maxillofacial Plates 9. Small Fragment Plates 10. Large Fragment Plates 	CLAA/MD/03/2015-DC	<p>Mfg Licno No./SC/P of 2018 Validity</p> <p>(28-12-2017)</p> <p>(Approval for Products)</p>

		11.Angle blade plates for Intertrochanteric femoral ostrotomies in Adults 12.Distal Femur Bone Plate 13.Locking plate system 14.LockingLateral Proximal tibia plate left/Right 15.Locking Rigid Plate 16.Intramedullary Interlocking Nails 17.Unreamed Tibial Nails (Stainless Steel) 18.Multiangle Tibia Nail 19.Tibial interlocking nails distal 20.Femoral Nails 21.Unreamed Femoral Nail 22.Supra Condylar Nail/IMSC Nail 23.Humerus Nail Cannulated 24.Reconstruction Nail (Right) 25.Titanium Elastic Nail 26.V-Nails 27.Square Nail for Ulna 28.Pins,Wires,Staples 29.Arthroplasty-Femoral Heads 30.Single Lock Mono Axial Screw Pedicle 31.Harrington Rod 32.Harrington Hook 33.Expandable Cage 34.Plif Cage 35.Tlif Cage		
175.	M/s Bio Med Helath Care Products Pvt.Ltd.,Plot No:-62-65,Setcor-6A,SIDCUL,Ranipur,Haridwar-249403,Uttarakhand	1.Bio Flon-Polyproylyene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D 2.Neo Kath-Polypropylyene,Polyethylene, SS Needle,Catheter,Poly Acetal- Sterile-5- E &D 3.Maxima- Polyproylyene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D 4.Agary- Polyproylyene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D 5.Vasocath- Polyproylyene,Polyethylene,SS Needle, Catheter,Poly	CLAA/MD/UA /02/2007-DC	Mfg Licno No.123/UA/SC/P-2007 Validity from 02/11/2017 to 01/11/2022 (29-12-2017) (Approval for Additonal Products)

	<p>acetal,Silicon Tube-Sterile-5- E &D</p> <p>6.Vasocath Advantive- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>7.Vagon- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>8.Plus Med- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>9.Romed- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>10.Union- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>11.Irismed- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>12.Medicare- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>13.Cadiflo- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>14.Medicann- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>15.Heur- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>16.Bio Pen- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal-Sterile-5- E &D</p> <p>17.Bio Neo- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal-Sterile-5- E &D</p> <p>18.Medi Life- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>19.Solidor- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal-Sterile-5- E &D</p> <p>20.Bio Sis- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>21.Demotek- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>22.Mediplast- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p>		
--	---	--	--

		<p>23.Holden- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>24.SafeMed- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>25.ANGELTOUCH- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>26.Rosemed- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>27.Beflon- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>28.Kareflon- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>29.Bulife- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>30.Unimed- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>31.IEmedical- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>32.HM Healthcare- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>33.Asco- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>34.Lifecare- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>35.Heflon- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>36.AG Cath- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>37.Anarmed- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>38.Flawless- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>39.Ritemed- Polyproylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>40.Manaflo- Polyproylene,Polyethylene,SS Needle, Catheter,Poly</p>		
--	--	---	--	--

		<p>acetal,Silicon Tube-Sterile-5- E &D</p> <p>41.Kings- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>42.Health- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p> <p>43.Top Care- Polypropylene,Polyethylene,SS Needle, Catheter,Poly acetal,Silicon Tube-Sterile-5- E &D</p>		
176.	M/s V.R.Ortho,Plot No:- 3,Anand Nagar,Behind Prakash Talkies,Daimond City,Navsari- 369445,Gujarat	<p>1.Cortex Screw-SS316L and titanium-Non Sterile-5- E&D</p> <p>2.Cancellous Screw- SS316L and titanium-Non Sterile-5- E&D</p> <p>3.Cannulated cancellous screw- SS316L and titanium-Non Sterile-5- E&D</p>	CLAA/MD/GU J/07/2017-DC	<p>Mfg Licno No.G/28/1630</p> <p>Validity from 23/08/2017 to 22/08/2022</p> <p>(29-12-2017)</p> <p>(Approval for Products)</p>
177.	M/s Indius Medical Pvt.Ltd., Unit No:- 12,electrnoci Co-Operation Estate,Pune-Satara Road,Pune-411009	<p>1.Cervical spine interbody spacer-Titanium alloy(tiaLV) and polyether ether ketone (PEEK)-Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>2.Anterior Cervical spine plate system-Titanium- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>3.Anterior cervical integrated plate and spacer system- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>4.Lumbar Spine Interbody spacer- Polyether ether ketone (PEEK) and Titanium alloy (tiaLV)- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>5.Spine Interbody spacer corpectomy spacer- Polyether ether ketone (PEEK) and Titanium alloy (tiaLV)- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>6.Lumbar pedicle screw system-Titanium Alloy (tiaLV)- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p> <p>7.Cervical spine lateral mass/pedicle screw system-Titanium alloy (TiALV) cobalt chrome alloy (CoCr Mo)- Non sterile-Not labeled (GSR690 (E) 25/09/2014)- E&D</p>	CLAA/MD/MS /01/2017-DC	<p>Mfg Licno No.MH/102470</p> <p>Validity from 21/12/2017 to 20/12/2022</p> <p>(29-12-2017)</p> <p>(Approval for Products)</p>

178.	M/s Shreji Medicine,142/161,Platinum Industrial Park, Opp.Jaimini Park Society,SLM Mil compound,Vatva,Ahmedabad-382440,Gujarat	1.Bone screw-SS316L and titanium-5-Non sterile-E &D 2.Bone Plate locking- SS316L and titanium-5-Non sterile-E &d 3. Nail- SS316L and titanium-5-Non sterile-E &d 4.Spine- SS316L and titanium-5-Non sterile-E &d 5.Arthoscopic implants- SS316L and titanium-5-Non sterile-E &d 6.Cancellous screw (low profile)- SS316L and titanium-5-Non sterile-E &d	CLAA/MD/GU J/16/2017-DC	Mfg Licno No.G/28/1637 Validity from 03/10/2017 to 02/10/2022 (29-12-2017) (Approval for Products)
179.	M/s Angiplast Pvt.Ltd.,Plot No:-4803 and 4804,Phase- IV, GIDC,Vatva,Ahmedabad- 382445,Gujarat	1.Yankauer suction set-Yankauer handle-Rigid PVC, Connector- PVC,Tubing-PVC-Sterile-5 years from date of manufacturing – Domestic 2.Suction Catheter-Catheter-PVC,Thumb Control-Rigid PVC- Sterile-5 years from date of manufacturing –Domestic 3.Thoracic catheter with trocar-Thoracic catheter-PVC,Trocar- Aluminium,Connector for thoracic catheter-ABS,Cap for Trocar-ABS- Sterile-5 years from date of manufacturing –Domestic 4.Closed wound drainage unit-Bellow container-LDPE,Perforated catheter-PVC, Curved Needle-SS,Connecting Tube-PVC, C Clamp- HDPE, Hanger-HDPE, Below container cap-HDPE,Inner Cap-PVC- Sterile-5 years from date of manufacturing –Domestic 5.Abdominal drain kit-Bag-soft PVC Connector,Catheter-PVC,Tubing- PVC,Connector-PVC,Hanger-PVC- Sterile-5 years from date of manufacturing –Domestic	CLAA/MD/GU J/17/2010-DC	Mfg Licno No.G/28/1357 Validity from 25/11/2015 to 24/11/2020 (29-12-2017) (Approval for Additional Products)
180.	M/s Meril Life Sciences Pvt.Ltd., Survey No:- 135/2a and 135/2/B,174/2 (M1-M2,Meril Park,4 th floor),Near G.M.Bilakhia stadium,Muktanand Marg,Chala,Vapi- 396191,Valsad,Gujarat	1.Intramedullary nails-SS-Non Sterile-E &D 2. Intramedullary nails- Ti and SS-Non Sterile-E &D 3.Bone Plates-TAN and TiCP-Non Sterile-E &D 4.Bone Screws--TAN and TiCP-Non Sterile-E &D 5.Hip Acetabular System-UHMWPE,HXLPE,E-CiMa-Sterile-5-E &D 6.Bone Plates-SS316LVM-Non Sterile-E &D 7.Bone Screws- SS316LVM-Non Sterile-E &D	CLAA/MD/GU J/05/2013-DC	Mfg Licno No.G/28/1465 Validity from 01/02/2014 to 31/01/2019 (29-12-2017) (Approval for Additional Products)

181.	M/s Sahajanand LaserTechnology Ltd., A-8, GIDC Electronic Estate, Sector-25,Gandhi Nagar-382016	1.Cobalt chromium coronary stent(with delivery system)-Cobalt Chromium-Sterile and Non sterile-24 months-E &D 2. Cobalt chromium Rapamycin Eluting-Cobalt Chromium-Sterile and Non sterile-18 months-E &D 3.Coronary stent(with delivery system)- 4.Cobalt chromium rapamycin eluting coronary stent(with delivery system)-Cobalt Chromium-Sterile and Non Sterile-18 months-E &D 5.Vector-Nylon-Sterile and Non Sterile-24months- E &D	CLAA/MD/GU J/08/2009-DC	Mfg Licno No.G/28/1314 Validity from 27/08/2014 to 26/08/2019 (29-12-2017) (Approval for Additional Products)
------	---	---	-------------------------	--