

No. 6-4/2018-S/ 2152
OFFICE OF THE DIRECTOR
CENTRAL DRUGS LABORATORY,
GOVERNMENT OF INDIA
MIN. OF HEALTH & FAMILY WELFARE
3, KYD STREET, KOLKATA-16

Kolkata Dated the 25/01/19


NOTICE INVITING QUOTATION

The Director-in Charge, Central Drugs Laboratory, Kolkata invites quotation from the interested Manufacturer, Distributer or Dealer for providing and supplying of **Dehumidifier** at Central Drugs Laboratory, Kolkata.

For detailed Notice. please visit <http://www.cdsc.gov.in>

The detailed terms & condition of quotation for **Dehumidifier** is available on the above website. Filled in quotation, accepting the stipulated terms & condition may be dropped in tender box available at the wing of Jr. Administrative Officer, Central Drugs Laboratory, Kolkata. Tender / Quotation on any working day between 9:30 am. to 6:00 pm. and upto **3:00 pm on 27.02.2019**. The quotation will be opened on the same day at **3:30 pm**. in presence of interested Quotationers.

Yours faithfully


25/01/19

(Dr. Anjan Pal)
SSO (Microbiology)
for Director
Central Drugs Laboratory, Kolkata.

No. 6-4/2018-S/
OFFICE OF THE DIRECTOR
CENTRAL DRUGS LABORATORY,
GOVERNMENT OF INDIA
MIN. OF HEALTH & FAMILY WELFARE
3, KYD STREET, KOLKATA-16

Kolkata Dated the


NOTICE INVITING QUOTATION

Quotation is hereby invited from the interested and bona fide Dealer, Distributors, Retailers for providing and supplying of **DEHUMIDIFIER** under the following specification, terms and conditions:

Name of the Item	Qty.
<p><u>DEHUMIDIFIER</u></p> <p><u>Specification</u></p> <ol style="list-style-type: none"> 1. Type :Desiccant Type Table Top. 2. Capacity :600±5% CMH or equivalent CFM 3. Power :230 V/1Ø/ 50 Hz 4. Processor :Should have Microprocessor based controlling system with LED / LCD display 5. Should have automatic Humidistat with display 6. All the units preferably CNC fabricated with powder coated stainless steel body. 7. Rotor shall be non flammable and should have robust internal structures with stainless steel perimeter flange to extend media and self life and for durability and easy serviceability. 8. Rotor media shall have in-situ synthesized high performance metal silicates. 9. The system should have Differential Air Pressure switch for providing reactivation air flow. 10. Process and reactivation air flow sectors shall be insulated and have PTFE bonded seal design to minimize air leakage. 11. Dehumidifier should have independent blower and motor for both process and reactivation air flow system. 12. Dehumidifier should have electrical interlocking of fan motors, heaters and rotor drive. 13. Working temperature range from 10°C to 40°C. 14. Dehumidifier shall have reactivation cool down safety device. 15. Dehumidifier Processor and Reactivation unit should have high capacity metallic mesh filter preferably of EU5 grade. 16. Dehumidifier should have additional cooling thermostat and cut-out for safety measure. 17. Dimension of Dehumidifier preferably within : 900mm (L) x 650mm (D) x 500mm (H) 18. Dehumidifier shall have CE and ISO certifications. 19. Manufacturer must have adequate number of Service Engineers to provide prompt service after sales. 	<p>1</p> <p>No.</p>

N.B: Applicable Taxes and duties if any, which may be levied by the Government, should be shown separately.

Yours faithfully


 Dr. Anjan Pal, 25/01/19
 SSO (Microbiology)
 for Director

Central Drugs Laboratory, Kolkata.

No. 6-4/2018-S/
OFFICE OF THE DIRECTOR
CENTRAL DRUGS LABORATORY,
GOVERNMENT OF INDIA
MIN. OF HEALTH & FAMILY WELFARE
3, KYD STREET, KOLKATA-16

Kolkata Dated the

NOTICE INVITING QUOTATION

Quotation is hereby invited from the interested and bona fide Dealer, Distributors, Retailers for Providing and supplying of **DEHUMIDIFIER** under the following specification, terms and conditions:

Name of the Item	Qty.
<p><u>DEHUMIDIFIER</u></p> <p><u>Specification</u></p> <ol style="list-style-type: none"> 1. Type : Desiccant Type Bench Top. 2. Capacity : 175±5% CMH or equivalent CFM 3. Power : 230 V/1Ø/ 50 Hz 4. Processor : Should have Microprocessor based controlling system with LED / LCD display 5. Should have automatic Humidistat with display 6. All the units preferably CNC fabricated with powder coated stainless steel body. 7. Rotor shall be non flammable and should have robust internal structures with stainless steel perimeter flange to extend media and self life and for durability and easy serviceability. 8. Rotor media shall have in-situ synthesized high performance metal silicates. 9. The system should have Differential Air Pressure switch for providing reactivation air flow. 10. Process and reactivation air flow sectors shall be insulated and have PTFE bonded seal design to minimize air leakage. 11. Dehumidifier should have independent blower and motor for both process and reactivation air flow system. 12. Dehumidifier should have electrical interlocking of fan motors, heaters and rotor drive. 13. Working temperature range from 10°C to 40°C. 14. Dehumidifier shall have reactivation cool down safety device. 15. Dehumidifier Processor and Reactivation unit should have high capacity metallic mesh filter preferably of EU5 grade. 16. Dehumidifier should have additional cooling thermostat and cut-out for safety measure. 17. Dimension of Dehumidifier preferably within : 700mm (L) x 500mm (D) x 400mm (H) 18. Dehumidifier shall have CE and ISO certifications. 19. Manufacturer must have adequate number of Service Engineers to provide prompt service after sales. 	<p>2</p> <p>Nos.</p>

N.B: Applicable Taxes and duties if any, which may be levied by the Government, should be shown separately.

Yours faithfully

Anjan Pal
Dr. Anjan Pal, 25/01/19

SSO (Microbiology)


for Director

Central Drugs Laboratory, Kolkata.

Terms & Conditions. :

1. The quotation/tender should be submitted in sealed covers. Open quotation/tender will be summarily rejected.
2. The quotation should be valid for a period of 90 (ninety) days from the date of opening of the tender / quotation.
3. The Director-in-Charge, CDL, Kolkata reserves the right to reject the quotation and/or supply order at any stage without assigning any reasons thereof.
4. The quoted rate should be all inclusive and nothing will be paid extra for transportation, loading, unloading at Central Drugs Laboratory, Kolkata.
5. The Quotationer/supplier should supply of **Dehumidifier** in good condition & strictly as per specification within thirty (30) days from the date of issue of supply order.
6. The person duly authorized by the quotationer will be allowed to deliver **Dehumidifier** at Central Drugs Laboratory, Kolkata on any working day within the time schedule.
7. No "Way Bill" in connection with supply of **Dehumidifier** will be issued from Central Drugs Laboratory, Kolkata.
8. **The quotationer should clearly mention the name of manufacturer and Product No.**
9. That on successful completion of supply of **Dehumidifier** as per specification, make and Catalogue Number to Store Section, CDL, Kolkata, the quotationer/supplier should submit bill in triplicate along with mandate form to release e-payment through Pay & Accounts Office, Min. of Health & F.W. Kolkata.
10. **Only manufacturing firms / authorized distributors are eligible to participate.**
11. The Quotationer should clearly quote their rates for both in figures and in words **in net price (including any other taxes/discounts etc.) failing which your quotation will not be entertained.** Any overwriting in rates or any change and/or alteration without attestation is liable to culminate in the outright rejection of the quotation. In case there is any error between the amount specified in figures and words, the lower of the two will be taken for consideration.
12. The Quotation duly completed as per the instruction given herein and accepting the same may be dropped in the "Tender Box" on any working days between **11.00 AM to 3.00 P.M. up to 27.02.2019** at the wing of Jr. Administrative Officer, Central Drugs Laboratory, Government of India, Ministry of Health & Family Welfare, 3, Kyd Street, Kolkata-700016.
13. Late quotation i.e. quotation received after the specified date and time of receipt will not be considered in any circumstances.
14. The Quotation so received will be opened on the same day i.e. on **27.02.2019 at 3.30 P.M.** in presence of interested Quotationer(s) or their authorized representative.

Yours faithfully


Dr. Anjan Pal,
SSO (Microbiology)
for Director
Central Drugs Laboratory, Kolkata.